Welcome!				
This document contains information about all the classes in the Python-POD_API.	Each class has its own sheet which describes the class's origin, parent and child classes, imports, global and instance variables, and all methods with details about its parameters, returns, and exceptions.			
Classes				
POD_8206HR				
POD_8401HR				
POD_Basics				
POD_Packets				
POD_Commands				
COM_io				

Class						
Name	File	Description	Parent	Child	Author	
POD_8206HR	PodDevice_ 8206HR.py	Handles communication using an 8206HR POD device.	POD_Basics	N/A	Thresa Kelly	
Imports						
Name	Origin	Description	From			
POD_Basics	Local	For inheritance	BasicPodProtocol			
POD_Packets	Local	For handling POD packets	PodPacketHandling			
POD_Commands	Local	For command constants	PodCommands			
Variables						
Name	Scope	Description	Value			
B4LENGTH	Class	Constant containing the number of bytes for a full Binary4 packet	16			
B4BINARYLENGTH	Class	Constant containing the number of bytes for binary data in a Binary4 packet	8			
_preampGain	Instance	Preamplifier gain	10 or 100			
Methods						
Name	Туре	Description	Parameter Name	Parameter Purpose	Return	Exception
init	Dunder	Runs when an instance is constructed. It runs the parent's initialization. Then it updates the _commands	port	String of the serial port to be opened. Used when initializing the COM_io instance.	N/A	N/A
	Dunder	to contain the appropriate commands for an 8306HR POD device.	baudrate=9600	Integer baud rate of the opened serial port. Used when initializing the COM_io instance.		177
UnpackPODpacket_Binary	Static	Overwrites the parent's method. Separates the components of a binary4 packet into a dictionary.	msg	Bytes string containing a complete binary4 Pod packet: STX (1 byte) + command (4 bytes) + packet number (1 bytes) + TTL (1 byte) + ch0 (2 bytes) + ch1 (2 bytes) + ch2 (2 bytes) + checksum (2 bytes) + ETX (1 byte)	A dictionary containing 'Command Number', 'Packet #', 'TTL', 'Ch0', 'Ch1', and 'Ch2' in bytes.	An exception is raised if (1) the packet does not have the minimum number of bytes, (2) does not begin with STX, or (3) does not end with ETX.
TranslatePODpacket_Binary	Static	Overwrites the parent's method. Unpacks the binary4 POD packet and converts the values of the ASCII-encoded bytes into integer values and the values of binary-encoded bytes into integers. Channel values are given in Volts.	msg	Bytes string containing a complete binary4 Pod packet: STX (1 byte) + command (4 bytes) + packet number (1 bytes) + TTL (1 byte) + ch0 (2 bytes) + ch1 (2 bytes) + ch2 (2 bytes) + checksum (2 bytes) + ETX (1 byte)	A dictionary containing 'Command Number', 'Packet #, 'TTL', 'Ch0', 'Ch1', and 'Ch2' as numbers.	N/A
TranslatePODpacket	Instance	Overwrites the parent's method. Determines if the packet is standard or binary, and translates accordingly. Adds a check for the 'GET TTL PORT' command.		Bytes string containing either a standard or binary packet	A dictionary containing the unpacked message in numbers	N/A
_TranslateTTLbyte_ASCII	Static	Separates the bits of each TTL (0-3) from a byte.	ttlByte	One Byte string for the TTL (ASCII encoded)	Tuple of the integer TTLs (0-3). 1 when input, 0 when output.	N/A
_TranslateTTLbyte_Binary	Static	Separates the bits of each TTL (0-3) from a byte.	ttlByte	One Byte string for the TTL (binary encoded)	Tuple of the integer TTLs (0-3). 1 when input, 0 when output.	N/A
_BinaryBytesToVoltage	Instance	Converts a binary bytes value read from POD device and converts it to the real voltage value at the preamplifier input	value	Bytes string containing voltage measurement	A number containing the voltage in Volts [V].	N/A
Read Binary	Instance	After receiving the prePacket, it reads the 8 bytes(TTL+channels) and then reads to ETX	prePacket	Bytes string containing the beginning of a POD packet: STX (1 byte) + command number (4 bytes)	Bute string for a hinary4 POD nacket	N/Δ
	motanice	(checksum+ETX).	validateChecksum=True	Set to True to validate the checksum. Set to False to skip validation	Byte string for a binary4 POD packet.	N/A

Class						
Name	File	Description	Parent	Child	Author	
name	-	Description	Parent		Autnor	
POD_8401HR	8401HR.py	Handles communication using an 8401HR POD device.	POD_Basics	N/A	Thresa Kelly	
Imports						
Name	Origin	Description	From			
POD_Basics	Local	For inheritance	BasicPodProtocol			
POD_Packets	Local	For handling POD packets	PodPacketHandling			
POD_Commands	Local	For command constants	PodCommands			
Variables						
Name	Scope	Description	Value			
B5LENGTH	Class	number of bytes for a Binary 5 packet	31			
B5BINARYLENGTH	Class	number of binary bytes for a Binary 5 packet	23			
_channelMap	Instance	Dictionary of the channel lables	Set by _GetChannelMapping(). Dictionary keys are ['A','B','C','D']			
_ssGain	Instance	Dictionary of the second stage gain for all four channels	1, 5, or None. Dictionary keys are ['A','B','C','D']			
_preampGain	Instance	Dictionary of the preamplifier gain for all four channels.	10, 100, or None. Dictionary keys are ['A','B','C','D']			
Methods						
Name	Туре	Description	Parameter Name	Parameter Purpose	Return	Exception
		·	port	String of the serial port to be opened. Used when initializing the COM io instance.		·
		Runs when an instance is constructed. It runs the parent's initialization. Then it updates the _commands to contain the appropriate commands for an 8401HR POD device. Sets the _channelMap, _ssGain, and _preampGain.	deviceName	String of the corresponding device/sensor name		An exception is raised if (1) the ssGain or preampGain have improper keys, (2) the device/sensor does not exist, (3) the ssGain was given bad values, or (4) the preampGain was given bad values
	Dunder		ssGain={'A':None,'B':None,'C':None			
init			,'D':None}	1 - 1	N/A	
			preampGain={'A':None,'B':None,'C': None,'D':None}	Dictionary of the preamplifier gain		
			baudrate=9600	Integer baud rate of the opened serial port. Used when initializing the COM_io instance.		
UnpackPODpacket_Binary	Static	Overwrites the parent's method. Separates the components of a binary5 packet into a dictionary.	msg	Bytes string containing a complete binary5 Pod packet: STX (1 byte) + command (4) + packet number (1) + status (1) + channels (9) + analog inputs (12) + checksum (2) + ETX (1)	A dictionary containing 'Command Number', 'Packet #,' 'Status', 'Channels', 'Analog EXT0', 'Analog EXT1', 'Analog TTL1', 'Analog TTL2', 'Analog TTL3', 'Analog TTL4', in bytes.	An exception is raised if (1) the packet does not have the minimum number of bytes, (2) does not begin with STX, or (3) does not end with ETX.
TranslatePODpacket_Binary	Instance	Overwrites the parent's method. Unpacks the binary5 POD packet and converts the values of the ASCII-encoded bytes into integer values and the values of binary-encoded bytes into integers. The channels and analogs are converted to volts (V).	msg	Bytes string containing a complete binary5 Pod packet: STX (1 byte) + command (4) + packet number (1) + status (1) + channels (9) + analog inputs (12) + checksum (2) + ETX (1)	A dictionary containing 'Command Number', 'Packet #', 'Status', 'CH3', 'CH2', 'CH1', 'CH0', 'Analog EXT0', 'Analog EXT1', 'Analog TTL1', 'Analog TTL2', 'Analog TTL3', 'Analog TTL4', as numbers.	N/A
GetChannelMapping	Static	Get the channel mapping (channel labels for A,B,C,D) for a given device.	deviceName	String for the device/sensor name.	Dictionary with keys A,B,C,D with values of the channel names. Returns None if the device name does not exist.	N/A
			value	Value to be converted to voltage		
Voltage PrimaryChannels	Static	Converts a value to a voltage for a primary channel.	ssGain=None	Second stage gain	Number of the voltage in volts [V]. Returns value if	N/A
		and the second s	PreampGain=None	Preamplifier gain	no gain is given (no-connect).	
			value	Value to be converted to voltage		
		Converts a value to a voltage for an EEC/EMC	ssGain	Second stage gain		
_Voltage_PrimaryChannels_EEGEMG	Static	Converts a value to a voltage for an EEG/EMG primary channel.	PreampGain	Preamplifier gain	Number of the voltage in volts [V].	N/A
		Converts a value to a voltage for a biosensor primary	value	Value to be converted to voltage		
_Voltage_PrimaryChannels_Biosensor	Static	channel.	ssGain	Second stage gain	Number of the voltage in volts [V].	N/A
_Voltage_SecondaryChannels	Static	Converts a value to a voltage for a secondary channel.	value	Value to be converted to voltage	Number of the voltage in volts [V].	N/A
		After receiving the prePacket, it reads the 23 bytes	prePacket	Bytes string containing the beginning of a POD packet: STX (1 byte) + command number (4 bytes)		
_Read_Binary	Instance	(binary data) and then reads to ETX (checksum+ETX).	validateChecksum=True	Set to True to validate the checksum. Set to False to skip validation	Byte string for a binary5 POD packet.	N/A

Class							
Name	File	Description	Parent	Child	Author		
POD_Packets	PodPacket Handling.py	Collection of methods for creating and interpreting POD packets	N/A	N/A	Thresa Kelly		
Imports							
Name	Origin	Description	From				
N/A	N/A	N/A	N/A				
Methods							
Name	Туре	Description	Parameter Name	Parameter Purpose	Return	Exception	
STX	Static	Get STX in bytes. STX marks the starting byte of a POD Packet	N/A	N/A	Bytes for STX (0x02)	N/A	
ETX	Static	Get ETXin bytes. ETX marks the end byte of a POD Packet	N/A	N/A	Bytes for ETX(0x03)	N/A	
IntToAsciiBytes	Static	0:	Converts an integer value into ASCII-encoded bytes. First, it converts the integer value into a usable uppercase hexadecimal string. Then it converts the ASCII code for each character into bytes. Lastly, it	value	Integer value to be converted into ASCII-encoded bytes	Bytes that are ASCII-encoded conversions of the	N/A
	J. J	ensures that the final message is the desired length. Example: if value=2 and numBytes=4, the returned ASCII will show b'0002', which is '0x30 0x30 0x30 0x30 0x32' in bytes.	numBytes	Number bytes to be the length of the ASCII-encoded message.	value parameter.	N/A	
AsciiBytesToInt	Static	Converts a ASCII-encoded bytes message into an integer. It does this using a base-16 conversion.	msg_b	Bytes message to be converted to an integer. The bytes must be base-16 or the conversion will fail.	Integer result from the ASCII-encoded byte conversion.	N/A	
BinaryBytesToInt Static		: Converts binary-encoded bytes into an integer	msg	Bytes message holding binary information to be converted into an integer.	Integer result from the binary-encoded bytes message.	N/A	
	Static		byteorder='big'	Ordering of bytes. 'big' for big endian and 'little' for little endian.			
			signed=False	Boolean flag to mark if the msg is signed (True) or unsigned (False)			
		Converts a specific bit range in an ASCII-encoded bytes object to an integer.	msg	Bytes message holding binary information to be converted into an integer.	Integer result from the ASCII-encoded bytes message in a given bit range.	e N/A	
ASCIIbytesToInt_Split	Static		keepTopBits	Integer position of the msb of desired bit range			
			cutBottomBits	Integer number of lsb to remove			
			msg	Bytes message holding binary information to be converted into an integer.		9 N/A	
			keepTopBits	Integer position of the msb of desired bit range			
BinaryBytesToInt_Split	Static	Converts a specific bit range in a binary-encoded	cutBottomBits	Integer number of lsb to remove	Integer result from the binary-encoded bytes message		
binary bytes rount_opin	Otatio	bytes object to an integer	byteorder='big'	Ordering of bytes. 'big' for big endian and 'little' for little endian.	in a given bit range.	IWA	
			signed=False	Boolean flag to mark if the msg is signed (True) or unsigned (False)			
Checksum		Calculates the checksum of a given bytes message. This is achieved by summing each byte in the message, inverting, and taking the last byte.	bytesIn	Bytes message containing POD packet data	Two ASCII-encoded bytes containing the checksum for bytesIn	N/A	
BuildPODpacket_Standard		Builds a standard POD packet STX (1 byte) + command number (4 bytes) + optional packet (? bytes) + checksum (2 bytes) + ETX (1 bytes) as	commandNumber	Integer representing the command number. This will be converted into a 4 byte long ASCII-encoded bytes string.	Bytes string of a complete standard POD packet	N/A	
		bytes.	payload=None	bytes string containing the payload			
PayloadToBytes	Static	Converts a payload into a bytes string	payload	Integer, bytes, or tuple containing the payload	Bytes string of the payload	Raises an Exception when the payload argument is	
rayioau iobytes Static		auc Converts a payload into a bytes string	argSizes	Tuple of the argument sizes		an incorrect type or formatted incorrectly.	

Class	FII.	Daniel dan	B	Object	A41	
Name	File	Description	Parent	Child	Author	
POD_Basics	BasicPodPr otocol.py	Handle basic communication with a POD device, including reading and writing packets and packet interpretation.	N/A	POD_8206HR	Thresa Kelly	
Imports						
Name	Origin	Description	From			
COM_io	Local	For opening and connecting serial COM ports	SerialCommunication			
POD_Packets	Local	For handling POD packets	PodPacketHandling			
POD_Commands	Local	Used to contain all POD commands in the class instance	PodCommands			
Variables						
Name	Scope	Description	Value			
numPod	Class	Integer equal to the number of POD_Basics class instances. Incremented on construction and decremented on destruction	0			
MINSTANDARDLENGTH	Class	integer minimum number of bytes in a standard POD packet	8			
MINBINARYLENGTH	Class	integer minimum number of bytes in a binary POD packet	15			
_port	Instance	Open serial port via COM_io class instance	COM_io			
_commands	Instance	Command handler POD_Commands class instance	POD_Commands			
Methods						
Name	Туре	Description	Parameter Name	Parameter Purpose	Return	Exception
init	Dunder	Runs when an instance of POD_Basics is constructed. It initializes the instance variable for the COM port communication (port) and for the	port	String of the serial port to be opened. Used when initializing the COM_io instance.	N/A	N/A
		command handler (_commands). It also increments the POD device counter (NUMPOD).	baudrate=9600	Integer baud rate of the opened serial port. Used when initializing the COM_io instance.		
del	Dunder	Runs when an instance is destructed. It decrements the POD device counter (NUMPOD)	N/A	N/A	N/A	N/A
GetNumberOfPODDevices	Static	Get the POD device counter	N/A	N/A	Integer of the number of class instances (NUMPOD).	N/A
UnpackPODpacket_Standard	Static	Converts a standard POD packet into a dictionary containing the command number and payload (if applicable) in bytes.	msg		A dictionary containing the POD packet's 'Command Number' and 'Payload' (if applicable) in bytes.	An exception is raised if (1) the msg does not have the minimum number of bytes in a standard pod packet, (2) does not begin with STX, and (3) does not end with ETX.
UnpackPODpacket_Binary	Static	Converts a variable-length binary packet into a dictionary containing the command number, binary packet length, and binary data in bytes.	msg	Bytes message containing a variable-length POD packet: STX (1 byte) + command number (4 bytes) + length of binary (4 bytes) + checksum (2 bytes) + ETX (1 bytes) + binary (LENGTH bytes) + checksum (2 bytes) + ETX (1 bytes)	A dictionary containing the 'Command Number', 'Binary Packet Length', and 'Binary Data' in bytes.	An exception is raised if (1) the msg does not have the minimum number of bytes in a standard pod packet, (2) does not begin with STX, (3) does not end with ETX, and (4) does not have an ETX after standard packet.
TranslatePODpacket_Standard	Instance	Unpacks the standard POD packet and converts the ASCII-encoded bytes values into integer values.	msg	Bytes message containing a standard POD packet	A dictionary containing the POD packet's 'Command Number' and 'Payload' (if applicable) in integers.	N/A
TranslatePODpacket_Binary	Static	Unpacks the variable-length binary POD packet and converts the values of the ASCII-encoded bytes into integer values and leaves the binary-encoded bytes as is.	msg	bytes message containing a variable-length POD	A dictionary containing the 'Command Number' and 'Binary Packet Length' in integers, and 'Binary Data' in bytes.	N/A
_ValidateChecksum	Static	Validates the checksum of a given POD packet. The checksum is valid if the calculated checksum from the data matches the checksum written in the packet.	msg	Bytes message containing a POD packet: STX (1 bytes) + data (? bytes) + checksum (2 bytes) + ETX (1 byte).	Returns True if the checksum is correct, false otherwise.	An exception is raised if the msg does not begin with STX or end with ETX.
GetDeviceCommands	Instance	Gets the dictionary containing the class instance's available POD commands.	N/A		Dictionary containing the available commands and their information. Formatted as key(command number): value([command name, number of argument ASCII bytes, number of return bytes, binary flag])	N/A
SetBaudrateOfDevice	Instance	If the port is open, it will change the baud rate to the parameter's value	baudrate	Integer baud rate to set for the open serial port.	True if successful at setting the baud rate, false otherwise	N/A
UnpackPODpacket	Static	Determines if the packet is standard or binary, and unpacks accordingly.	msg		A dictionary containing the unpacked message in bytes	N/A
TranslatePODpacket	Instance	Determines if the packet is standard or binary, and translates accordingly.	msg	packet	A dictionary containing the unpacked message in numbers	N/A
			cmd	An integer representing the command number.	Bytes string containing a POD packet beginning with	
WriteRead	Instance	Writes a command with optional payload to POD device, then reads (once) the device response.	payload=None	None when there is no payload. If there is a payload, set to an integer value or a bytes string.	STX and ending with ETX. This may be a standard packet, binary packet, or an unformatted packet	N/A
			validateChecksum=True	Set to True to validate the checksum. Set to False to skip validation	(STX+something+ETX).	
		Ruilds a POD nacket and writes it to a POD device via	cmd	An integer representing the command number.		An exception is raised if (1) the command does not

GetPODpacket	Instance	Dullus a P-U-Dacket and writes it to a P-U-Device vial COM port. If an integer payload is give, the method will convert it into a bytes string of the length expected by the command. If a bytes payload is given, it must be the correct length.		None when there is no payload. If there is a payload, set to an integer value, bytes string, or tuple	Returns the bytes string of the POD packet.	exist for the instance, (2) a payload is not given when the command expects one, (3) the payload (given in bytes) is the size not expected by the command, or (4) the payload is given as a type other than integer or bytes.
WritePacket	Instance	Builds a POD packet and writes it to the POD device.	cmd	An integer representing the command number.	Returns the bytes string that was written to the POD	N/A
	mounise	Dallac a 1 GS pasies and miles it is allo 1 GS assists.	payload=None	None when there is no payload. If there is a payload, set to an integer value, bytes string, or tuple	device	
ReadPODpacket	Instance	Reads a complete POD packet, either in standard or binary format, beginning with STX and ending with ETX. Reads first STX and then starts recursion.	validateChecksum=True	Set to True to validate the checksum. Set to False to skip validation	Bytes string containing a POD packet beginning with STX and ending with ETX. This may be a standard packet, binary packet, or an unformatted packet (STX+something+ETX).	N/A
_ReadPODpacket_Recursive	Instance	Reads the command number. If the command number ends in ETX, the packet is returned. Next, it checks if the command is allowed. Then, it checks if the command is standard or binary and reads accordingly, then returns the packet.	validateChecksum=True	Set to True to validate the checksum. Set to False to skip validation	Bytes string containing a POD packet beginning with STX and ending with ETX. This may be a standard packet, binary packet, or an unformatted packet (STX+something+ETX).	N/A
_Read_GetCommand	Instance	Reads one byte at a time up to 4 bytes to get the ASCII-encoded bytes command number. For each byte read, it can (1) start the recursion over if an STX is found, (2) returns if ETX is found, or (3) continue building the command number.	validateChecksum=True	Set to True to validate the checksum. Set to False to skip validation	4 byte long string containing the ASCII-encoded command number.	An exception is raised if the command number is not allowed for the POD device
_Read_ToETX	Instance	Reads one byte at a time until an ETX is found. It will restart the recursive read if an STX is found anywhere.	validateChecksum=True	Set to True to validate the checksum. Set to False to skip validation	Bytes string ending with ETX	N/A
Read Standard	Instance	Instance Reads the payload, checksum, and ETX. Then it builds the complete standard POD packet in bytes.	prePacket	Bytes string containing the beginning of a POD packet: STX (1 byte) + command number (4 bytes)	Dutan atrian for a complete standard DOD market	An exception is raised if the checksum is invalid (only
_Read_Standard	instance		validateChecksum=True	Set to True to validate the checksum. Set to False to skip validation	Bytes string for a complete standard POD packet	if validateChecksum=True)
		(prePacket+payload+checksum+ETX). Then it Instance determines how long the binary packet is from the payload of the standard POD packet and reads that	prePacket	Bytes string containing the beginning of a POD packet: STX (1 byte) + command number (4 bytes)	Bytes string for a variable-length binary POD packet	An exception is raised if the checksum is invalid (only if validateChecksum=True)
_Read_Binary	modrice		validateChecksum=True	Set to True to validate the checksum. Set to False to skip validation		

Class						
Class	Eu.	D	B4	Object	A45	
Name	File	Description	Parent	Child	Author	
POD_Commands	PodComma nds.py	Manages a dictionary containing available commands for a POD device.	N/A	N/A	Thresa Kelly	
Imports						
Name	Origin	Description	From			
N/A	N/A	N/A	N/A			
Variables						
Name	Scope	Description	Value			
NAME	Class	index key for the command name forcommands list values	0			
ARGUMENTS	Class	index key for the number of bytes in an argument forcommands list values	1			
RETURNS	Class	index key for the number of bytes in the return forcommands list values	2			
BINARY	Class	index key for the binary flag forcommands list values	3			
NOVALUE	Class	Integer used to mark when a list item incommands means 'no value' or undefined.	-1			
_U8	Class		2			
U16	Class	Number of bytes for an unsigned 16-bit value	4			
commands	Instance	Dictionary containing the available commands for a POD device. Each entry is formatted as { key(command number) : value([command name, number of argument ASCII bytes, number of return bytes, binary flag) }	{ key : value }			
Methods						
Name	Туре	Description	Parameter Name	Parameter Purpose	Return	Exception
init	Dunder	Runs whan an instance is constructed. It sents the commands dictionary to the basic command set.	N/A	N/A	N/A	N/A
NoValue	Static	Gets value ofNOVALUE	N/A	N/A	Value ofNOVALUE	N/A
U8	Static	Gets value ofU8	N/A	N/A	Value ofU8	N/A
U16	Static	Gets value ofU16	N/A	N/A	Value ofU16	N/A
GetBasicCommands	Static	Creates a dictionary containing the basic POD command set (0,1,2,3,4,5,6,7,8,9,10,11,12)	N/A	N/A	N/A	N/A
GetBasicCommands GetCommands	Static	command set (0,1,2,3,4,5,6,7,8,9,10,11,12)	N/A	N/A N/A		N/A N/A
		command set (0,1,2,3,4,5,6,7,8,9,10,11,12) Gets the contents of the current command dictionary			N/A	
GetCommands	Instance	command set (0,1,2,3,4,5,6,7,8,9,10,11,12) Gets the contents of the current command dictionary (_commands) Sets the current commands (_commands) to the	N/A	N/A	N/A	N/A
GetCommands	Instance	command set (0,1,2,3,4,5,6,7,8,9,10,11,12) Gets the contents of the current command dictionary (_commands) Sets the current commands (_commands) to the basic POD command set.	N/A N/A	N/A N/A	N/A N/A	N/A
GetCommands RestoreBasicCommands	Instance	command set (0,1,2,3,4,5,6,7,8,9,10,11,12) Gets the contents of the current command dictionary (_commands) Sets the current commands (_commands) to the basic POD command set. Adds a command entry to the current commands	N/A N/A commandNumber	N/A N/A Integer of the command number	N/A N/A True if the command was successfully added, False if	N/A N/A
GetCommands	Instance	command set (0,1,2,3,4,5,6,7,8,9,10,11,12) Gets the contents of the current command dictionary (_commands) Sets the current commands (_commands) to the basic POD command set.	N/A N/A commandNumber commandName	N/A N/A Integer of the command number String of the command's name	N/A N/A	N/A N/A
GetCommands RestoreBasicCommands	Instance	command set (0,1,2,3,4,5,6,7,8,9,10,11,12) Gets the contents of the current command dictionary (_commands) Sets the current commands (_commands) to the basic POD command set. Adds a command entry to the current commands dictionary (_commands) if the command does not	N/A N/A commandNumber commandName argumentBytes	N/A N/A Integer of the command number String of the command's name Integer of the number of bytes in the argument	N/A N/A True if the command was successfully added, False if the command could not be added because t already exists.	N/A N/A
GetCommands RestoreBasicCommands	Instance	command set (0,1,2,3,4,5,6,7,8,9,10,11,12) Gets the contents of the current command dictionary (_commands) Sets the current commands (_commands) to the basic POD command set. Adds a command entry to the current commands dictionary (_commands) if the command does not	N/A N/A commandNumber commandName argumentBytes returnBytes	N/A N/A Integer of the command number String of the command's name Integer of the number of bytes in the argument Integer of the number of bytes in the return Boolean flag to mark if the command is binary (True)	N/A N/A True if the command was successfully added, False if the command could not be added because t already	N/A N/A
GetCommands RestoreBasicCommands AddCommand	Instance Instance Instance	command set (0,1,2,3,4,5,6,7,8,9,10,11,12) Gets the contents of the current command dictionary (_commands) Sets the current commands (_commands) to the basic POD command set. Adds a command entry to the current commands dictionary (_commands) if the command does not exist Removes the entry for a given command in	N/A N/A commandNumber commandName argumentBytes returnBytes isBinary	N/A N/A Integer of the command number String of the command's name Integer of the number of bytes in the argument Integer of the number of bytes in the return Boolean flag to mark if the command is binary (True) or standard (False)	N/A N/A True if the command was successfully added, False if the command could not be added because t already exists. True if the command was successfully removed, False if the command does not exist.	N/A N/A
GetCommands RestoreBasicCommands AddCommand RemoveCommand	Instance Instance Instance Instance	command set (0,1,2,3,4,5,6,7,8,9,10,11,12) Gets the contents of the current command dictionary (_commands) Sets the current commands (_commands) to the basic POD command set. Adds a command entry to the current commands dictionary (_commands) if the command does not exist Removes the entry for a given command in _commands dictionary. Gets the command number from the command	N/A N/A commandNumber commandName argumentBytes returnBytes isBinary cmd name	N/A N/A Integer of the command number String of the command's name Integer of the number of bytes in the argument Integer of the number of bytes in the return Boolean flag to mark if the command is binary (True) or standard (False) integer command number or string command name.	N/A N/A True if the command was successfully added, False if the command could not be added because t already exists. True if the command was successfully removed, False if the command does not exist. Integer representing the command number. If the	N/A N/A N/A N/A N/A
GetCommands RestoreBasicCommands AddCommand RemoveCommand CommandNumberFromName	Instance Instance Instance Instance Instance	command set (0,1,2,3,4,5,6,7,8,9,10,11,12) Gets the contents of the current command dictionary (_commands) Sets the current commands (_commands) to the basic POD command set. Adds a command entry to the current commands dictionary (_commands) if the command does not exist Removes the entry for a given command incommands dictionary. Gets the command number from the command dictionary using the command's name Gets the tuple for the number of bytes in the argument	N/A N/A commandNumber commandName argumentBytes returnBytes isBinary cmd name cmd	N/A N/A Integer of the command number String of the command's name Integer of the number of bytes in the argument Integer of the number of bytes in the return Boolean flag to mark if the command is binary (True) or standard (False) integer command number or string command name. string of the command's name	N/A N/A True if the command was successfully added, False if the command could not be added because t already exists. True if the command was successfully removed, False if the command does not exist. Integer representing the command number. If the command could not be found, return None. Tuple representing the number of bytes in the argument for cmd. If the command could not be found, return None. Tuple representing the number of bytes in the return None.	N/A N/A N/A N/A N/A
GetCommands RestoreBasicCommands AddCommand RemoveCommand CommandNumberFromName ArgumentBytes	Instance Instance Instance Instance Instance Instance	command set (0,1,2,3,4,5,6,7,8,9,10,11,12) Gets the contents of the current command dictionary (_commands) Sets the current commands (_commands) to the basic POD command set. Adds a command entry to the current commands dictionary (_commands) if the command does not exist Removes the entry for a given command incommands dictionary. Gets the command number from the command dictionary using the command's name Gets the tuple for the number of bytes in the argument for a given command. Gets the tuple for the number of bytes in the return for	N/A N/A commandNumber commandName argumentBytes returnBytes isBinary cmd name cmd	N/A N/A Integer of the command number String of the command's name Integer of the number of bytes in the argument Integer of the number of bytes in the return Boolean flag to mark if the command is binary (True) or standard (False) integer command number or string command name. string of the command's name integer command number or string command name.	N/A N/A True if the command was successfully added, False if the command could not be added because t already exists. True if the command was successfully removed, False if the command does not exist. Integer representing the command number. If the command could not be found, return None. Tuple representing the number of bytes in the argument for cmd. If the command could not be found, return None. Tuple representing the number of bytes in the return for cmd. If the command could not be found, return None.	N/A N/A N/A N/A N/A N/A N/A N/A
GetCommands RestoreBasicCommands AddCommand RemoveCommand CommandNumberFromName ArgumentBytes ReturnBytes	Instance Instance Instance Instance Instance Instance Instance	command set (0,1,2,3,4,5,6,7,8,9,10,11,12) Gets the contents of the current command dictionary (_commands) Sets the current commands (_commands) to the basic POD command set. Adds a command entry to the current commands dictionary (_commands) if the command does not exist Removes the entry for a given command in _commands dictionary. Gets the command number from the command dictionary using the command's name Gets the tuple for the number of bytes in the argument for a given command. Gets the tuple for the number of bytes in the return for a given command.	N/A N/A commandNumber commandName argumentBytes returnBytes isBinary cmd name cmd	N/A N/A Integer of the command number String of the command's name Integer of the number of bytes in the argument Integer of the number of bytes in the return Boolean flag to mark if the command is binary (True) or standard (False) integer command number or string command name. string of the command's name integer command number or string command name. integer command number or string command name.	N/A N/A True if the command was successfully added, False if the command could not be added because t already exists. True if the command was successfully removed, False if the command does not exist. Integer representing the command number. If the command could not be found, return None. Tuple representing the number of bytes in the argument for cmd. If the command could not be found, return None. Tuple representing the number of bytes in the argument for cmd. If the command could not be found, return None. Boolean flag that is True if the command is binary and False if standard. If the command could not be found, return None.	N/A N/A N/A N/A N/A N/A N/A N/A

Class							
Name	File	Description	Parent	Child	Author		
COM_io	SerialCommu nication.py	Handle serial communication (read/write) using COM ports.	N/A	N/A	Thresa Kelly		
Imports							
Name	Origin	Description	From				
serial.tools.list_ports	Enviornment	For accessing the COM ports on the computer	N/A				
Variables							
Name	Scope	Description	Value				
serialInst	Instance	Serial object to set the port and baud rate to. It can be opened or closed.	Serial				
Methods							
Name	Туре	Description	Parameter Name	Parameter Purpose	Return	Exception	
GetCOMportsList	Static	Finds all the available COM ports on the user's computer and appends them to an accessible list.	N/A	N/A	List containing the names of available COM ports	N/A	
init	Dunder		port	String of the serial port to be opened.	N/A	N/A	
init	Dunder	serialInst to a given COM port with a set baudrate.	baudrate=9600	Integer baud rate of the opened serial port.	IN/A	IV/A	
del	Dunder	Runs when the object is destructed. It closes the serial port, if open.	N/A	N/A	N/A	N/A	
BuildPortName	Instance	Converts the port parameter into the "COM"+ <number> format</number>	port	String name of a COM port. Can be an integer or string.	N/A	N/A	
IsSerialOpen	Instance	Returns True if the serial instance port is open, false otherwise	N/A	N/A	N/A	N/A	
IsSerialClosed	Instance	Returns False if the serial instance port is open, True otherwise	N/A	N/A	N/A	N/A	
CloseSerialPort	Instance	Closes the instance serial port if it is open.	N/A	N/A	N/A	N/A	
OnenSerialPort	Instance	First, it closes the serial port if it is open. Then, it	port	String of the serial port to be opened.	N/A	Raises an exception if the given port does not exist.	
OpenSerialPort	mstance	opens a serial port with a set baud rate.	baudrate=9600	Integer baud rate of the opened serial port.	IN/A	raises an exception if the given port does not exist.	
SetBaudrate	Instance	If the port is open, it will change the baud rate to the parameter's value	baudrate	Integer baud rate to set for the open serial port.	True if successful at setting the baud rate, false otherwise	N/A	
GetPortName	Instance	Gets the name of the open port.	N/A	N/A	If the serial port is open, it will return a string of the port's name. If the port is closed, it will return None.	N/A	
Read	Instance	Reads a specified number of bytes from the open serial port.	numBytes	Integer number of bytes to read	If the serial port is open, it will return a set number of read bytes. If it is closed, it will return None.	N/A	
ReadLine	Instance	Reads until a new line ('\n') from the open serial port.	N/A	N/A	If the serial port is open, it will return a complete read line. If closed, it will return None.	N/A	
ReadUntil	Instance	Reads until a set character from the open serial port.	eol	end-of-line character	If the serial port is open, it will return a read line ending in eol. If closed, it will return None.	N/A	
Write	Instance	Write a set message to the open serial port.	message	byte string containing the message to write	N/A	N/A	