Python-POD-API_Documentation Welcome

Welcome!				
This document contains information about all the classes in the Python-POD_API.	Each class has its own sheet which describes the class's origin, parent and child classes, imports, global and instance variables, and all methods with details about its parameters, returns, and exceptions.			
Classes				
Setup_8206HR				
POD_8206HR				
POD_Basics				
POD_Packets				
POD_Commands				
COM_io				

4/4/2023

2

Class						
Name	File	Description	Parent	Child	Author	
Setup_8206HR	CodeSetup_ 8206HR.py	Allows a user to set up and stream from any number of 8206HR POD devices. The streamed data is saved to a file.	N/A	N/A	Thresa Kelly	
Imports						
Name	Origin	Description	From			
COM_io	Local	Used to get the COM ports list	SerialCommunication			
POD_8206HR	Local	For managing active 8206HR POD devices	PodDevice_8206HR			
texttable	Enviornment		N/A			
threading	Enviornment	Used to stream from multiple POD devices and ask for user input concurrently.	N/A			
time	Enviornment	For timing the duration of methods	N/A			
numpy	Enviornment	For arrays and data handling	N/A			
path	Enviornment	Used for file handling	os			
EdfWriter		For writing to EDF files	pyedflib			
Variables			F)			
Name	Scope	Description	Value			
_PHYSICAL_BOUND_uV	Class	Physical min/max for EDF file	4069			
podDevices		Dictionary containing all POD 8206HR class objects	Set by user			
_podDevices _podParametersDict	Instance	Dictionary containing all POD_6206HR class objects Dictionary containing the setup parameters for all POD 8206HR objects	Set by user			
saveFileName	Instance	String containing the path, filename, and file extension to a file to save streaming data to. The filename will be extended with " <device number="">" for each device.</device>	Set by user			
_options	Instance	Dictionary listing the different options for the user to complete	{1: 'Start Streaming.', 2: 'Show current settings.', 3: 'Edit save file path.', 4: 'Edit POD device parameters.', 5: 'Connect a new POD device.', 6: 'Reconnect current POD devices.', 7: 'Generate initialization code.', 8: 'Quit.'}			
Methods						
Name	Type	Description	Parameter Name	Parameter Purpose	Return	Exception
Name SetupPODparameters		Description Prompts user for the POD device parameters if none are provided. Used in initialization.	Parameter Name podParametersDict=None	Parameter Purpose Dictionary containing parameters to setup a POD device	Return N/A	Exception N/A
	Instance	Prompts user for the POD device parameters if none are		Dictionary containing parameters to setup a POD		
SetupPODparameters	Instance	Prompts user for the POD device parameters if none are provided. Used in initialization. Asks user for a file path and name to save the POD	podParametersDict=None	Dictionary containing parameters to setup a POD device String describing the file path and file name with	N/A N/A	N/A
SetupPODparameters SetupSaveFile	Instance	Prompts user for the POD device parameters if none are provided. Used in initialization. Asks user for a file path and name to save the POD stream data to if not provided. Used in initialization.	podParametersDict=None saveFile=None	Dictionary containing parameters to setup a POD device String describing the file path and file name with extension	N/A N/A	N/A N/A
SetupPODparameters SetupSaveFile Run	Instance Instance Instance Static	Prompts user for the POD device parameters if none are provided. Used in initialization. Asks user for a file path and name to save the POD stream data to if not provided. Used in initialization. Gives the user options until quitting.	podParametersDict=None saveFile=None N/A	Dictionary containing parameters to setup a POD device String describing the file path and file name with extension N/A	N/A N/A	N/A N/A N/A
SetupPODparameters SetupSaveFile Run _SetNumberOfDevices	Instance Instance Static Instance	Prompts user for the POD device parameters if none are provided. Used in initialization. Asks user for a file path and name to save the POD stream data to if not provided. Used in initialization. Gives the user options until quitting. Asks the user how many POD devices they want.	podParametersDict=None saveFile=None N/A N/A	Dictionary containing parameters to setup a POD device String describing the file path and file name with extension N/A N/A	N/A N/A N/A Integer number for the number of devices	N/A N/A N/A The number must be an integer greater than zero.
SetupPODparameters SetupSaveFile Run _SetNumberOfDevices _DisconnectAllPODdevices	Instance Instance Instance Static Instance Instance	Prompts user for the POD device parameters if none are provided. Used in initialization. Asks user for a file path and name to save the POD stream data to if not provided. Used in initialization. Gives the user options until quitting. Asks the user how many POD devices they want. deletes all POD devices Connects all POD devices Connects all pod devices after disconnecting previous	podParametersDict=None saveFile=None N/A N/A N/A	Dictionary containing parameters to setup a POD device String describing the file path and file name with extension N/A N/A N/A Unique dictionary key representing the POD device dictionary containing parameters to setup a POD	N/A N/A Integer number for the number of devices N/A N/A	N/A N/A N/A The number must be an integer greater than zero. N/A
SetupPODparameters SetupSaveFile Run _SetNumberOfDevices _DisconnectAllPODdevices _ConnectAllPODdevices	Instance Instance Instance Static Instance Instance	Prompts user for the POD device parameters if none are provided. Used in initialization. Asks user for a file path and name to save the POD stream data to if not provided. Used in initialization. Gives the user options until quitting. Asks the user how many POD devices they want. deletes all POD devices Connects all pod devices after disconnecting previous devices Creates a POD_8206HR object and write the POD	podParametersDict=None saveFile=None N/A N/A N/A N/A N/A deviceNum	Dictionary containing parameters to setup a POD device String describing the file path and file name with extension N/A N/A N/A N/A Unique dictionary key representing the POD device	N/A N/A N/A Integer number for the number of devices N/A N/A True when the device is successfully connected, false	N/A N/A N/A The number must be an integer greater than zero. N/A N/A
SetupPODparameters SetupSaveFile Run _SetNumberOfDevices _DisconnectAllPODdevices _ConnectAllPODdevices _ConnectPODdevice	Instance Instance Instance Static Instance Instance Instance Instance	Prompts user for the POD device parameters if none are provided. Used in initialization. Asks user for a file path and name to save the POD stream data to if not provided. Used in initialization. Gives the user options until quitting. Asks the user how many POD devices they want. deletes all POD devices Connects all pod devices after disconnecting previous devices Creates a POD_8206HR object and write the POD parameters to the device. Generates new device number and asks the user for the	podParametersDict=None saveFile=None N/A N/A N/A N/A N/A deviceNum deviceParams	Dictionary containing parameters to setup a POD device String describing the file path and file name with extension N/A N/A N/A Unique dictionary key representing the POD device dictionary containing parameters to setup a POD device	N/A N/A N/A Integer number for the number of devices N/A N/A True when the device is successfully connected, false otherwise.	N/A N/A N/A The number must be an integer greater than zero. N/A N/A N/A
SetupPODparameters SetupSaveFile Run SetNumberOfDevices _DisconnectAllPODdevices _ConnectAllPODdevices _ConnectPODdevice	Instance Instance Instance Static Instance Instance Instance Instance Instance Instance Instance Instance	Prompts user for the POD device parameters if none are provided. Used in initialization. Asks user for a file path and name to save the POD stream data to if not provided. Used in initialization. Gives the user options until quitting. Asks the user how many POD devices they want. deletes all POD devices Connects all pOd devices after disconnecting previous devices Creates a POD_8206HR object and write the POD parameters to the device. Generates new device number and asks the user for the POD device parameters. Asks user for the number of POD devices and requests	podParametersDict=None saveFile=None N/A N/A N/A N/A N/A deviceNum deviceParams N/A	Dictionary containing parameters to setup a POD device String describing the file path and file name with extension N/A N/A N/A N/A Unique dictionary key representing the POD device dictionary containing parameters to setup a POD device N/A	N/A N/A N/A Integer number for the number of devices N/A N/A True when the device is successfully connected, false otherwise. N/A N/A N/A Dictionary containing the POD parameters for one	N/A N/A N/A N/A The number must be an integer greater than zero. N/A N/A N/A N/A
SetupPODparameters SetupSaveFile Run _SetNumberOfDevices _DisconnectAllPODdevices _ConnectAllPODdevices _ConnectPODdevice _AddPODdevice _SetParam_allPODdevices	Instance	Prompts user for the POD device parameters if none are provided. Used in initialization. Asks user for a file path and name to save the POD stream data to if not provided. Used in initialization. Gives the user options until quitting. Asks the user how many POD devices they want. deletes all POD devices Connects all pod devices after disconnecting previous devices Creates a POD_8206HR object and write the POD parameters to the device. Generates new device number and asks the user for the POD device parameters. Asks user for the number of POD devices and requests the parameters for all devices. Creates dictionary from POD device parameters user	podParametersDict=None saveFile=None N/A N/A N/A N/A N/A deviceNum deviceParams N/A N/A	Dictionary containing parameters to setup a POD device String describing the file path and file name with extension N/A N/A N/A Unique dictionary key representing the POD device dictionary containing parameters to setup a POD device N/A N/A	N/A N/A N/A Integer number for the number of devices N/A N/A True when the device is successfully connected, false otherwise. N/A N/A N/A Dictionary containing the POD parameters for one device (port, sample rate, preamplifier gain, and low	NI/A NI/A NI/A The number must be an integer greater than zero. NI/A NI/A NI/A NI/A NI/A NI/A NI/A
SetupPODparameters SetupSaveFile Run _SetNumberOfDevices _DisconnectAllPODdevices _ConnectAllPODdevices _ConnectPODdevice _AddPODdevice _SetParam_allPODdevices _GetParam_onePODdevice	Instance Static	Prompts user for the POD device parameters if none are provided. Used in initialization. Asks user for a file path and name to save the POD stream data to if not provided. Used in initialization. Gives the user options until quitting. Asks the user how many POD devices they want. deletes all POD devices Connects all pod devices after disconnecting previous devices Creates a POD_8206HR object and write the POD parameters to the device. Generates new device number and asks the user for the POD device parameters. Asks user for the number of POD devices and requests the parameters for all devices. Creates dictionary from POD device parameters user input Gets a list of available ports and asks user which port to	podParametersDict=None saveFile=None N/A N/A N/A N/A N/A deviceNum deviceParams N/A N/A N/A forbidden=[]	Dictionary containing parameters to setup a POD device String describing the file path and file name with extension N/A N/A N/A Unique dictionary key representing the POD device dictionary containing parameters to setup a POD device N/A N/A N/A N/A List of String containing the names of ports not	N/A N/A N/A Integer number for the number of devices N/A N/A True when the device is successfully connected, false otherwise. N/A N/A Dictionary containing the POD parameters for one device (port, sample rate, preamplifier gain, and low pass)	N/A N/A N/A The number must be an integer greater than zero. N/A N/A N/A N/A N/A N/A N/A N/A
SetupPODparameters SetupSaveFile Run SetNumberOfDevices DisconnectAllPODdevices _ConnectAllPODdevices _ConnectPODdevice _AddPODdevice _SetParam_allPODdevices _GetParam_onePODdevice _ChoosePort	Instance Static	Prompts user for the POD device parameters if none are provided. Used in initialization. Asks user for a file path and name to save the POD stream data to if not provided. Used in initialization. Gives the user options until quitting. Asks the user how many POD devices they want. deletes all POD devices Connects all pod devices after disconnecting previous devices Creates a POD_8206HR object and write the POD parameters to the device. Generates new device number and asks the user for the POD device parameters. Asks user for the number of POD devices and requests the parameters for all devices. Creates dictionary from POD device parameters user input Gets a list of available ports and asks user which port to connect to the current POD device. Gets a list of available port names. If there are no ports, the program will wait until another POD device is plugged	podParametersDict=None saveFile=None N/A N/A N/A N/A N/A deviceNum deviceParams N/A N/A N/A forbidden=[]	Dictionary containing parameters to setup a POD device String describing the file path and file name with extension N/A N/A N/A N/A Unique dictionary key representing the POD device dictionary containing parameters to setup a POD device N/A N/A N/A List of String containing the names of ports not allowed to be connected to List of String containing the names of ports not	N/A N/A N/A Integer number for the number of devices N/A N/A True when the device is successfully connected, false otherwise. N/A N/A Dictionary containing the POD parameters for one device (port, sample rate, preamplifier gain, and low pass) String name of the chosen port	N/A N/A N/A The number must be an integer greater than zero. N/A N/A N/A N/A N/A N/A N/A The user input must match an allowed COM port
SetupPODparameters SetupSaveFile RunSetNumberOfDevicesDisconnectAllPODdevicesConnectAllPODdevicesConnectPODdeviceAddPODdeviceSetParamallPODdevicesGetParam_onePODdeviceChoosePortGetPortsList	Instance Static Static Static	Prompts user for the POD device parameters if none are provided. Used in initialization. Asks user for a file path and name to save the POD stream data to if not provided. Used in initialization. Gives the user options until quitting. Asks the user how many POD devices they want. deletes all POD devices Connects all pod devices after disconnecting previous devices Creates a POD_8206HR object and write the POD parameters to the device. Generates new device number and asks the user for the POD device parameters. Asks user for the number of POD devices and requests the parameters for all devices. Creates dictionary from POD device parameters user input Gets a list of available ports and asks user which port to connect to the current POD device. Gets a list of available port names. If there are no ports, the program will wait until another POD device is plugged in. Asks user for the sample rate.	podParametersDict=None saveFile=None N/A N/A N/A N/A N/A deviceNum deviceParams N/A N/A N/A forbidden=[]	Dictionary containing parameters to setup a POD device String describing the file path and file name with extension N/A N/A N/A N/A Unique dictionary key representing the POD device dictionary containing parameters to setup a POD device N/A N/A N/A N/A List of String containing the names of ports not allowed to be connected to	N/A N/A N/A Integer number for the number of devices N/A N/A True when the device is successfully connected, false otherwise. N/A N/A N/A Dictionary containing the POD parameters for one device (port, sample rate, preamplifier gain, and low pass) String name of the chosen port List of strings containing the name of all COM ports. Integer number between 100-2000 Hz for the sample rate	NI/A NI/A NI/A The number must be an integer greater than zero. NI/A NI/A NI/A NI/A NI/A NI/A The user input must match an allowed COM port There must be COM ports in use Sample rate must be an integer between 100-2000
SetupPODparameters SetupSaveFile Run _SetNumberOfDevices _DisconnectAllPODdevices _ConnectAllPODdevices _ConnectPODdevice _AddPODdevice _SetParam_allPODdevices _GetParam_onePODdevice _ChoosePort _GetPortsList _ChooseSampleRate	Instance Static Static Static	Prompts user for the POD device parameters if none are provided. Used in initialization. Asks user for a file path and name to save the POD stream data to if not provided. Used in initialization. Gives the user options until quitting. Asks the user how many POD devices they want. deletes all POD devices Connects all pod devices after disconnecting previous devices Creates a POD _8206HR object and write the POD parameters to the device. Generates new device number and asks the user for the POD device parameters. Asks user for the number of POD devices and requests the parameters dictionary from POD device parameters user input Gets a list of available ports and asks user which port to connect to the current POD device. Gets a list of available port names. If there are no ports, the program will wait until another POD device is plugged in.	podParametersDict=None saveFile=None N/A N/A N/A N/A N/A deviceNum deviceParams N/A N/A N/A forbidden=[] forbidden=[]	Dictionary containing parameters to setup a POD device String describing the file path and file name with extension N/A N/A N/A N/A Unique dictionary key representing the POD device dictionary containing parameters to setup a POD device N/A N/A N/A List of String containing the names of ports not allowed to be connected to List of String containing the names of ports not allowed to be connected to N/A	N/A N/A N/A Integer number for the number of devices N/A N/A True when the device is successfully connected, false otherwise. N/A N/A Dictionary containing the POD parameters for one device (port, sample rate, preamplifier gain, and low pass) String name of the chosen port List of strings containing the name of all COM ports. Integer number between 100-2000 Hz for the sample	NI/A NI/A NI/A The number must be an integer greater than zero. NI/A NI/A NI/A NI/A NI/A NI/A NI/A The user input must match an allowed COM port There must be COM ports in use
SetupPODparameters SetupSaveFile Run _SetNumberOfDevices _DisconnectAllPODdevices _ConnectAllPODdevices _ConnectPODdevice _AddPODdevice _SetParam_allPODdevices _GetParam_onePODdevice _ChoosePort _GetPortsList _ChooseSampleRate _ChoosePreampGain	Instance Instance Instance Static Instance Instance Instance Instance Instance Instance Instance Instance Static Static Static Static Static Static Static Static	Prompts user for the POD device parameters if none are provided. Used in initialization. Asks user for a file path and name to save the POD stream data to if not provided. Used in initialization. Gives the user options until quitting. Asks the user how many POD devices they want. deletes all POD devices Connects all pod devices after disconnecting previous devices Connects all pod devices after disconnecting previous devices Creates a POD_8206HR object and write the POD parameters to the device. Generates new device number and asks the user for the POD device parameters. Asks user for the number of POD devices and requests the parameters for all devices. Creates dictionary from POD device parameters user input Gets a list of available ports and asks user which port to connect to the current POD device. Gets a list of available port names. If there are no ports, the program will wait until another POD device is plugged in. Asks user for the sample rate. Asks user for the preamplifier gain of their POD device	podParametersDict=None saveFile=None N/A N/A N/A N/A N/A N/A deviceNum deviceParams N/A N/A N/A forbidden=[] forbidden=[] N/A N/A	Dictionary containing parameters to setup a POD device String describing the file path and file name with extension N/A N/A N/A N/A Unique dictionary key representing the POD device dictionary containing parameters to setup a POD device N/A N/A N/A List of String containing the names of ports not allowed to be connected to List of String containing the names of ports not allowed to be connected to N/A N/A	N/A N/A N/A Integer number for the number of devices N/A N/A True when the device is successfully connected, false otherwise. N/A N/A Dictionary containing the POD parameters for one device (port, sample rate, preamplifier gain, and low pass) String name of the chosen port List of strings containing the name of all COM ports. Integer number between 100-2000 Hz for the sample rate Integer 10 or 100 for the preamplifier gain Dictionary containing lowpass filters (EEG1, EEG2,	NI/A NI/A NI/A The number must be an integer greater than zero. NI/A NI/A NI/A NI/A NI/A NI/A NI/A The user input must match an allowed COM port There must be COM ports in use Sample rate must be an integer between 100-2000 Gain must be an integer value of 10 or 100

4/4/2023

_EditParams	Instance	Asks the user which device to edit, then asks for the new parameters	N/A	N/A	N/A	N/A
_SelectPODdeviceFromDictToEdit	Instance	Asks user for the pod device number to edit	N/A	N/A	Integer number for the pod device key	User must input an integer number that exists in the POD parameters dict keys
_GetForbiddenNames	Instance	Gets a list of all port names used by the POD devices. One string name can be excluded.	exclude=None	String of the name to exclude in the returned list of names	List of strings containing the port names for all POD devices	N/A
_PrintDeviceNumber	Static	Prints a string wrapper containing the device number	num	Integer representing the device number	N/A	N/A
_PrintPODdeviceParamDict	Instance	Prints string of the POD parameters dictionary. This can be used as a parameter to SetupPODparameters()	N/A	N/A	N/A	N/A
_DisplayPODdeviceParameters	Instance	Prints a table containing the parameters for all POD devices	N/A	N/A	N/A	N/A
BuildFileName	Static	Adds the device number to the file name with path	fileName devNum	String of the file path and name with extension Integer number of the device	String of the filename with added device number: path\filename <device#>.ext</device#>	N/A
_PrintSaveFile	Instance	Prints the path/file name of the file that streaming data is saved to	N/A	N/A	N/A	N/A
_CheckFileExt	Static	Checks for valid file extension	f flsExt=True goodExt=['.csv','.txt','.edf'] printErr=True	file name or extension Boolean flag that is true if f is an extension, false otherwise List of valid file extensions Boolean flag that, when true, will print an error statement	True if extension is in goodExt list, False otherwise	N/A
_GetFilePath	Static	Asks user for a path and filename to save streaming data to.	N/A	N/A	String of the file path, name, and extension.	Filename must end in .csv, .txt, or .edf
_GetFileName	Static	Asks the user for a filename	N/A	N/A	String of the file name and extension	Filename must end in .csv, .txt, or .edf
	Instance	Opens a file using the user's "path/filename_ <device#>.extension"</device#>	devNum	Integer number for the DEVICE#	Opened write file	N/A
_OpenSaveFile_TXT	Static	Opens a text file and write the column names	fname	String filename	Opened file	N/A
_OpenSaveFile_EDF	Instance	Opens EDF file and write header	fname devNum	String filename Integer device number	Opened file	N/A
_WriteDataToFile_TXT	Static	Writes data to an open text file	file data sampleRate t	opened write file List of 3 items, one for each channel Integer sample rate in Hz integer time (in seconds) corresponding to the data	N/A	N/A
_WriteDataToFile_EDF	Static	Writes data to an open EDF file	file data	opened EDF file List of 3 items, one for each channel	N/A	N/A
_AskToStopStream	Instance	Writes command to stop streaming to all POD devices after user clicks enter	N/A	N/A	N/A	N/A
_StreamUntilStop	Static	Continually reads and translates streaming data from POD device until end stream command is read. Streaming data is written to a file every second.	pod file sampleRate	POD_8206HR device to read data from Open file to write data to Sample rate for the pod device	N/A	N/A
_StreamThreading	Instance	Creates a thread for each POD device to stream and a thread for user input to stop streaming.	N/A	N/A	N/A	N/A
_Stream	Instance	Tests that all devices are connected then starts streaming. It also prints the streaming execution time	N/A	N/A	N/A	N/A
_PrintOptions	Instance	Prints options available for user	N/A	N/A	N/A	N/A
_AskOption	Instance	Asks user which option to do	N/A	N/A	Integer number representing an option key	User input must be an integer that is a key in the options dictionary.
_DoOption	Instance	Performs the methods associated with the user selected option	choice	Integer number representing an option key	N/A	N/A
_AskYN	Static	Asks a user a yes or no question.	question	String of the question to ask the user	True when user inputs yes, False when user inputs no	User must input y/yes or n/no (case insentive)
_TimeFunc	Static	Runs a function and gets the calculated execultion time	func	function/method name	Float of the execution time in seconds rounded to 3 decimal places	N/A
_TestDeviceConnection	Static	Writes a PING command to the POD device and reads the response.	pod	POD_8206HR device to read data from	True if PING command is read from the device	N/A
_TestDeviceConnection_All	Instance	Checks the connection of all current POD devices	N/A	N/A	True if all devices are connected, false otherwise	N/A
_PrintlnitCode	Instance	Prints code that can be used to initialize and run Setup_8206HR with the current parameters.	N/A	N/A	N/A	N/A
_uV	Static	Converts volts to microVolts, rounded to 6 decimal places	voltage	number of volts	number of uV	N/A

Class						
Name	File	Description	Parent	Child	Author	
POD_8206HR	PodDevice_ 8206HR.py	Handles communication using an 8206HR POD	POD_Basics	N/A	Thresa Kelly	
Imports						
Name	Origin	Description	From			
POD_Basics	Local	For inheritance	BasicPodProtocol			
POD_Packets	Local	For handling POD packets	PodPacketHandling			
POD_Commands	Local	For command constants	PodCommands			
Variables						
Name	Scope	Description	Value			
B4LENGTH	Class	Constant containing the number of bytes for a full Binary4 packet	16			
_preampGain	Instance	Preamplifier gain	10 or 100			
Methods						
Name	Туре	Description	Parameter Name	Parameter Purpose	Return	Exception
init	Dunder	Runs when an instance is constructed. It runs the parent's initialization. Then it updates the _commands to contain the appropriate commands for an 8306HR	port baudrate=9600	String of the serial port to be opened. Used when initializing the COM_io instance. Integer baud rate of the opened serial port. Used	N/A	N/A
UnpackPODpacket_Binary	Static	POD device. Overwrites the parent's method. Separates the components of a binary4 packet into a dictionary.	msg	when initializing the COM_io instance. Bytes string containing a complete binary4 Pod packet: STX (1 byte) + command (4 bytes) + packet number (1 bytes) + TTL (1 byte) + ch0 (2 bytes) + ch1 (2 bytes) + ch2 (2 bytes) + checksum (2 bytes) + ETX (1 byte)	A dictionary containing 'Command Number', 'Packet #', 'TTL', 'Ch0', 'Ch1', and 'Ch2' in bytes.	An exception is raised if (1) the packet does not have the minimum number of bytes, (2) does not begin with STX, or (3) does not end with ETX.
TranslatePODpacket_Binary	Static	Overwrites the parent's method. Unpacks the binary4 POD packet and converts the values of the ASCII-encoded bytes into integer values and the values of binary-encoded bytes into integers. Channel values are given in Volts.	msg	Bytes string containing a complete binary4 Pod packet: STX (1 byte) + command (4 bytes) + packet number (1 bytes) + TTL (1 byte) + ch0 (2 bytes) + ch1 (2 bytes) + ch2 (2 bytes) + checksum (2 bytes) + ETX (1 byte)	A dictionary containing 'Command Number', 'Packet #, 'TTL', 'Ch0', 'Ch1', and 'Ch2' as numbers.	N/A
UnpackPODpacket	Static	Determines if the packet is standard or binary, and unpacks accordingly.	msg	Bytes string containing either a standard or binary4 packet	A dictionary containing the unpacked message in bytes	N/A
TranslatePODpacket	Static	Determines if the packet is standard or binary, and translates accordingly.	msg	Bytes string containing either a standard or binary4 packet	A dictionary containing the unpacked message in integers	N/A
BinaryBytesToVoltage	Instance	Converts a binary bytes value read from POD device and converts it to the real voltage value at the preamplifier input	value	Bytes string containing voltage measurement	A number containing the voltage in Volts [V].	N/A
Road Binary	Instance	After receiving the prePacket, it reads the 8	prePacket	Bytes string containing the beginning of a POD packet: STX (1 byte) + command number (4 bytes)	Byte string for a binary4 POD packet.	N/A
_Read_Binary Instance	instance		validateChecksum=True	Set to True to validate the checksum. Set to False to skip validation	byte string for a binary4 POD packet.	IVA

Class						
Name	File	Description	Parent	Child	Author	
POD_Packets	PodPacket	Collection of methods for creating and interpreting POD packets	N/A	N/A	Thresa Kelly	
Imports						
Name	Origin	Description	From			
N/A	N/A	N/A	N/A			
Methods						
Name	Type	Description	Parameter Name	Parameter Purpose	Return	Exception
STX	Static	Get STX in bytes. STX marks the starting byte of a POD Packet	N/A	N/A	Bytes for STX (0x02)	N/A
ETX	Static	Get ETXin bytes. ETX marks the end byte of a POD Packet	N/A	N/A	Bytes for ETX(0x03)	N/A
IntToAsciiBytes	Static	uppercase hexadecimal string. Then it converts the	value	Integer value to be converted into ASCII-encoded bytes	Bytes that are ASCII-encoded conversions of the value parameter.	N/A
Introduction	Static	ensures that the final message is the desired length. Example: if value=2 and numBytes=4, the returned ASCII will show b'0002', which is '0x30 0x30 0x30 0x32' in bytes.	numBytes	Number bytes to be the length of the ASCII-encoded message.		
AsciiBytesToInt	Static	Converts a ASCII-encoded bytes message into an integer. It does this using a base-16 conversion.	msg_b	Bytes message to be converted to an integer. The bytes must be base-16 or the conversion will fail.	Integer result from the ASCII-encoded byte conversion.	N/A
			msg	Bytes message holding binary information to be converted into an integer.		
BinaryBytesToInt	Static	Converts binary-encoded bytes into an integer	byteorder='big'	Ordering of bytes. 'big' for big endian and 'little' for little endian.	Integer result from the binary-encoded bytes message.	N/A
			signed=False	Boolean flag to mark if the msg is signed (True) or unsigned (False)		
Checksum	Static	Calculates the checksum of a given bytes message. This is achieved by summing each byte in the message, inverting, and taking the last byte.	bytesIn	Bytes message containing POD packet data	Two ASCII-encoded bytes containing the checksum for bytesIn	N/A
BuildPODpacket_Standard	Static	Builds a standard POD packet STX (1 byte) + command number (4 bytes) + optional packet (? bytes) + checksum (2 bytes) + ETX (1 bytes) as	commandNumber	Integer representing the command number. This will be converted into a 4 byte long ASCII-encoded bytes string.	Bytes string of a complete standard POD packet	N/A
		bytes.	payload=None	bytes string containing the payload		
PayloadToBytes	Static	Converts a payload into a bytes string	payload	Integer, bytes, or tuple containing the payload	Bytes string of the payload	Raises an Exception when the payload argument is
r ayload lobytes	Static	Converts a payload lifto a bytes stillig	argSizes	Tuple of the argument sizes	Dytes string of the payload	an incorrect type or formatted incorrectly.

Class						
Name	File	Description	Parent	Child	Author	
rume		Handle basic communication with a POD device,	I WIGHT	- Ciliu	Autivi	
POD_Basics	BasicPodPr otocol.py	including reading and writing packets and packet interpretation.	N/A	POD_8206HR	Thresa Kelly	
Imports						
Name	Origin	Description	From			
COM_io	Local	For opening and connecting serial COM ports	SerialCommunication			
POD_Packets	Local	For handling POD packets	PodPacketHandling			
POD_Commands	Local	Used to contain all POD commands in the class instance	PodCommands			
Variables						
Name	Scope	Description	Value			
NUMPOD	Class	Integer equal to the number of POD_Basics class instances. Incremented on construction and decremented on destruction	0			
_port	Instance	Open serial port via COM_io class instance	COM_io			
_commands	Instance	Command handler POD_Commands class instance	POD_Commands			
Methods						
Name	Туре	Description	Parameter Name	Parameter Purpose	Return	Exception
		Runs when an instance of POD Basics is		String of the serial port to be opened. Used when		
init	Dunder	constructed. It initializes the instance variable for the COM port communication (_port) and for the	port	initializing the COM_io instance. Integer baud rate of the opened serial port. Used	N/A	N/A
		command handler (_commands). It also increments the POD device counter (NUMPOD).	baudrate=9600	when initializing the COM_io instance.		
del	Dunder	Runs when an instance is destructed. It decrements the POD device counter (NUMPOD)	N/A	N/A	N/A	N/A
GetNumberOfPODDevices	Static	Get the POD device counter	N/A	N/A	Integer of the number of class instances (NUMPOD).	N/A
UnpackPODpacket_Standard	Static	Converts a standard POD packet into a dictionary containing the command number and payload (if applicable) in bytes.	msg	Bytes message containing a standard POD packet: STX (1 byte) + command number (4 bytes) + optional packet (? bytes) + checksum (2 bytes) + ETX (1 bytes)	A dictionary containing the POD packet's 'Command Number' and 'Payload' (if applicable) in bytes.	An exception is raised if (1) the msg does not have the minimum number of bytes in a standard pod packet, (2) does not begin with STX, and (3) does not end with ETX.
UnpackPODpacket_Binary	Static	Converts a variable-length binary packet into a dictionary containing the command number, binary packet length, and binary data in bytes.	msg	Bytes message containing a variable-length POD packet: STX (1 byte) + command number (4 bytes) + length of binary (4 bytes) + checksum (2 bytes) + ETX (1 bytes) + binary (LENGTH bytes) + checksum (2 bytes) + ETX (1 bytes)	A dictionary containing the 'Command Number', 'Binary Packet Length', and 'Binary Data' in bytes.	An exception is raised if (1) the msg does not have the minimum number of bytes in a standard pod packet, (2) does not begin with STX, (3) does not end with ETX, and (4) does not have an ETX after standard packet.
TranslatePODpacket_Standard	Instance	Unpacks the standard POD packet and converts the ASCII-encoded bytes values into integer values.	msg	Bytes message containing a standard POD packet	A dictionary containing the POD packet's 'Command Number' and 'Payload' (if applicable) in integers.	N/A
TranslatePODpacket_Binary	Static	Unpacks the variable-length binary POD packet and converts the values of the ASCII-encoded bytes into integer values and leaves the binary-encoded bytes as is.	msg	Bytes message containing a variable-length POD packet	A dictionary containing the 'Command Number' and 'Binary Packet Length' in integers, and 'Binary Data' in bytes.	N/A
_ValidateChecksum	Static	Validates the checksum of a given POD packet. The checksum is valid if the calculated checksum from the data matches the checksum written in the packet.	msg	Bytes message containing a POD packet: STX (1 bytes) + data (? bytes) + checksum (2 bytes) + ETX (1 byte).	Returns True if the checksum is correct, false otherwise.	An exception is raised if the msg does not begin with STX or end with ETX.
GetDeviceCommands	Instance	Gets the dictionary containing the class instance's available POD commands.	N/A	N/A	Dictionary containing the available commands and their information. Formatted as key(command number): value([command name, number of argument ASCII bytes, number of return bytes, binary flag])	N/A
SetBaudrateOfDevice	Instance	If the port is open, it will change the baud rate to the parameter's value	baudrate	Integer baud rate to set for the open serial port.	True if successful at setting the baud rate, false otherwise	N/A
			cmd	An integer representing the command number.	Distance at single containing a DOD and the stand	
WriteRead	Instance	Writes a command with optional payload to POD device, then reads (once) the device response.	payload=None	None when there is no payload. If there is a payload, set to an integer value or a bytes string.	Bytes string containing a POD packet beginning with STX and ending with ETX. This may be a standard packet, binary packet, or an unformatted packet	N/A
		, , , , , , , , , , , , , , , , , , , ,	validateChecksum=True	Set to True to validate the checksum. Set to False to skip validation	(STX+something+ETX).	
		Builds a POD packet and writes it to a POD device via	cmd	An integer representing the command number.		An exception is raised if (1) the command does not
GetPODpacket	Instance	COM port. If an integer payload is give, the method will convert it into a bytes string of the length expected by the command. If a bytes payload is given, it must be the correct length.	payload=None	None when there is no payload. If there is a payload, set to an integer value, bytes string, or tuple	Returns the bytes string of the POD packet.	exist for the instance, (2) a payload is not given when the command expects one. (3) the payload (given in bytes) is the size not expected by the command, or (4) the payload is given as a type other than integer or bytes.
WritePacket	Instance	Ruilds a POD nacket and writes it to the POD device	cmd	An integer representing the command number.	Returns the bytes string that was written to the POD	N/Δ

Python-POD-API_Documentation POD_Basics

www.	Ingranac	Number of the first property and wireless in the tree				. PVO
vvriteracket	instance	Bullos a POD packet and writes it to the POD device.	payload=None	None when there is no payload. If there is a payload, set to an integer value, bytes string, or tuple	device	N/A
ReadPODpacket	Instance	Reads a complete POD packet, either in standard or binary format, beginning with STX and ending with ETX. Reads first STX and then starts recursion.	validateChecksum=True	Set to True to validate the checksum. Set to False to skip validation	Bytes string containing a POD packet beginning with STX and ending with ETX. This may be a standard packet, binary packet, or an unformatted packet (STX+something+ETX).	N/A
_ReadPODpacket_Recursive	Instance	Reads the command number. If the command number ends in ETX, the packet is returned. Next, it checks if the command is allowed. Then, it checks if the command is standard or binary and reads accordingly, then returns the packet.	validateChecksum=True	Set to True to validate the checksum. Set to False to skip validation	Bytes string containing a POD packet beginning with STX and ending with ETX. This may be a standard packet, binary packet, or an unformatted packet (STX+something+ETX).	N/A
_Read_GetCommand	Instance	Reads one byte at a time up to 4 bytes to get the ASCII-encoded bytes command number. For each byte read, it can (1) start the recursion over if an STX is found, (2) returns if ETX is found, or (3) continue building the command number.	validateChecksum=True	Set to True to validate the checksum. Set to False to skip validation	4 byte long string containing the ASCII-encoded command number.	An exception is raised if the command number is not allowed for the POD device
_Read_ToETX	Instance	Reads one byte at a time until an ETX is found. It will restart the recursive read if an STX is found anywhere.	validateChecksum=True	Set to True to validate the checksum. Set to False to skip validation	Bytes string ending with ETX	N/A
Read Standard	Instance	Reads the payload, checksum, and ETX. Then it	prePacket	Bytes string containing the beginning of a POD packet: STX (1 byte) + command number (4 bytes)	Bytes string for a complete standard POD packet	An exception is raised if the checksum is invalid (only
_Read_Standard	instance	builds the complete standard POD packet in bytes.	validateChecksum=True	Set to True to validate the checksum. Set to False to skip validation	bytes string for a complete standard POD packet	if validateChecksum=True)
	Instance	Reads the remaining part of the variable-length binary packet. It first reads the standard packet (prePacket+payload+checksum+ETX). Then it determines how long the binary packet is from the payload of the standard POD packet and reads that many bytes. It then reads to ETX to get the checksum+ETX.	prePacket	Bytes string containing the beginning of a POD packet: STX (1 byte) + command number (4 bytes)	Bytes string for a variable-length binary POD packet	An exception is raised if the checksum is invalid (only if validateChecksum=True)
_Read_Binary			validateChecksum=True	Set to True to validate the checksum. Set to False to skip validation		

Class						
Name	File	Description	Parent	Child	Author	
		Manages a dictionary containing available commands		-		
POD_Commands	nds.py	for a POD device.	N/A	N/A	Thresa Kelly	
Imports						
Name	Origin	Description	From			
N/A	N/A	N/A	N/A			
Variables						
Name	Scope	Description	Value			
NAME	Class	index key for the command name forcommands list values	0			
ARGUMENTS	Class	index key for the number of bytes in an argument forcommands list values	1			
RETURNS	Class	index key for the number of bytes in the return forcommands list values	2			
BINARY	Class	index key for the binary flag forcommands list values	3			
NOVALUE	Class	Integer used to mark when a list item incommands means 'no value' or undefined.	-1			
U8	Class	Number of bytes for an unsigned 8-bit value	2			
U16	Class	Number of bytes for an unsigned 16-bit value	4			
_commands	Instance	Dictionary containing the available commands for a POD device. Each entry is formatted as { key(command number) : value([command name, number of argument ASCII bytes, number of return bytes, binary flag) }	{ key : value }			
Methods						
Name	Туре	Description	Parameter Name	Parameter Purpose	Return	Exception
init	Dunder	Runs whan an instance is constructed. It sents the commands dictionary to the basic command set.	N/A	N/A	N/A	N/A
NoValue	Static	Gets value ofNOVALUE	N/A	N/A	Value ofNOVALUE	N/A
U8	Static	Gets value ofU8	N/A	N/A	Value ofU8	N/A
U16	Static	Gets value ofU16	N/A	N/A	Value ofU16	N/A
GetBasicCommands	Static	Creates a dictionary containing the basic POD command set (0,1,2,3,4,5,6,7,8,9,10,11,12)	N/A	N/A	N/A	N/A
GetCommands	Instance	Gets the contents of the current command dictionary (_commands)	N/A	N/A	N/A	N/A
RestoreBasicCommands	Instance	Sets the current commands (commands) to the basic POD command set.	N/A	N/A	N/A	N/A
			commandNumber	Integer of the command number		
		Adds a command entry to the current commands	commandName	String of the command's name	True if the command was successfully added, False if	
AddCommand	Instance	dictionary (commands) if the command does not	argumentBytes	Integer of the number of bytes in the argument	the command could not be added because t already	N/A
		exist	returnBytes	Integer of the number of bytes in the return	exists.	
			isBinary	Boolean flag to mark if the command is binary (True) or standard (False)		
RemoveCommand	Instance	Removes the entry for a given command incommands dictionary.	cmd	integer command number or string command name.	True if the command was successfully removed, False if the command does not exist.	N/A
CommandNumberFromName	Instance	Gets the command number from the command dictionary using the command's name	name	string of the command's name	Integer representing the command number. If the command could not be found, return None.	N/A
ArgumentBytes	Instance	Gets the tuple for the number of bytes in the argument for a given command.	cmd	integer command number or string command name.	Tuple representing the number of bytes in the argument for cmd. If the command could not be found, return None.	N/A
ReturnBytes	Instance	Gets the tuple for the number of bytes in the return for a given command.	cmd	integer command number or string command name.	Tuple representing the number of bytes in the return for cmd. If the command could not be found, return None.	N/A
IsCommandBinary	Instance	Gets the binary flag for a given command	cmd	integer command number or string command name.	Boolean flag that is True if the command is binary and False if standard. If the command could not be found, return None.	N/A
		Checks if a command exists in thecommands				

Class						
Name	File	Description	Parent	Child	Author	
Name			Parent	Child	Author	
COM_io		Handle serial communication (read/write) using COM ports.	N/A	N/A	Thresa Kelly	
Imports						
Name	Origin	Description	From			
serial.tools.list_ports	Enviornment	For accessing the COM ports on the computer	N/A			
Variables						
Name	Scope	Description	Value			
serialInst	Instance	Serial object to set the port and baud rate to. It can be opened or closed.	Serial			
Methods						
Name	Туре	Description	Parameter Name	Parameter Purpose	Return	Exception
GetCOMportsList	Static	Finds all the available COM ports on the user's computer and appends them to an accessible list.	N/A	N/A	List containing the names of available COM ports	N/A
1-14	D	Runs when the object is constructed. It initialized the	port	String of the serial port to be opened.	N/A	N/A
init	Dunder	serialInst to a given COM port with a set baudrate.	baudrate=9600	Integer baud rate of the opened serial port.		
del	Dunder	Runs when the object is destructed. It closes the serial port, if open.	N/A	N/A	N/A	N/A
BuildPortName	Instance	Converts the port parameter into the "COM"+ <number> format</number>	port	String name of a COM port. Can be an integer or string.	N/A	N/A
IsSerialOpen	Instance	Returns True if the serial instance port is open, false otherwise	N/A	N/A	N/A	N/A
IsSerialClosed	Instance	Returns False if the serial instance port is open, True otherwise	N/A	N/A	N/A	N/A
CloseSerialPort	Instance	Closes the instance serial port if it is open.	N/A	N/A	N/A	N/A
OnenCarialDart	Instance	First, it closes the serial port if it is open. Then, it	port	String of the serial port to be opened.	N/A	Raises an exception if the given port does not exist.
OpenSerialPort	Instance	opens a serial port with a set baud rate.	baudrate=9600	Integer baud rate of the opened serial port.	N/A	Raises an exception if the given port does not exist.
SetBaudrate	Instance	If the port is open, it will change the baud rate to the parameter's value	baudrate	Integer baud rate to set for the open serial port.	True if successful at setting the baud rate, false otherwise	N/A
GetPortName	Instance	Gets the name of the open port.	N/A	N/A	If the serial port is open, it will return a string of the port's name. If the port is closed, it will return None.	N/A
Read	Instance	Reads a specified number of bytes from the open serial port.	numBytes	Integer number of bytes to read	If the serial port is open, it will return a set number of read bytes. If it is closed, it will return None.	N/A
ReadLine	Instance	Reads until a new line ('\n') from the open serial port.	N/A	N/A	If the serial port is open, it will return a complete read line. If closed, it will return None.	N/A
ReadUntil	Instance	Reads until a set character from the open serial port.	eol	end-of-line character	If the serial port is open, it will return a read line ending in eol. If closed, it will return None.	N/A
Write	Instance	Write a set message to the open serial port.	message	byte string containing the message to write	N/A	N/A

COM_io

4/4/2023