



Class Name	File	Description	Parent	Child	Author	
Setup_PodDevices	Setup_PodDevice.py	Setup_PodDevice allows a user to set up and stream data from any number of POD devices. The streamed data is saved to a file.	N/A		Thresa Kelly	
Imports Name	Origin	Description	From			
<code>time</code>	Environment	For timing the duration of methods				
<code>os</code>	Environment	Used for file handling				
<code>Thread</code>	Environment	Used to stream from multiple POD devices and ask for user input concurrently.	threading			
<code>floor</code>	Environment	For rounding numbers	math			
<code>Setup_8206HR</code>	Local	For managing active 8206-HR POD devices	Setup_8206HR			
Variables Name	Scope	Description	Value	Type		
<code>_setupPodDevices</code>	Instance	Dictionary containing the Setup_Interface subbases for each POD device.	{ '8206-HR' : Setup_8206HR }	dict<str,Setup_Interface>		
<code>_saveFileName</code>	Instance	String containing the path, filename, and file extension to a file to save streaming data to. The filename will be extended with "<DEVICE NAME>_<NUMBER>" for each device.	Set by user	str		
<code>_options</code>	Instance	Dictionary listing the different options for the user to complete	[{ '1' : Start Streaming, '2' : Show current settings, '3' : Edit save file path, '4' : Edit POD device parameters, '5' : Connect a new POD device, '6' : Reconnect current POD devices, '7' : Generate initialization code, '8' : 'Quit' }]	dict[int,str]		
Methods Name	Type	Description	Parameter Name	Parameter Purpose	Return	Exception
<code>__init__</code>	Dunder	Initializes the class. Sets the default values of the class instance variables. Calls functions to complete the class setup.	saveFile:str None=None	String describing the directory path and filename with an extension	N/A	N/A
<code>del__</code>	Dunder	Deletes all POD device setup objects.	podParametersDict:dict[str,N/A None]={'8206-HR':None}	Dictionary of POD devices and their respective initialization dictionaries.	N/A	N/A
<code>GetPODParametersDict</code>	Instance	Gets the POD device initialization dictionaries for all device types	N/A	N/A	Dictionary whose keys are the POD device name, and N/A	
<code>GetSaveFileName</code>	Instance	Gets the name of the class object's save file	N/A	N/A	Value the setup dictionary.	
<code>GetOptions</code>	Instance	Gets the dictionary of setup options	N/A	N/A	String of the save file name and path (saveFile:filename)	
<code>SetupPODParameters</code>	Instance	Sets up each POD device type. Used in initialization.	podParametersDict:dict[str,N/A None]={'8206-HR':None}	Dictionary of all POD devices initialization. The keys are the device name and the entities are the initialization dictionaries.	N/A	N/A
<code>SetupSavefile</code>	Instance	Gets the path/file name from the user and stores it. Used in initialization.	saveFile:str None=None	String of the save file which includes the directory path, filename, and file extension	N/A	N/A
<code>Run</code>	Instance	Prints the options and asks the user what to do. Loops until 'Quit' is chosen.	N/A	N/A	N/A	N/A
<code>PrintOptions</code>	Instance	Prints options available for user	N/A	N/A	N/A	N/A
<code>_AskOption</code>	Instance	Asks user which option to do	choice:int	Integer number representing an option key	N/A	User input must be an integer that is a key in the options dictionary.
<code>_DoOption</code>	Instance	Performs the methods associated with the user selected option	N/A	N/A	Float of the execution time in seconds	N/A
<code>_Stream</code>	Instance	Streams data from all POD devices and prints the execution time.	N/A	N/A	N/A	N/A
<code>_ShowCurrentSettings</code>	Instance	Displays the POD device settings for all devices, and then prints the save file name	N/A	N/A	N/A	N/A
<code>_EditSaveFilePath</code>	Instance	Asks the user for a new file name and path, then sets the value to the POD devices.	N/A	N/A	N/A	N/A
<code>_EditCheckConnect</code>	Instance	Displays the POD devices parameters, asks the user to edit the device, and then reconnects the device for each POD device type.	N/A	N/A	N/A	N/A
<code>_ConnectNewDevice</code>	Instance	Asks the user for the POD device type, then it sets up that device	N/A	N/A	N/A	N/A
<code>_Reconnect</code>	Instance	Reconnects all POD devices	N/A	N/A	Bool that is true if all devices were successfully connected. False otherwise	N/A
<code>_PrintInitCode</code>	Instance	Prints code that can be used to initialize and run SetupPodDevices with the current parameters.	N/A	N/A	N/A	N/A

_PrintSaveFile	Instance	Prints the file path and name that data is saved to. Note that the device name and number will be appended to the end of the filename.	N/A	N/A	N/A	N/A	N/A
_CheckFfileExt	Static	Checks for valid file extension	f	file name or extension Booleam flag that is true if f is an extension, false otherwise	True if extension is in goodExt list, False otherwise	N/A	N/A
_GetFilePath	Static	Asks user for a path and filename to save streaming data to.	N/A	String of the file path, name, and extension.	Filname must end in .csv, .txt, or .edf		
_GetFileName	Static	Asks the user for a filename	N/A	String of the file name and extension	Filname must end in .csv, .txt, or .edf		
_SetfilenameToDevices	Instance	Sets the filename to each POD device type	N/A	N/A	N/A	N/A	N/A
_StreamAllDevices	Instance	Streams data from all the devices. User is asked to click enter to stop streaming. Data is saved to file. Uses threading.	N/A	N/A	N/A	N/A	N/A
_AskToStopStream	Instance	Asks user to press enter to stop streaming. The program will then prompt all POD devices to end stream.	N/A	N/A	N/A	N/A	N/A
_TimeFunc	Static	Runs a function and gets the calculated execution time	'func: 'function'	function/method name	Float of the execution time in seconds rounded to 3 decimal places	N/A	N/A

Class Name	File	Description	Parent	Child	Author
Setup_Interface	Setup_PodInterface.py	Setup_Interface provides the basic interface of required methods for subclasses to implement. SetupPodDevices.py is designed to handle any of these children.	N/A	Setup_E206HR	Thresa Kelly
Imports Name	Origin	Description	From		
os	Environment	For file path handling.	pyedflib		
os EdfWriter	Environment	For writing to EDF files.	threading		
Thread	Environment	For streaming from multiple POD devices.	io		
IOBase	Environment	For return annotations for text file operations.	SerialCommunication		
COM_io	Local	For getting available COM ports.	BasicPodProtocol		
POD_Basics	Local	For annotating POD devices as function parameters.			
Variables Name	Scope	Description	Type		
_NAME	Class	Device name, should be overwritten by child subclasses.	'GENERIC'		
PORTKEY	Class	Dictionary key for the COM port.	'Port'		
_podDevices	Instance	Dict of pod device objects. MUST have keys as device#	dict[int:POD_Basics]		
_podParametersDict	Instance	dictionary of device information. MUST have keys as device#, and each value must have [_PORTKEY , str...other values...].	dict[int:dict]		
_saveFileName	Instance	string filename: <path>/file.txt. The device name and number will be appended to the filename	str		
Methods Name	Type	Description	Parameter Name	Parameter Purpose	Return
_GetParam_onePODdevice	Instance	(Interface) Prompts the user to input all device setup parameters	forbiddenNames: list[str]	List of port names that are already used.	Dictionary of the device parameters.
_GetPODdeviceParameterTable	Instance	(Interface) get a text table that displays the parameters of all POD devices.	N/A	N/A	N/A
_ConnectPODdevice	Instance	(Interface) Write setup commands to initialize the POD device with the user's parameters	deviceNum: int deviceParams: dict	Integer key for the device# dictionary of the device parameters.	True for successful connection, false otherwise
_StreamThreading	Instance	(Interface) Stream data and save data to a file. Each POD device has its own thread	N/A	N/A	dictionary with the key as the device# and value as the thread object
_StopStream	Instance	(Interface) Tell POD devices to stop streaming	N/A	N/A	N/A
_OpenSaveFile_TXT	Static	(Interface) Open a text file and write column names	frame: str	String file name	opened file object IOBase
_OpenSaveFile_EDF	Instance	(Interface) Create an EDF file and write all channel information.	frame: str devNum: int	String file name Integer of the device#	EdfWriter file object
__init__	Dunder	Initializes the class instance variables	N/A	N/A	N/A
__del__	Dunder	Disconnects all POD devices.	N/A	N/A	N/A
SetFileName	Instance	Sets the filenames to save data to. Note that the device name and number will be appended to the end.	fileName: str	String file name	N/A
GetPODparametersDict	Instance	Gets a dictionary whose keys are the device number and the value is the device parameters dict.	N/A	N/A	N/A
SetupPODparameters	Instance	Sets the parameters for the POD devices.	podParametersDict: dict[None=None devices, name: str]	dictionary of the device parameters for all devices, Name or the POD device type.	N/A
_SetNumberOfDevices	Instance	Asks the user for how many devices they want to setup	N/A	N/A	True if all devices are successfully connected, false otherwise.
_ConnectAllPODDevices	Instance	Connects all POD devices by deleted all POD objects.	N/A	N/A	N/A
_DisconnectAllPODdevices	Instance	Disconnects all POD devices by deleted all POD objects.	N/A	N/A	N/A
_AddPODdevice	Instance	Asks the user for the parameters for the new device. A new device# is generated.	N/A	N/A	N/A
_SetParam_allPODdevices	Instance	First gets the number of POD devices, then asks the user for the information for each device.	N/A	N/A	N/A
_ChoosePort	Static	Asks the user to select a COM port.	forbiddenNames: list[] = []	List of port names that are already used.	String name of the port.
_GetPortsList	Static	Gets the names of all available ports.	forbiddenNames: list[] = []	List of port names that are already used.	List of port names
_ValidateParams	Instance	Displays a table of the parameters of all devices, then asks the user if everything is correct. The user can then edit the parameters of a device.	N/A	N/A	N/A
_EditParams	Instance	Asks the user which device to edit, and then asks them to re-input the device parameters	N/A	N/A	N/A
_SelectPODdeviceFromDictToEdit	Instance	Asks the user to select a valid device number. The input must be an integer number of an existing device.	key:str=Port	String key to access the _podParametersDict	Integer for the device#
_GetForbiddenNames	Instance	Generates a list of port names used by the active pod devices. There is an option to exclude an additional name from the list.	excludeStr=None e	String port name to exclude from the returned list	List of string names of ports in use.
_PrintDeviceNumber	Instance	Prints a title with the device#	num: int	Integer of the device#	N/A

<code>_DisplayPODDeviceParameters</code>	Instance	Display all the pod device parameters in a table	N/A	N/A	N/A	N/A
<code>_OpenSaveFile</code>	Instance	Opens a save file for a given device	devNum: int	Integer of the device#	Open OBase for a text file, or EdtWriter for EDt file.	N/A
<code>_BuildFileName</code>	Instance	Appends the device name and number to the end of the file name,	devNum: int	Integer of the device#	String file name.	N/A
<code>_Stream</code>	Instance	Tests that all devices are connected then starts streaming data	N/A	N/A	Dictionary with integer device# keys and Thread values.	Test connection failed.
<code>_TestDeviceConnection</code>	Instance	Writes a PING packet, then reads the response. A connection is successful if PING is read back	pod: POD_Basics	POD device	True for successful connection, false otherwise	N/A
<code>_TestDeviceConnection_All</code>	Instance	Tests the connection of all POD devices	N/A	N/A	True when all devices are successfully connected, false otherwise	N/A
<code>_AskYN</code>	Static	Asks the user a yes or no question	question: str	String containing the question	True for yes, False otherwise.	N/A

Class Name	File	Description	Parent	Child	Author		
Setup_8206HR	Setup_8206HR RPy	Setup_8206HR provides the setup functions for an 8206-HR POD device.	Setup_Interface	N/A	Theresa Kelly		
Name	Origin	Description	From	As			
os	Environment	For displaying the parameters in a table.					
numpy	Environment	For name handling.					
Thread	Environment	For streaming from multiple devices simultaneously.	threading				
IOWriter	Environment	For writing to EDF files.	pyedf				
IOWriter	Environment	For return annotations for text files.	io				
Setup_Interface	Local	For inheritance.	Setup_PodInterface				
POD_8206HR	Local	For communicating with 8206-HR POD devices.	PodDevice_8206HR				
datetime	Local	For implementing the current date and the time.	datetime				
time	Local	For implementing the current time of the execution.	datetime				
date	Local	For implementing the current date.	datetime				
gmtime	Local	For implementing the GMT time.	time				
Variables	Name	Scope	Description	Type			
-PARAMKEYS	Class	List of dictionary keys for device parameters	[Setup_Interface, PORTKEY, Sample Rate, Preamplifier Gain, Low Pass]	list[str]			
-LOWPASSKEYS	Class	List of dictionary keys for the Low Pass parameter.	[EEG1, EEG2, EEG3, EMG1]	list[str]			
-PHYSICAL_BOUND_UV	Class	Physical max-min stream value in uV	4069	int			
-NAME	Class	Name of the POD device. Overwritten from Parent	8206-HR	str			
Methods	Name	Type	Description	Parameter Name	Parameter Purpose	Return	Exception
_ConnectPODdevice	Instance	Creates a POD_8206HR object and write the setup parameters to it.	deviceNum: int deviceParams: dict[str,int dict[str,int]] forbiddenNames: list[str]	Integer of the device#	Dictionary of the device#s parameters	True of connection was successful, false otherwise.	N/A
_GetParam_onePODdevice	Instance	Asks the user to input all the device parameters			List of port names already used by other devices	Dictionary of device parameters	N/A
_ChooseSampleRate	Static	Asks user for the sample rate.	N/A			Integer number between 100-2000 Hz for the sample rate	Sample rate must be an integer between 100-2000
_ChoosePreampGain	Static	Asks user for the preamplifier gain of their POD device	N/A			Integer 10 or 100 for the preamp/gain	Gain must be an integer value of 10 or 100
_ChooseLowpass	Static	Builts dictionary of all lowpass filters	N/A			Dictionary containing lowpass filters (EEG1, EEG2, EEG3, EMG)	N/A
_ChooseLowpassForEEG	Static	Asks user for lowpass value for a given EEG	eeg: str			Integer number between 11-500 Hz for EEG	User input must be an integer between 11-500
_DisplayPODdeviceParameters	Instance	Prints a table containing the parameters for all POD devices	N/A			N/A	N/A
_OpenSaveFile_TXT	Static	Opens a text file and write the column names. Writes Currentdate and current time at the top of the txt file. Shows GMT.	frame: str	String filename	Opened file		N/A
_OpenSaveFile_EDF	Instance	Opens EDF file and write header	frameNum: int	String filename Integer device number	Opened file		N/A
_WriteDataToFile_TXT	Static	Writes data to an open text file	file: IOBase data: list[Pod.IODarray]	opened write file List of 3 items, one for each channel		N/A	N/A
_WriteDataToFile_EDF	Static	Writes data to an open EDF file	t: np.ndarray file: EdfWriter	Integer sample rate in Hz list with the time stamps (in seconds) opened EDF file		N/A	N/A
_StreamThreading	Instance	Opens a save file, then creates a thread for each device to stream and write data from.	data: list[Pod.IODarray]	List of 3 items, one for each channel		Dictionary with keys as the device# and values as the started Thread.	N/A
-StreamUnitStop	Instance	Streams data from a POD device and saves data to file. Stops looking when a stop stream command is read. Calculates average time difference across multiple packets to collect a continuous time series data.	pod: POD_8206HR file: IOBaseEdWriter	POD device open file		N/A	N/A
-StopStream	Instance	Write a command to stop streaming data to all POD devices	sampleRate: int	Integer sample rate in Hz		N/A	N/A
-UV	Static	Converts volts to microVolts, rounded to 6 decimal places	voltage: float/int	number of uVs		number of uVs	N/A

Class Name	File	Description	Parent	Child	Author
POD_8206HR	PodDevice_8206HR.py	Handles communication using an 8206HR POD device.	POD_Basics	N/A	Theresa Kelly
Imports Name			From		
POD_Basics	Local	For inheritance.	BasicPodProtocol		
POD_Packets	Local	For handling POD packets	PodPacketHardcoding		
POD_Commands	Local	For command constants	PodCommands		
Variables Name	Scope	Description	Value	Type	
_B4LENGTH	Class	Constant containing the number of bytes for a full Binary4 packet	16	int	
_B4BINARYLENGTH	Class	Constant containing the number of bytes for binary data in a Binary4 packet	8	int	
preampGain	Instance	Preamplifier gain	10 or 100	int	
Methods Name	Type	Description	Parameter Name	Parameter Purpose	Return
__init__	Dunder	Runs when an instance is constructed. It runs the parent's initialization. Then it updates the __commands to contain the appropriate commands for an 8306HR POD device.	port: str/int preampGain: int baudrate: int=9600	String of the serial port to be opened. Used when initializing the CON_Lo instance. Preamplifier gain. Must be 10 or 100. Integer baud rate of the opened serial port. Used when initializing the COM_M_lo instance.	N/A
UnpackPODpacket_Binary	Static	Overwrites the parent's method. Separates the components of a binary4 packet into a dictionary.	msg: bytes	Bytes string containing a complete binary4 Pod packet. STX (1 byte) + command (4 bytes) + packet number (1 bytes) + TTL (1 byte) + ch0 (2 bytes) + ch1 (2 bytes) + ch2 (2 bytes) + checksum (2 bytes) + ETX (1 byte)	An exception is raised if (1) the packet does not have the minimum number of bytes, (2) does not begin with STX, or (3) does not end with ETX.
TranslatePODpacket_Binary	Static	Overwrites the parent's method. Unpacks the binary4 POD packet and converts the values of the ASCII-encoded bytes into integer values and the values of binary-encoded bytes into integers. Channel values are given in Volts.	msg: bytes	Bytes string containing a complete binary4 Pod packet. STX (1 byte) + command (4 bytes) + packet number (1 bytes) + TTL (1 byte) + ch0 (2 bytes) + ch1 (2 bytes) + ch2 (2 bytes) + checksum (2 bytes) + ETX (1 byte)	A dictionary containing 'Command Number', 'Packet #', 'TTL', 'Ch0', 'Ch1', and 'Ch2' in bytes.
TranslatePODpacket	Instance	Overwrites the parent's method. Determines if the packet is standard or binary, and translates accordingly. Adds a check for the 'GET TTL PORT' command.	msg: bytes	Bytes string containing either a standard or binary packet	A dictionary containing the unpacked message in numbers
TranslateTTLbyte_ASCII	Static	Separates the bits of each TTL (0-3) from a byte.	ttlByte: bytes	One Byte string for the TTL (ASCII encoded)	N/A
TranslateTTLbyte_Binary	Static	Separates the bits of each TTL (0-3) from a byte.	ttlByte: bytes	One Byte string for the TTL (binary encoded)	N/A
_BinaryBytesToVoltage	Instance	Converts a binary bytes value read from POD device and converts it to the real voltage value at the preamplifier input	value: bytes	Bytes string containing voltage measurement	N/A
Read_Binary	Instance	After receiving the prePacket, it reads the 8 bytes(TTL-channels) and then reads to ETX (checksum+ETX).	prePacket: bytes validateChecksum:bool=True	Bytes string containing the beginning of a POD packet. STX (1 byte) + command number (4 bytes) Set to True to validate the checksum. Set to False to skip validation	Bytes string for a binary4 POD packet.

Class Name	File	Description	Parent	Child	Author
POD_8401HR	PodDevice_8401HR.py	Handles communication using an 8401HR POD device.	POD_Basics	N/A	Thresa Kelly
<i>Imports</i>					
Name	Origin	Description	From		
POD_Basics	Local	For inheritance	BasicPodProtocol		
POD_Packets	Local	For handling POD packets	PodPacketHandling		
POD_Commands	Local	For command constants	PodCommands		
<i>Variables</i>					
Name	Scope	Description	Value	Type	
_BSLENGTH	Class	number of bytes for a Binary 5 packet	31	int	
_BSBNARYLENGTH	Class	number of binary bytes for a binary 5 packet	23	int	
_channelMap	Instance	Dictionary of the channel tables	Set by GetChannelMapping(). Dictionary keys are [A,B,C,D].	dict[str,int]None	An exception is raised if (1) the ssGain or preampGain have improper keys, (2) the device/sensor does not exist, (3) the ssGain was given bad values, or (4) the preampGain was given bad values
_ssGain	Instance	Dictionary of the second stage gain for all four channels	1, 5, or None. Dictionary keys are [A,B,C,D].	dict[str,int]None	
_preampGain	Instance	Dictionary of the preamplifier gain for all four channels.	10, 100, or None. Dictionary keys are [A,B,C,D].	dict[str,int]None	
<i>Methods</i>					
Name	Type	Description	Parameter Name	Parameter Purpose	Exception
__init__	Dunder	Runs when an instance is constructed. It runs the parent's initialization. Then it updates the _commands to contain the appropriate commands for an 8401HR POD device. Sets the _channelMap, _ssGain, and _preampGain.	port: str/int deviceName: str ssGain:dict[int,int]=[A:None, B:None,C:None,D:None] preampGain:dict[int,int]=[A:None, B:None,C:None,D:None] baudrate:int=9600	String of the serial port to be opened. Used when initializing the COM_Io instance. String of the corresponding device/sensor name Dictionary of the secondary stage gain Dictionary of the preamplifier gain Integer baud rate of the opened serial port. Used when initializing the COM_Io instance.	
UnpackPODpacket_Binary	Static	Overwrites the parent's method. Separates the components of a binary 5 packet into a dictionary.	msg: bytes	Bytes string containing a complete binary5 Pod packet: STX (1 byte) + command (4) + packet number (1) + status (1) + channels (9) + analog inputs (12) + checksum (2) + ETX (1)	A dictionary containing 'Command Number', 'Packet #', 'Status', 'Channels', 'Analog EX10', 'Analog EX11', 'Analog TTL1', 'Analog TTL2', 'Analog EXT1', 'Analog TTL3', 'Analog TTL4'.
TranslatePODpacket_Binary	Instance	Overwrites the parent's method. Unpacks the binary5 POD packet and converts the values of the ASCII encoded bytes into integer values and the values of binary-encoded bytes into integers. The channels and analogs are converted to volts (V).	msg: bytes	Bytes string containing a complete binary5 Pod packet: STX (1 byte) + command (4) + packet number (1) + status (1) + channels (9) + analog inputs (12) + checksum (2) + ETX (1)	A dictionary containing 'Command Number', 'Packet #', 'Status', 'CH3', 'CH2', 'CH1', 'CHO', 'Analog EX1', 'Analog EXT1', 'Analog TTL1', 'Analog TTL2', 'Analog TTL3', 'Analog TTL4'. N/A
TranslatePODpacket	Instance	Overwrites the parent's method. Determines the packet is standard or binary and translates the accordingly. Specially handles TTL packet payloads.	msg: bytes	Bytes string containing either a standard or binary packet.	A dictionary containing the unpacked message in numbers
GetChannelMapping	Static	Get the channel mapping (channel labels for A,B,C,D) for a given device.	device: str	String for the device/sensor name.	Dictionary with keys A,B,C,D with values of the channel names. Returns None if the device name does not exist.
_Voltage_PrimaryChannels	Static	Builds an integer, which represents a binary mask, that can be used for TTL command arguments.	ext0:bool=0 ext1:bool=0 ttl4:bool=0 ttl3:bool=0 ttl2:bool=0	boolean bit boolean bit boolean bit boolean bit boolean bit	Integer number to be used as a bit mask
GetTTIbitmask_Int			value: int	Value to be converted to voltage	Number of the voltage in volts [V]. Returns value if no gain is given (mc-connect).
_Voltage_PrimaryChannels_EEGEMG	Static	Converts a value to a voltage for a primary channel.	ssGain:int[None=None PreampGain:int[None=None value:int	Second stage gain Preamplifier gain Value to be converted to voltage	N/A
_Voltage_PrimaryChannels_Biosensor	Static	Converts a value to a voltage for a biosensor primary channel.	value:int ssGain: int	Second stage gain Preamplifier gain Value to be converted to voltage	N/A
_Voltage_SecondaryChannels	Static	Converts a value to a voltage for a secondary channel.	value:int	Second stage gain Value to be converted to voltage	N/A
_Read_Binary	Instance	After receiving the prePacket, it reads the 23 bytes (binary data) and then reads to ETX (checksum+ETX).	prePacket: bytes validateChecksum:bool=True	Bytes string containing the beginning of a POD packet: STX (1 byte) + command number (4 bytes) Set to True to validate the checksum. Set to False to skip validation	N/A

Class Name	File	Description	Parent	Child	Author
POD_Packets	PodPacket	Collection of methods for creating and interpreting POD packets	N/A	N/A	Thresa Kelly
Imports Name		PodPacket Handling API			
Methods Name	Origin	Description	From		
N/A	N/A		N/A		
Methods Name	Type	Description	Parameter Name	Parameter Purpose	Return
STX	Static	Get STX in bytes. STX marks the starting byte of a POD Packet	N/A	N/A	Bytes or STX(0x02)
ETX	Static	Get ETX in bytes. ETX marks the end byte of a POD Packet	N/A	N/A	Bytes or ETX(0x03)
		Converts an integer value into ASCII-encoded bytes.		Integer value to be converted into ASCII-encoded bytes	
		First, it converts the integer value into a usable uppercase hexadecimal string. Then it converts the ASCII code for each character into bytes. Lastly, it ensures that the final message is the desired length.			
		Example: if value=2 and numBytes=4, the returned ASCII will show b'0002', which is 0x30 0x30 0x30 0x32' in bytes.		Number bytes to be the length of the ASCII-encoded message.	N/A
IntToAsciiBytes	Static	Converts a POD-encoded bytes message into an integer. It does this using a base-16 conversion.	msg: bytes	Bytes message to be converted to an integer. The bytes must be base-16 or the conversion will fail.	Integer result from the ASCII-encoded bytes conversion.
AsciiBytesToInt	Static	Converts binary-encoded bytes into an integer	msg: bytes	Bytes message holding binary information to be converted into an integer.	N/A
BinaryBytesToInt	Static	Converts a specific bit range in an ASCII-encoded bytes object to an integer.	byteorder:str='big'	Ordering of bytes, 'big' for big endian and 'little' for little endian.	Integer result from the binary-encoded bytes message.
ASCIIBytesToIntInt_Split	Static	Converts a specific bit range in an ASCII-encoded bytes object to an integer.	signedBool=False	Boolean flag to mark if the msg is signed (True) or unsigned (False)	Integer result from the ASCII-encoded bytes message in a given bit range.
			msg: bytes	Bytes message holding binary information to be converted into an integer.	
			keepTopBits: int	Integer position of the msb of desired bit range	
			cutBottomBits: int	Integer number of lsbs to remove	
			msg: bytes	Bytes message holding binary information to be converted into an integer.	
			keepTopBits: int	Integer position of the msb of desired bit range	
			cutBottomBits: int	Integer number of lsbs to remove	
			byteorder:str='big'	Ordering of bytes, 'big' for big endian and 'little' for little endian.	
			signedBool=False	Boolean flag to mark if the msg is signed (True) or unsigned (False)	
Checksum		Calculates the checksum of a given bytes message. This is achieved by summing each byte in the message, inverting, and taking the last byte.	bytesIn: bytes	Bytes message containing POD packet data	Two ASCII-encoded bytes containing the checksum for bytesIn
BuildPODpacket_Standard		Builds a standard POD packet – STX (1 byte) + command number (4 bytes) + optional packet (?) bytes) + checksum (2 bytes) + ETX (1 bytes) – as bytes.	commandNumber: int payload: bytes None=None payload: intbytes tuple[intbytes] argsizes: tuple[int]	Integer representing the command number. This will be converted into a 4 byte long ASCII-encoded bytes string. bytes string containing the payload	Bytes string of a complete standard POD packet
PayloadToIntBytes	Static	Converts a payload into a bytes string	integer: bytes, or tuple containing the payload	Bytes string of the payload	Bytes string of the argument sizes

Class Name	File	Description	Parent	Child	Author
POD_Basics	BasicPodPr	Handle basic communication with a POD device, including reading and writing packets and packet interpretation.	N/A	POD_8206HR POD_8401THR	Theresa Kelly
Imports Name					
COM_Io	Origin	Description	From		
POD_Packets	Local	For opening and connecting serial COM ports	SerialCommunication		
POD_Commands	Local	For handling POD packets	PodPacketHandling		
Variables Name	Scope	Description	Value	Type	
_numPod	Class	Integer equal to the number of POD_Basics class instances, incremented on construction and decremented on destruction	0	int	
_MINSTANDARDLENGTH	Class	Integer minimum number of bytes in a standard POD packet	8	int	
_MINBINARYLENGTH	Class	Integer minimum number of bytes in a binary POD packet	15	int	
port	Instance	Open serial port via COM_Io class instance	COM_Io	COM_Io	
commands	Instance	Command handler POD_Commands class instance	POD_Commands	POD_Commands	
Methods Name	Type	Description	Parameter Name	Parameter Purpose	Return
init	Dunder	Runs when an instance of POD_Basics is constructed. It initializes the instance variable for the COM port communication, _port, and for the command handler, _commands. It also increments the POD device counter, _NUMPOD.	port: str/int baudrate:int=9600	String of the serial port to be opened. Used when initializing the COM_Io instance. Integer baud rate of the opened serial port. Used when initializing the COM_Io instance.	N/A
del	Dunder	Runs when an instance is destroyed. It decrements the POD device counter, _NUMPOD.	N/A	N/A	N/A
GetNumberOfPODDevices	Static	Get the POD device counter	N/A	N/A	Integer of the number of class instances (_NUMPOD).
UnpackPODpacket_Standard	Static	Converts a standard POD packet into a dictionary containing the command number and payload (if applicable) in bytes.	msg: bytes	Bytes message containing a standard POD packet: STX (1 byte) + command number (4 bytes) + optional packet (7 bytes) + checksum (2 bytes) + ETX (1 bytes)	A dictionary containing the POD packet's 'Command Number' and 'Payload' (if applicable) in bytes.
UnpackPODpacket_Binary	Static	Converts a variable-length binary packet into a dictionary containing the command number, binary packet length, and binary data in bytes.	msg: bytes	Bytes message containing a variable-length POD packet: STX (1 byte) + command number (4 bytes) + length of binary (4 bytes) + checksum (2 bytes) + ETX (1 bytes) + binary (LENGTH bytes) + checksum (2 bytes) + ETX (1 bytes)	A dictionary containing the 'Command Number', Binary Packet Length, and Binary Data in bytes.
TranslatePODpacket_Standard	Instance	Unpacks the standard POD packet and converts the ASCII-encoded bytes values into integer values.	msg: bytes	Bytes message containing a standard POD packet	A dictionary containing the POD packet's 'Command Number' and 'Payload' (if applicable) in integers.
TranslatePODpacket_Binary	Static	Unpacks the variable-length binary POD packet and converts the values of the ASCII-encoded bytes into integer values and leaves the binary-encoded bytes as is.	msg: bytes	Bytes message containing a variable-length POD packet	A dictionary containing the 'Command Number' and 'Binary Data' in bytes.
_ValidateChecksum	Static	Validates the checksum of a given POD packet. The checksum is valid if the calculated checksum from the data matches the checksum written in the packet.	msg: bytes	Bytes message containing a POD packet: STX (1 bytes) + data (?) bytes + checksum (2 bytes) + ETX (1 bytes).	Returns True if the checksum is correct, False otherwise.
GetDeviceCommands	Instance	Gets the dictionary containing the class instances available POD commands.	N/A	N/A	An exception is raised if the msg does not begin with STX or end with ETX.
SetBaudratePODdevice	Instance	If the port is open, it will change the baud rate to the parameter's value	baudrate: int	Integer baud rate to set for the open serial port.	An exception is raised if (1) the command does not have the minimum number of bytes in a standard pod packet, (2) does not begin with STX, and (3) does not end with ETX.
UnpackPODpacket	Static	Determines if the packet is standard or binary, and unpacks accordingly.	msg: bytes	Bytes string containing either a standard or binary packet	An exception is raised if (1) the msg does not have the minimum number of bytes in a standard pod packet, (2) does not begin with STX, (3) does not end with ETX, and (4) does not have an ETX after standard packet.
TranslatePODpacket	Instance	Determines if the packet is standard or binary, and translates accordingly.	msg: bytes	Bytes string containing either a standard or binary packet	An exception is raised if (1) the command does not have the minimum number of bytes in a standard pod packet, (2) does not begin with STX, and (3) does not end with ETX.
WriteRead	Instance	Writes a command with optional payload to POD device, then reads (once) the device response.	cnd: str/int validateChecksum:bool=True	Bytes string containing a POD packet beginning with STX and ending with ETX. This may be a standard packet, binary packet, or an unformatted packet (STX+something+ETX).	Set to True to validate the checksum. Set to False to skip validation
		Builds a POD packet and writes it to a POD device via	cnd: str/int	An integer representing the command number.	An exception is raised if (1) the command does not

GetPODpacket	Instance	Buils a POD packet and writes it to a FUD uerface via COM port. If an integer payload is give, the method will convert it into a bytes string of the length expected by the command. If a bytes payload is given, it must be the correct length.	cmd: str/int payload:int/bytestuple[int [bytes]=None]	None when there is no payload. If there is a payload, set to an integer value, bytes string, or tuple	Returns the bytes string of the POD packet.
WritePacket	Instance	Buils a POD packet and writes it to the POD device.	cmd: str/int payload:int/bytestuple[int [bytes]=None]	None when there is no payload. If there is a payload, set to an integer value, bytes string, or tuple	Returns the bytes string that was written to the POD device
ReadPODpacket	Instance	Reads a complete POD packet, either in standard or binary format, beginning with STX and ending with ETX. Reads first STX and then starts recursion.	validateChecksum:bool= True	Set to True to validate the checksum. Set to False to skip validation	Bytes string containing a POD packet beginning with STX and ending with ETX. This may be a standard packet, binary packet, or an unformatted packet (STX+something+ETX).
_ReadPODpacket_Recursive	Instance	Reads the command number. If the command number ends in ETX, the packet is returned. Next, it checks if the command is allowed. Then, it checks if the command is standard or binary and reads accordingly, then returns the packet.	validateChecksum:bool= True	Set to True to validate the checksum. Set to False to skip validation	Bytes string containing a POD packet beginning with STX and ending with ETX. This may be a standard packet, binary packet, or an unformatted packet (STX+something+ETX).
_Read_GetCommand	Instance	Reads one byte at a time up to 4 bytes to get the ASCII-encoded bytes command number. For each byte read, it can (1) start the recursion over if an STX is found, (2) returns if ETX is found, or (3) continue building the command number.	validateChecksum:bool= True	Set to True to validate the checksum. Set to False to skip validation	4 byte long string containing the ASCII-encoded command number.
_Read_ToETX	Instance	Reads one byte at a time until an ETX is found. It will restart the recursive read if an STX is found anywhere.	validateChecksum:bool= True	Set to True to validate the checksum. Set to False to skip validation	Bytes string ending with ETX
_Read_Standard	Instance	Reads the payload, checksum, and ETX. Then it builds the complete standard POD packet in bytes.	prePacket: bytes validateChecksum:bool= True	Bytes string containing the beginning of a POD packet. STX (1 byte) + command number (4 bytes)	An exception is raised if the checksum is invalid (only if validateChecksum=True)
_Read_Binary	Instance	Reads the remaining part of the variable-length binary packet. It first reads the standard packet (prePacket+payload+checksum+ETX). Then it determines how long the binary packet is from the payload of the standard POD packet and reads that many bytes. It then reads to ETX to get the checksum+ETX.	prePacket: bytes validateChecksum:bool= True	Bytes string containing the beginning of a POD packet. STX (1 byte) + command number (4 bytes)	An exception is raised if the checksum is invalid (only if validateChecksum=True)

Class Name	File	Description	Parent	Child	Author
<code>POD_Commands</code>	<code>PodCommands.py</code>	Manages a dictionary containing available commands for a POD device.	N/A	N/A	Tressa Kelly
Imports	Name	Description	From		
<code>__NAME</code>	Class	index key for the command name for __commands dict	int		
<code>__ARGUMENTS</code>	Class	index key for the number of bytes in an argument for __commands dict values	int		
<code>__RETURNS</code>	Class	index key for the number of bytes in the return for __commands dict values	int		
<code>__BINARY</code>	Class	index key for the binary flag for __commands dict	int		
<code>__NOVALUE</code>	Class	Integer used to mark when a list item in __commands means no value or undefined.	int		
<code>U8</code>	Class	Number of bytes for an unsigned 8-bit value	int		
<code>U16</code>	Class	Dictionary containing the available commands for a POD device. Each entry is formatted as {key(command number) : value(command name, number of argument ASCII bytes, number of return bytes, binary flag)}	POD_Commands.GetBasicCommands()	dict[int,list[tuple[int,bool]]]	
Methods	Name	Type	Description	Parameter Name	Parameter Purpose
<code>__init__</code>	Dunder	Runs when an instance is constructed. It sends the commands dictionary to the basic command set.	N/A	N/A	N/A
<code>NoValue</code>	Static	Gets value of __NOVALUE	N/A	N/A	Value of __NOVALUE
<code>U8</code>	Static	Gets value of __U8	N/A	N/A	Value of __U8
<code>U16</code>	Static	Gets value of __U16	N/A	N/A	Value of __U16
<code>GetBasicCommands</code>	Static	Creates a dictionary containing the basic POD command set [0,1,2,3,4,5,6,7,8,9,10,11,12]	N/A	N/A	N/A
<code>GetCommands</code>	Instance	Gets the contents of the current command dictionary (<code>__commands</code>)	N/A	N/A	N/A
<code>RestoreBasicCommands</code>	Instance	Sets the current command set (<code>__commands</code>) to the basic POD command set.	N/A	N/A	N/A
<code>AddCommand</code>	Instance	Adds a command entry to the current commands dictionary (<code>__commands</code>) if the command does not exist.	commandNumber: int commandName: str argumentBytes: tuple[int]	integer of the command number String of the command's name Integer of the number of bytes in the argument	True if the command was successfully added. False if the command could not be added because it already exists.
<code>RemoveCommand</code>	Instance	Removes the entry for a given command in <code>__commands</code> dictionary.	returnBytes: tuple[int] isBinary: bool cmd: int,str	Integer of the number of bytes in the return Boolean flag to mark if the command is binary (True) integer command number or string command name. string of the command's name	True if the command was successfully removed. False if the command does not exist.
<code>CommandNumberFromName</code>	Instance	Gets the command number from the command dictionary using the command's name	name: str	integer representing the command number. If the command could not be found, return None.	N/A
<code>ArgumentBytes</code>	Instance	Gets the tuple for the number of bytes in the argument for a given command.	cmd: int,str	Tuple representing the number of bytes in the argument for cmd. If the command could not be found, return None.	N/A
<code>ReturnBytes</code>	Instance	Gets the tuple for the number of bytes in the return for a given command.	cmd: int,str	Tuple representing the number of bytes in the return for cmd. If the command could not be found, return None.	N/A
<code>IsCommandBinary</code>	Instance	Gets the binary flag for a given command	cmd: int,str	Boolean flag that is True if the command is binary and False if standard. If the command could not be found, return None.	N/A
<code>DoesCommandExist</code>	Instance	Checks if a command exists in the __commands dictionary	cmd: int,str	True if the command exists, false otherwise.	N/A

Class Name	File	Description	Parent	Child	Author
COM_io	SerialCommunication.py	Handle serial communication (read/write) using COM ports.	N/A	N/A	Thresa Kelly
Imports Name	Origin	Description	From		
serial.tools.list_ports	Environment	For accessing the COM ports on the computer	N/A		
Variables Name	Scope	Description	Value	Type	
_serialInst	Instance	Serial object to set the port and baud rate to. It can be opened or closed.	Serial	serial.Serial	
Methods Name	Type	Description	Parameter Name	Parameter Purpose	Return
GetCOMPortsList	Static	Finds all the available COM ports on the user's computer and appends them to an accessible list.	N/A	N/A	List containing the names of available COM ports
init	Dunder	Runs when the object is initialized. It initializes the serialinst to a given COM port with a set baudrate.	port: str/int baudrate:int=9600	String of the serial port to be opened. Integer baud rate of the opened serial port.	N/A
del	Dunder	Runs when the object is destructed. It closes the serial port, if open.	N/A	N/A	N/A
_BuildPortName	Instance	Converts the port parameter into the "COM<#>" format	port: str/int	Name of a COM port. Can be an integer or string.	N/A
IsSerialOpen	Instance	Returns True if the serial instance port is open, false otherwise	N/A	N/A	N/A
IsSerialClosed	Instance	Returns False if the serial instance port is open, True otherwise	N/A	N/A	N/A
CloseSerialPort	Instance	Closes the instance serial port if it is open.	N/A	N/A	N/A
OpenSerialPort	Instance	First, it closes the serial port if it is open. Then, it opens a serial port with a set baud rate.	port: str/int baudrate:int=9600	String of the serial port to be opened. Integer baud rate of the opened serial port.	Raises an exception if the given port does not exist.
SetBaudrate	Instance	If the port is open, it will change the baud rate to the parameter's value	baudrate: int	Integer baud rate to set for the open serial port.	True if successful at setting the baud rate, false otherwise
GetPortName	Instance	Gets the name of the open port.	N/A	N/A	If the serial port is open, it will return a string of the port's name. If the port is closed, it will return None.
Read	Instance	Reads a specified number of bytes from the open serial port.	numBytes: int	Integer number of bytes to read	If the serial port is open, it will return a set number of read bytes. If it is closed, it will return None.
ReadLine	Instance	Reads until a new line ('\n') from the open serial port.	N/A	N/A	If the serial port is open, it will return a complete read line. If closed, it will return None.
ReadUntil	Instance	Reads until a set character from the open serial port.	eof: bytes	end-of-line character	If the serial port is open, it will return a read line ending in eof. If closed, it will return None.
Write	Instance	Write a set message to the open serial port.	message: bytes	byte string containing the message to write	N/A