

---

# Python POD API

*Release v1.3.2*

**Thresa Kelly**

**Aug 03, 2023**



**CONTENTS:**

|          |   |           |
|----------|---|-----------|
| <b>1</b> | <b>Python POD API Project</b>           | <b>1</b>  |
| 1.1      | PodApi package . . . . .                | 1         |
| 1.1.1    | Subpackages . . . . .                   | 1         |
| 1.1.1.1  | PodApi.Commands package . . . . .       | 1         |
| 1.1.1.2  | PodApi.Devices package . . . . .        | 5         |
| 1.1.1.3  | PodApi.Packets package . . . . .        | 26        |
| 1.1.1.4  | PodApi.Parameters package . . . . .     | 44        |
| 1.1.2    | Module contents . . . . .               | 51        |
| 1.2      | Setup package . . . . .                 | 51        |
| 1.2.1    | Subpackages . . . . .                   | 51        |
| 1.2.1.1  | Setup.Inputs package . . . . .          | 51        |
| 1.2.1.2  | Setup.SetupOneDevice package . . . . .  | 57        |
| 1.2.2    | Submodules . . . . .                    | 74        |
| 1.2.3    | Setup.Setup_PodDevices module . . . . . | 74        |
| 1.2.4    | Module contents . . . . .               | 77        |
| <b>2</b> | <b>Indices and tables</b>               | <b>79</b> |
|          | <b>Python Module Index</b>              | <b>81</b> |
|          | <b>Index</b>                            | <b>83</b> |



## PYTHON POD API PROJECT

### 1.1 PodApi package

#### 1.1.1 Subpackages

##### 1.1.1.1 PodApi.Commands package

##### Submodules

##### PodApi.Commands.PodCommands module

**class** PodApi.Commands.PodCommands.CommandSet

Bases: object

POD\_Commands manages a dictionary containing available commands for a POD device.

**\_\_commands**

Dictionary containing the available commands for a POD device. Each entry is formatted as { key(command number) : value([command name, number of argument ASCII bytes, number of return bytes, binary flag ) }.

**Type**

dict[int,list[str|tuple[int]]|bool]]

**AddCommand**(*commandNumber: int, commandName: str, argumentBytes: tuple[int], returnBytes: tuple[int], isBinary: bool, description: str*) → bool

Adds a command entry to the current commands dictionary (\_\_commands) if the command does not exist.

**Parameters**

- **commandNumber** (*int*) – Integer of the command number.
- **commandName** (*str*) – String of the command's name.
- **argumentBytes** (*tuple[int]*) – Integer of the number of bytes in the argument.
- **returnBytes** (*tuple[int]*) – Integer of the number of bytes in the return.
- **isBinary** (*bool*) – Boolean flag to mark if the command is binary (True) or standard (False).
- **description** (*str*) – String description of the command.

**Returns**

True if the command was successfully added, False if the command could not be added because it already exists.

**Return type**

bool

**ArgumentHexChar**(*cmd: int | str*) → tuple[int] | None

Gets the tuple for the number of hex characters in the argument for a given command.

**Parameters**

**cmd** (*int | str*) – Integer command number or string command name.

**Returns**

Tuple representing the number of bytes in the argument for cmd. If the command could not be found, return None.

**Return type**

tuple[int]|None

**CommandNumberFromName**(*name: str*) → int | None

Gets the command number from the command dictionary using the command's name.

**Parameters**

**name** (*str*) – String of the command's name.

**Returns**

Integer representing the command number. If the command could not be found, return None.

**Return type**

int|None

**Description**(*cmd: int | str*) → str | None

Gets the description for a given command.

**Parameters**

**cmd** (*int | str*) – Integer command number or string command name.

**Returns**

String description for the command. If the command could not be found, return None.

**Return type**

str|None

**DoesCommandExist**(*cmd: int | str*) → bool

Checks if a command exists in the \_\_commands dictionary.

**Parameters**

**cmd** (*int | str*) – Integer command number or string command name.

**Returns**

True if the command exists, false otherwise.

**Return type**

bool

**static GetBasicCommands**() → dict[int, list[str | tuple[int] | bool]]

Creates a dictionary containing the basic POD command set (0,1,2,3,4,5,6,7,8,9,10,11,12).

**Returns**

Dictionary containing the available commands for this POD device. Each entry is formatted as

{ key(command number) : value([command name, number of argument ASCII bytes, number of return bytes, binary flag, description) }.

**Return type**

dict[int, list[str|tuple[int]]|bool|str]]

**GetCommands()** → dict[int, list[str | tuple[int] | bool]]

Gets the contents of the current command dictionary (\_\_commands).

**Returns**

Dictionary containing the available commands for a POD device. Each entry is formatted as { key(command number) : value([command name, number of argument ASCII bytes, number of return bytes, binary flag, description) }.

**Return type**

dict[int, list[str|tuple[int]]|bool|str]]

**IsCommandBinary(cmd: int | str)** → bool | None

Gets the binary flag for a given command.

**Parameters**

**cmd** (int | str) – Integer command number or string command name.

**Returns**

Boolean flag that is True if the command is binary and False if standard. If the command could not be found, return None.

**Return type**

bool|None

**static NoValue()** → int

Gets value of \_\_NOVALUE.

**Returns**

Value of \_\_NOVALUE.

**Return type**

int

**RemoveCommand(cmd: int | str)** → bool

Removes the entry for a given command in \_\_commands dictionary.

**Parameters**

**cmd** (int | str) – Integer command number or string command name.

**Returns**

True if the command was successfully removed, False if the command does not exist.

**Return type**

bool

**RestoreBasicCommands()** → None

Sets the current commands (\_\_commands) to the basic POD command set.

**ReturnHexChar(cmd: int | str)** → tuple[int] | None

Gets the tuple for the number of hex characters in the return for a given command.

**Parameters**

**cmd** (int | str) – Integer command number or string command name.

**Returns**

Tuple representing the number of hex characters in the return for cmd. If the command could not be found, return None.

**Return type**

tuple[int]|None

**Search**(cmd: int | str, idx: int | None = None) → str | tuple[int] | bool | None

Searches the \_\_commands dictionary for the command.

**Parameters**

- **cmd** (int | str) – Integer command number or string command name.
- **idx** (int, optional) – Index for the desired value in the command information list. Defaults to None.

**Returns**

If an idx was given, this returns the idx value of the command information list if the command was found (None otherwise). If no idx is given, this returns true if the command is found (False otherwise).

**Return type**

str|tuple[int]|bool|str|None

**static U16**() → int

Gets value of \_\_U16.

**Returns**

Value of \_\_U16.

**Return type**

int

**static U32**() → int

Gets value of \_\_U32.

**Returns**

Value of \_\_U32.

**Return type**

int

**static U8**() → int

Gets value of \_\_U8.

**Returns**

Value of \_\_U8.

**Return type**

int

**\_\_ARGUMENTS: int = 1**

Class-level integer representing the index key for the number of bytes in an argument for \_\_commands list values.

**\_\_BINARY: int = 3**

Class-level integer representing the index key for the binary flag for \_\_commands list values.

**\_\_DESCRIPTION: int = 4**

Class-level integer representing the index key for the description for \_\_commands list values.



**\_\_NAME: int = 0**

Class-level integer representing the index key for the command name for \_\_commands list values.

**\_\_NOVALUE: int = -1**

Class-level integer used to mark when a list item in \_\_commands means ‘no value’ or is undefined.

**\_\_RETURNS: int = 2**

Class-level integer representing the index key for the number of bytes in the return for \_\_commands list values.

**\_\_U16: int = 4**

Class-level integer representing the number of hexadecimal characters for an unsigned 16-bit value.

**\_\_U32: int = 8**

Class-level integer representing the number of hexadecimal characters for an unsigned 32-bit value.

**\_\_U8: int = 2**

Class-level integer representing the number of hexadecimal characters for an unsigned 8-bit value.

## Module contents

### 1.1.1.2 PodApi.Devices package

#### Subpackages

#### PodApi.Devices.SerialPorts package

#### Submodules

#### PodApi.Devices.SerialPorts.PortAccess module

**class PodApi.Devices.SerialPorts.PortAccess.FindPorts**

Bases: object

**static ChoosePort**(*forbidden: list[str] = []*) → str

Systems checks user’s Operating System, and chooses ports accordingly.

#### Parameters

**forbidden** (*list[str], optional*) – List of port names that the user should not use. This may be because these ports are already in use or that the port is not a POD device. Defaults to [].

#### Returns

String name of the port.

#### Return type

str

**static GetAllPortNames**() → list[str]

Finds all the available COM ports on the user’s computer and appends them to an accessible list.

#### Returns

List containing the names of available COM ports.

#### Return type

list[str]

**static GetSelectPortNames**(*forbidden: list[str] = []*) → list[str]

Gets the names of all available ports.

**Parameters**

**forbidden**(*list[str], optional*) – List of port names that the user should not use. This may be because these ports are already in use or that the port is not a POD device. Defaults to [].

**Returns**

List of port names.

**Return type**

list[str]

**static \_ChoosePortLinux**(*forbidden: list[str] = []*) → str

User picks Serial port in Linux.

**Parameters**

**forbidden**(*list[str], optional*) – List of port names that the user should not use. This may be because these ports are already in use or that the port is not a POD device. Defaults to [].

**Returns**

String name of the port.

**Return type**

str

**static \_ChoosePortWindows**(*forbidden: list[str] = []*) → str

User picks COM port in Windows.

**Parameters**

**forbidden**(*list[str], optional*) – List of port names that the user should not use. This may be because these ports are already in use or that the port is not a POD device. Defaults to [].

**Returns**

String name of the port.

**Return type**

str

## PodApi.Devices.SerialPorts.SerialComm module

**class PodApi.Devices.SerialPorts.SerialComm.PortIO**(*port: str | int, baudrate: int = 9600*)

Bases: object

COM\_io handles serial communication (read/write) using COM ports.

**\_\_serialInst**

Instance-level serial COM port.

**Type**

Serial

**CloseSerialPort**() → None

Closes the instance serial port if it is open.

**Flush()** → bool

Reset the input and output serial buffer.

**Returns**

True if the buffers are flushed, False otherwise.

**Return type**

bool

**GetPortName()** → str | None

Gets the name of the open port.

**Returns**

If the serial port is open, it will return a string of the port's name. If the port is closed, it will return None.

**Return type**

str|None

**IsSerialClosed()** → bool

Returns False if the serial instance port is open, True otherwise.

**Returns**

True if the COM port is closed, False otherwise.

**Return type**

bool

**IsSerialOpen()** → bool

Returns True if the serial instance port is open, false otherwise.

**Returns**

True if the COM port is open, False otherwise.

**Return type**

bool

**OpenSerialPort**(*port: str | int, baudrate: int = 9600*) → None

First, it closes the serial port if it is open. Then, it opens a serial port with a set baud rate.

**Parameters**

- **port** (*str | int*) – String of the serial port to be opened.
- **baudrate** (*int, optional*) – Integer baud rate of the opened serial port. Defaults to 9600.

**Raises**

**Exception** – Port does not exist.

**Read**(*numBytes: int, timeout\_sec: int | float = 5*) → bytes | None

Reads a specified number of bytes from the open serial port.

**Parameters**

- **numBytes** (*int*) – Integer number of bytes to read.
- **timeout\_sec** (*int | float, optional*) – Time in seconds to wait for serial data. Defaults to 5.

**Raises**

**Exception** – Timeout for serial read.

**Returns**

If the serial port is open, it will return a set number of read bytes. If it is closed, it will return None.

**Return type**

bytes|None

**ReadLine()** → bytes | None

Reads until a new line is read from the open serial port.

**Returns**

If the serial port is open, it will return a complete read line. If closed, it will return None.

**Return type**

bytes|None

**ReadUntil(*eol: bytes*)** → bytes | None

Reads until a set character from the open serial port.

**Parameters**

**eol** (*bytes*) – end-of-line character.

**Returns**

If the serial port is open, it will return a read line ending in eol. If closed, it will return None.

**Return type**

bytes|None

**SetBaudrate(*baudrate: int*)** → bool

Sets the baud rate of the serial port

**Parameters**

**baudrate** (*int*) – Baud rate, or signals per second.

**Returns**

True if the baudrate was set, False otherwise.

**Return type**

bool

**Write(*message: bytes*)** → None

Write a set message to the open serial port.

**Parameters**

**message** (*bytes*) – byte string containing the message to write.

**\_\_BuildPortName(*port: str | int*)** → str

Converts the port parameter into the “COM”+<number> format for Windows or “/dev/tty...”+<number> for Linux.

**Parameters**

**port** (*str | int*) – Name of a COM port. Can be an integer or string.

**Returns**

Name of the COM port.

**Return type**

str

## Module contents

### Submodules

#### PodApi.Devices.BasicPodProtocol module

**class** PodApi.Devices.BasicPodProtocol.**Pod**(port: str | int, baudrate: int = 9600)

Bases: object

POD\_Basics handles basic communication with a generic POD device, including reading and writing packets and packet interpretation.

**\_port**

Instance-level COM\_io object, which handles the COM port

**Type**

COM\_io

**\_commands**

Instance-level POD\_Commands object, which stores information about the commands available to this POD device.

**Type**

POD\_Commands

**static ChoosePort**(forbidden: list[str] = []) → str

Systems checks user's Operating System, and chooses ports accordingly.

**Parameters**

**forbidden** (list[str], optional) – List of port names that are already used. Defaults to [].

**Returns**

String name of the port.

**Return type**

str

**FlushPort**() → bool

Reset the input and output serial port buffer.

**Returns**

True of the buffers are flushed, False otherwise.

**Return type**

bool

**GetDeviceCommands**() → dict[int, list[str | tuple[int] | bool]]

Gets the dictionary containing the class instance's available POD commands.

**Returns**

Dictionary containing the available commands and their information. Formatted as key(command number) : value([command name, number of argument ASCII bytes, number of return bytes, binary flag ])

**Return type**

dict[int, list[str|tuple[int]]|bool]]

**GetPODpacket**(*cmd: str | int, payload: int | bytes | tuple[int | bytes] | None = None*) → bytes

Builds a POD packet and writes it to a POD device via COM port. If an integer payload is give, the method will convert it into a bytes string of the length expected by the command. If a bytes payload is given, it must be the correct length.

**Parameters**

- **cmd** (*str | int*) – Command number.
- **payload** (*int | bytes | tuple[int | bytes], optional*) – None when there is no payload. If there is a payload, set to an integer value, bytes string, or tuple. Defaults to None.

**Raises**

- **Exception** – POD command does not exist.
- **Exception** – POD command requires a payload.

**Returns**

Bytes string of the POD packet.

**Return type**

bytes

**static GetU**(*u: int*) → int

number of hexadecimal characters for an unsigned u-bit value.

**Parameters**

**u** (*int*) – 8, 16, or 32 bits. Enter any other number for NOVALUE.

**Returns**

number of hexadecimal characters for an unsigned u-bit value.

**Return type**

int

**ReadPODpacket**(*validateChecksum: bool = True, timeout\_sec: int | float = 5*) → *Packet*

Reads a complete POD packet, either in standard or binary format, beginning with STX and ending with ETX. Reads first STX and then starts recursion.

**Parameters**

- **validateChecksum** (*bool, optional*) – Set to True to validate the checksum. Set to False to skip validation. Defaults to True.
- **timeout\_sec** (*int | float, optional*) – Time in seconds to wait for serial data. Defaults to 5.

**Returns**

POD packet beginning with STX and ending with ETX. This may be a standard packet, binary packet, or an unformatted packet (STX+something+ETX).

**Return type**

*Packet*

**SetBaudrateOfDevice**(*baudrate: int*) → bool

If the port is open, it will change the baud rate to the parameter's value.

**Parameters**

**baudrate** (*int*) – Baud rate to set for the open serial port.

**Returns**

True if successful at setting the baud rate, false otherwise.

**Return type**

bool

**WritePacket**(*cmd: str | int, payload: int | bytes | tuple[int | bytes] | None = None*) → *PacketStandard*

Builds a POD packet and writes it to the POD device.

**Parameters**

- **cmd** (*str | int*) – Command number.
- **payload** (*int | bytes | tuple[int | bytes], optional*) – None when there is no payload. If there is a payload, set to an integer value, bytes string, or tuple. Defaults to None.

**Returns**

Packet that was written to the POD device.

**Return type**

Packet\_Standard

**WriteRead**(*cmd: str | int, payload: int | bytes | tuple[int | bytes] | None = None, validateChecksum: bool = True*) → bytes

Writes a command with optional payload to POD device, then reads (once) the device response.

**Parameters**

- **cmd** (*str | int*) – Command number.
- **payload** (*int | bytes | tuple[int|bytes], optional*) – None when there is no payload. If there is a payload, set to an integer value or a bytes string. Defaults to None.
- **validateChecksum** (*bool, optional*) – Set to True to validate the checksum. Set to False to skip validation. Defaults to True.

**Returns**

Bytes string containing a POD packet beginning with STX and ending with ETX. This may be a standard packet, binary packet, or an unformatted packet (STX+something+ETX).

**Return type**

bytes

**\_ReadPODpacket\_Recursive**(*validateChecksum: bool = True*) → *Packet*

Reads the command number. If the command number ends in ETX, the packet is returned. Next, it checks if the command is allowed. Then, it checks if the command is standard or binary and reads accordingly, then returns the packet.

**Parameters****validateChecksum** (*bool, optional*) – Set to True to validate the checksum. Set to False to skip validation. Defaults to True.**Raises****Exception** – Cannot read an invalid command.**Returns**

POD packet beginning with STX and ending with ETX. This may be a standard packet, binary packet, or an unformatted packet (STX+something+ETX).

**Return type***Packet*|Packet\_Standard|Packet\_BinaryStandard

**`_Read_Binary`**(*prePacket: bytes, validateChecksum: bool = True*) → *PacketBinary*

Reads the remaining part of the variable-length binary packet. It first reads the standard packet (*prePacket*+payload+checksum+ETX). Then it determines how long the binary packet is from the payload of the standard POD packet and reads that many bytes. It then reads to ETX to get the checksum+ETX.

**Parameters**

- **`prePacket`** (*bytes*) – Bytes string containing the beginning of a POD packet: STX (1 byte) + command number (4 bytes)
- **`validateChecksum`** (*bool, optional*) – Set to True to validate the checksum. Set to False to skip validation. Defaults to True.

**Raises**

**Exception** – An exception is raised if the checksum is invalid (only if *validateChecksum=True*).

**Returns**

Variable-length binary POD packet.

**Return type**

`Packet_BinaryStandard`

**`_Read_GetCommand`**(*validateChecksum: bool = True*) → bytes

Reads one byte at a time up to 4 bytes to get the ASCII-encoded bytes command number. For each byte read, it can (1) start the recursion over if an STX is found, (2) returns if ETX is found, or (3) continue building the command number.

**Parameters**

**`validateChecksum`** (*bool, optional*) – Set to True to validate the checksum. Set to False to skip validation. Defaults to True.

**Returns**

4 byte long string containing the ASCII-encoded command number.

**Return type**

bytes

**`_Read_Standard`**(*prePacket: bytes, validateChecksum: bool = True*) → *PacketStandard*

Reads the payload, checksum, and ETX. Then it builds the complete standard POD packet in bytes.

**Parameters**

- **`prePacket`** (*bytes*) – Bytes string containing the beginning of a POD packet: STX (1 byte) + command number (4 bytes).
- **`validateChecksum`** (*bool, optional*) – Set to True to validate the checksum. Set to False to skip validation. Defaults to True.

**Raises**

**Exception** – An exception is raised if the checksum is invalid (only if *validateChecksum=True*).

**Returns**

Complete standard POD packet.

**Return type**

`Packet_Standard`

**`_Read_ToETX`**(*validateChecksum: bool = True*) → bytes

Reads one byte at a time until an ETX is found. It will restart the recursive read if an STX is found anywhere.



**Parameters**

**validateChecksum** (*bool*, *optional*) – Set to True to validate the checksum. Set to False to skip validation. Defaults to True.

**Returns**

Bytes string ending with ETX.

**Return type**

bytes

**static** **\_ValidateChecksum** (*msg: bytes*) → bool

Validates the checksum of a given POD packet. The checksum is valid if the calculated checksum from the data matches the checksum written in the packet.

**Parameters**

**msg** (*bytes*) – Bytes message containing a POD packet: STX (1 bytes) + data (? bytes) + checksum (2 bytes) + ETX (1 byte).

**Returns**

True if the checksum is correct, false otherwise.

**Return type**

bool

**Raises**

**Exception** – msg does not begin with STX or end with ETX.

**PodApi.Devices.PodDevice\_8206HR module**

**class** PodApi.Devices.PodDevice\_8206HR.**Pod8206HR** (*port: str | int, preampGain: int, baudrate: int = 9600*)

Bases: *Pod*

POD\_8206HR handles communication using an 8206HR POD device.

**\_preampGain**

Instance-level integer (10 or 100) preamplifier gain.

**Type**

int

**ReadPODpacket** (*validateChecksum: bool = True, timeout\_sec: int | float = 5*) → *Packet*

Reads a complete POD packet, either in standard or binary format, beginning with STX and ending with ETX. Reads first STX and then starts recursion.

**Parameters**

- **validateChecksum** (*bool*, *optional*) – Set to True to validate the checksum. Set to False to skip validation. Defaults to True.
- **timeout\_sec** (*int | float*, *optional*) – Time in seconds to wait for serial data. Defaults to 5.

**Returns**

POD packet beginning with STX and ending with ETX. This may be a standard packet, binary packet, or an unformatted packet (STX+something+ETX).

**Return type**

*Packet*

**`_Read_Binary`**(*prePacket: bytes, validateChecksum: bool = True*) → *PacketBinary4*

After receiving the prePacket, it reads the 8 bytes(TTL+channels) and then reads to ETX (checksum+ETX).

**Parameters**

- **prePacket** (*bytes*) – Bytes string containing the beginning of a POD packet: STX (1 byte) + command number (4 bytes).
- **validateChecksum** (*bool, optional*) – Set to True to validate the checksum. Set to False to skip validation. Defaults to True.

**Raises**

**Exception** – Bad checksum for binary POD packet read.

**Returns**

Binary4 POD packet.

**Return type**

Packet\_Binary4

**`static _TranslateTTLbyte_ASCII`**(*tTlByte: bytes*) → dict[str, int]

Separates the bits of each TTL (0-3) from a ASCII encoded byte.

**Parameters**

**tTlByte** (*bytes*) – One byte string for the TTL (ASCII encoded).

**Returns**

Dictionary of the TTLs. Values are 1 when input, 0 when output.

**Return type**

dict[str,int]

## PodApi.Devices.PodDevice\_8229 module

**`class PodApi.Devices.PodDevice_8229.Pod8229`**(*port: str | int, baudrate: int = 19200*)

Bases: *Pod*

POD\_8229 handles communication using an 8229 POD device.

**`static BuildSetDayScheduleArgument`**(*day: str | int, hours: list | tuple[bool | int], speed: int | list | tuple[int]*) → tuple[int]

Appends the day of the week code to the front of the encoded hourly schedule. this tuple is formatted to be used as the #141 'SET DAY SCHEDULE' argument.

**Parameters**

- **day** (*str | int*) – Day of the week. Can be either the name of the day (i.e. Sunday, Monday, etc.) or the 0-6 day code (0 for Sunday increacing to 6 for Saturday).
- **hours** (*list | tuple[bool | int]*) – Array of 24 items. The value is 1 for motor on and 0 for motor off.
- **speed** (*int | list | tuple[int]*) – Speed of the motor (0-100). This is an integer of all hours are the same speed. If there are multiple speeds, this should be an array of 24 items.

**Returns**

`_description_`

**Return type**

tuple[int]

**static CodeDayOfWeek**(*day: str*) → int

Converts the day of the week to an integer code understandable by the POD device. The day is determined by the first 1-2 characters of the string, which supports multiple abbreviations for days of the week.

**Parameters**

**day** (*str*) – Day of the week.

**Raises**

**Exception** – Invalid day of the week.

**Returns**

Code for the day of the week. Values are 0-6, with 0 for Saturday, 1 for Monday, ..., and 6 for Saturday.

**Return type**

int

**static CodeDaySchedule**(*hours: list | tuple[bool | int], speed: int | list | tuple[int]*) → list[int]

Bitmasks the day schedule to encode the motor on/off status and the motor speed. Use this for getting the command #141 'SET DAY SCHEDULE' U8x24 argument component.

**Parameters**

- **hours** (*list | tuple[bool | int]*) – Array of 24 items. The value is 1 for motor on and 0 for motor off.
- **speed** (*int | list | tuple[int]*) – Speed of the motor (0-100). This is an integer of all hours are the same speed. If there are multiple speeds, this should be an array of 24 items.

**Returns**

List of 24 integer items. The msb is the motor on/off flag and the remaining 7 bits are the speed.

**Return type**

list[int]

**static DecodeDayAndSchedule**(*dayschedule: bytes*)

**static DecodeDayOfWeek**(*day: int*) → str

Converts the integer code for a day of the week to a human-readable string.

**Parameters**

**day** (*int*) – Day of the week code must be 0-6.

**Returns**

Day of the week ('Sunday', 'Monday', etc.).

**Return type**

str

**static DecodeDaySchedule**(*schedule: bytes*) → dict[str, int | tuple[int]]

Interprets the return bytes from the command #142 'GET DAY SCHEDULE'.

**Parameters**

**schedule** (*bytes*) – 24 byte long bitstring with one U8 per hour in a day.

**Returns**

Dictionary with 'Hour' as a tuple of 24 0/1 values (0 is motor off and 1 is motor on) and 'Speed' as the motor speed (0-100). If the motor speed is the same every hour, 'Speed' will be an integer; otherwise, 'Speed' will be a tuple of 24 items.

**Return type**

dict[str,int|tuple[int]]

**static DecodeLCDSchedule**(*schedule: bytes*) → dict[str, str | list[int]]

Interprets the return bytes from the command #202 'LCD SET DAY SCHEDULE'.

**Parameters****schedule** (*bytes*) – 4 Byte long bitstring. Byte 3 is weekday, Byte 2 is hours 0-7, Byte 1 is hours 8-15, and byte 0 is hours 16-23.**Returns**

Dictionary with Day as the day of the week, and Hours containing a list of 24 0/1 values (one for each hour). Each bit represents the motor state in that hour, 1 for on and 0 for off.

**Return type**

dict[str,int|list[int]]

**static GetCurrentTime**() → tuple[int]

Gets a tuple to use as the argument for command #140 SET TIME containing values for the current time.

**Returns**

Tuple of 7 integer values. The format is (Seconds, Minutes, Hours, Day, Month, Year (without century, so 23 for 2023), Weekday (0 for Sunday))

**Return type**

tuple[int]

**ReadPODpacket**(*validateChecksum: bool = True, timeout\_sec: int | float = 5*) → *PacketStandard*

Reads a complete POD packet, either in standard or binary format, beginning with STX and ending with ETX. Reads first STX and then starts recursion.

**Parameters**

- **validateChecksum** (*bool, optional*) – Set to True to validate the checksum. Set to False to skip validation. Defaults to True.
- **timeout\_sec** (*int | float, optional*) – Time in seconds to wait for serial data. Defaults to 5.

**Returns**

POD packet beginning with STX and ending with ETX. This may be a standard packet, binary packet, or an unformatted packet (STX+something+ETX).

**Return type***Packet***WritePacket**(*cmd: str | int, payload: int | bytes | tuple[int | bytes] | None = None*) → bytes

Builds a POD packet and writes it to the POD device.

**Parameters**

- **cmd** (*str | int*) – Command number.
- **payload** (*int | bytes | tuple[int | bytes], optional*) – None when there is no payload. If there is a payload, set to an integer value, bytes string, or tuple. Defaults to None.

**Returns**

Bytes string that was written to the POD device.

**Return type**

bytes

**static** `_CodeDecimalAsHex(val: int) → int`

Builds an integer that equals the val argument when converted into hexadecimal. All integers are converted to hexadecimal ASCII encoded bytes. Some commands (i.e. 8229 #140) need decimal ASCII encoded bytes; to do this, give the return value of `_CodeDecimalAsHex()` as the payload. Example: I want a number that is equal to 16 in hex.  $1*16^1 + 6*16^0 = 22$ . Calling `_CodeDecimalAsHex(16)` will return 22.

**Parameters**

**val** (*int*) – Unsigned integer number.

**Returns**

integer that equals the val argument when converted into hexadecimal.

**Return type**

int

**static** `_Custom140SETTIME(payload: tuple[int]) → tuple[int]`

Custom function to translate the payload for command #140 SET TIME.

**Parameters**

**payload** (*tuple[int]*) – Default translated payload.

**Returns**

Tuple of times.

**Return type**

tuple[int]

**static** `_DecodeDecimalAsHex(val: int) → int`

Interprets an integer that was converted to a hexadecimal representation of a decimal value. In other words, this method reverses `_CodeDecimalAsHex()`.

**Parameters**

**val** (*int*) – Unsigned integer that was converted to a hexadecimal representation of a decimal value.

**Returns**

Unsigned integer as a true decimal number.

**Return type**

int

**static** `_Validate_Day(day: str | int) → int`

Raises an exception if the day is incorrectly formatted. If the day is given as a string, it will be converted to its integer code.

**Parameters**

**day** (*str | int*) – String day of the week or its respective integer code.

**Raises**

- **Exception** – The day integer argument must be 0-6.
- **Exception** – The day argument must be a str or int.

**Returns**

Integer code representing a day of the week.

**Return type**

int

**static** `_Validate_Hours(hours: list | tuple[bool | int]) → list[bool | int]`

Raises an exception if the hours is incorrectly formatted. Converts the hours into a list before returning it.

**Parameters**

**hours** (*list* | *tuple[bool | int]*) – Array with 24 items with values of 1/0 representing each hour

**Raises**

- **Exception** – The hours argument must be a list or tuple.
- **Exception** – The hours argument must have exactly 24 items.
- **Exception** – The hours items must be 0 or 1.

**Returns**

List with 24 items for each hour. The values are 1/0.

**Return type**

list[bool|int]

**static \_Validate\_Schedule**(*schedule: bytes, size: int*) → bytes

Raises an exception if the schedule is incorrectly formatted

**Parameters**

- **schedule** (*bytes*) – Bytes string containing the day schedule.
- **size** (*int*) – Number of U8 bytes.

**Raises**

- **Exception** – The schedule must be bytes.
- **Exception** – The schedule is the incorrect size

**Returns**

Same as the schedule argument.

**Return type**

bytes

**static \_Validate\_Speed**(*speed: int | list | tuple[int]*) → list[int]

Raises an exception if the speed is incorrectly formatted. If an integer speed is given, it will convert it to a list.

**Parameters**

**speed** (*int | list | tuple[int]*) – Motor speed (0-100). This can either be an integer or a tuple/list of 24 speeds.

**Raises**

- **Exception** – The speed argument must be an int, list, or tuple.
- **Exception** – The speed must be between 0 and 100.
- **Exception** – The speed argument must have exactly 24 items as a list/tuple.
- **Exception** – The speed must be between 0 and 100 for every list/tuple item.

**Returns**

List containing 24 motor speeds.

**Return type**

list[int]

**PodApi.Devices.PodDevice\_8401HR module**

```
class PodApi.Devices.PodDevice_8401HR.Pod8401HR(port: str | int, ssGain: tuple | list | dict[str, int | None]
      = {'A': None, 'B': None, 'C': None, 'D': None},
      preampGain: tuple | list | dict[str, int | None] = {'A':
      None, 'B': None, 'C': None, 'D': None}, baudrate: int
      = 9600)
```

Bases: *Pod*

POD\_8401HR handles communication using an 8401-HR POD device.

**\_ssGain**

Instance-level dictionary storing the second-stage gain for all four channels.

**Type**

dict[str,int|None]

**\_preampGain**

Instance-level dictionary storing the pramplifier gain for all four channels.

**Type**

dict[str,int|None]

**static CalculateBiasDAC\_GetDACValue**(vout: int | float) → int

Calculates the DAC value given the output voltage. Used for ‘GET/SET BIAS’ commands.

**Parameters**

**vout** (int | float) – Output voltage (+/- 2.048 V).

**Returns**

Integer of the DAC value (16 bit 2’s complement).

**Return type**

int

**static CalculateBiasDAC\_GetVout**(value: int) → float

Calculates the output voltage given the DAC value. Used for ‘GET/SET BIAS’ commands.

**Parameters**

**value** (int) – DAC value (16 bit 2’s complement).

**Returns**

Float of the output bias voltage [V].

**Return type**

float

**static DecodeChannelBitmask**(channels: bytes) → dict[str, int]

Converts the channel bitmask byte to a dictionary with each channel value. Use for ‘GET INPUT GROUND’ command payloads.

**Parameters**

**channels** (bytes) – U8 byte containing the channel configuration.

**Returns**

Dictionary with the channels as keys and values as the state. 0=Grounded and 1=Connected to Preamp.

**Return type**

dict[str,int]

**static DecodeSSConfigBitmask**(*config: bytes*)

Converts the SS configuration byte to a dictionary with the high-pass and gain. Use for 'GET SS CONFIG' command payloads.

**Parameters**

**config** (*bytes*) – U8 byte containing the SS configuration. Bit 0 = 0 for 0.5Hz Highpass, 1 for DC Highpass. Bit 1 = 0 for 5x gain, 1 for 1x gain.

**static DecodeTTLByte**(*tTlByte: bytes*) → dict[str, int]

Converts the TTL bytes argument into a dictionary of integer TTL values.

**Parameters**

**tTlByte** (*bytes*) – U8 byte containing the TTL bitmask.

**Returns**

Dictionary with TTL name keys and integer TTL values.

**Return type**

dict[str,int]

**static DecodeTTLPayload**(*payload: bytes*) → tuple[dict[str, int]]

Decodes a payload with the two TTL bytes.

**Parameters**

**payload** (*bytes*) – Bytes string of the POD packet payload.

**Returns**

Tuple with two TTL dictionaries.

**Return type**

tuple[dict[str, int]]

**static GetChannelBitmask**(*a: bool, b: bool, c: bool, d: bool*) → int

Gets a bitmask, represented by an unsigned integer, used for 'SET INPUT GROUND' command.

**Parameters**

- **a** (*bool*) – State for channel A, 0=Grounded and 1=Connected to Preamp.
- **b** (*bool*) – State for channel B, 0=Grounded and 1=Connected to Preamp.
- **c** (*bool*) – State for channel C, 0=Grounded and 1=Connected to Preamp.
- **d** (*bool*) – State for channel D, 0=Grounded and 1=Connected to Preamp.

**Returns**

Integer representing a bitmask.

**Return type**

int

**static GetChannelMapForPreampDevice**(*preampName: str*) → dict[str, str] | None

Get the channel mapping (channel labels for A,B,C,D) for a given device.

**Parameters**

**preampName** (*str*) – String for the device/sensor name.

**Returns**

Dictionary with keys A,B,C,D with values of the channel names. Returns None if the device name does not exist.

**Return type**

dict[str,str]|None



**static GetSSConfigBitmask**(*gain: int, highpass: float*) → int

Gets a bitmask, represented by an unsigned integer, used for ‘SET SS CONFIG’ command.

**Parameters**

- **gain** (*int*) – 1 for 1x gain. else for 5x gain.
- **highpass** (*float*) – 0 for DC highpass, else for 0.5Hz highpass.

**Returns**

Integer representing a bitmask.

**Return type**

int

**static GetSupportedPreampDevices**() → list[str]

Gets a list of device/sensor names used for channel mapping.

**Returns**

List of string names of all supported sensors.

**Return type**

list[str]

**static GetTTLbitmask**(*ext0: bool = 0, ext1: bool = 0, ttl4: bool = 0, ttl3: bool = 0, ttl2: bool = 0, ttl1: bool = 0*) → int

Builds an integer, which represents a binary mask, that can be used for TTL command arguments.

**Parameters**

- **ext0** (*bool, optional*) – boolean bit for ext0. Defaults to 0.
- **ext1** (*bool, optional*) – boolean bit for ext1. Defaults to 0.
- **ttl4** (*bool, optional*) – boolean bit for ttl4. Defaults to 0.
- **ttl3** (*bool, optional*) – boolean bit for ttl3. Defaults to 0.
- **ttl2** (*bool, optional*) – boolean bit for ttl2. Defaults to 0.
- **ttl1** (*bool, optional*) – boolean bit for ttl1. Defaults to 0.

**Returns**

Integer number to be used as a bit mask.

**Return type**

int

**static IsPreampDeviceSupported**(*name: str*) → bool

Checks if the argument exists in channel map for all preamp sensors.

**Parameters**

**name** (*str*) – name of the device

**Returns**

True if the name exists in \_\_CHANNELMAPALL, false otherwise.

**Return type**

bool

**ReadPODpacket**(*validateChecksum: bool = True, timeout\_sec: int | float = 5*) → *Packet*

Reads a complete POD packet, either in standard or binary format, beginning with STX and ending with ETX. Reads first STX and then starts recursion.

**Parameters**

- **validateChecksum** (*bool, optional*) – Set to True to validate the checksum. Set to False to skip validation. Defaults to True.
- **timeout\_sec** (*int/float, optional*) – Time in seconds to wait for serial data. Defaults to 5.

**Returns**

POD packet beginning with STX and ending with ETX. This may be a standard packet, binary packet, or an unformatted packet (STX+something+ETX).

**Return type**

*Packet*

**WritePacket** (*cmd: str | int, payload: int | bytes | tuple[int | bytes] | None = None*) → *PacketStandard*

Builds a POD packet and writes it to the POD device.

**Parameters**

- **cmd** (*str | int*) – Command number.
- **payload** (*int | bytes | tuple[int | bytes], optional*) – None when there is no payload. If there is a payload, set to an integer value, bytes string, or tuple. Defaults to None.

**Returns**

Packet that was written to the POD device.

**Return type**

*Packet\_Standard*

**static \_FixABCDtype** (*info: tuple | list | dict, thisIs: str = ""*) → dict

Converts the info argument into a dictionary with A, B, C, and D as keys.

**Parameters**

- **info** (*tuple | list | dict*) – Variable to be converted into a dictionary.
- **thisIs** (*str, optional*) – Description of the info argument, which is used in Exception statements. Defaults to ‘’.

**Raises**

- **Exception** – The dictionary has improper keys; keys must be ['A','B','C','D'].
- **Exception** – The argument must have only four values.
- **Exception** – The argument must be a tuple, list, or dict.

**Returns**

The info argument converted to a dictionary with A, B, C, and D as keys.

**Return type**

dict

**\_Read\_Binary** (*prePacket: bytes, validateChecksum: bool = True*)

After receiving the prePacket, it reads the 23 bytes (binary data) and then reads to ETX.

**Parameters**

- **prePacket** (*bytes*) – Bytes string containing the beginning of a POD packet: STX (1 byte) + command number (4 bytes).
- **validateChecksum** (*bool, optional*) – Set to True to validate the checksum. Set to False to skip validation. Defaults to True.

**Raises**

**Exception** – Bad checksum for binary POD packet read.

**static** `_ValidatePreampGain(preampGain: dict) → None`

Checks that the preamplifier gain dictionary has proper values (10, 100, or None).

**Parameters**

**preampGain** (*dict*) – preamplifier gain dictionary.

**Raises**

**Exception** – EEG/EMG preampGain must be 10 or 100. For biosensors, the preampGain is None.

**static** `_ValidateSsGain(ssgain: dict)`

Checks that the second stage gain dictionary has proper values (1, 5, or None).

**Parameters**

**ssgain** (*dict*) – Second stage gain dictionary.

**Raises**

**Exception** – The ssGain must be 1 or 5; set ssGain to None if no-connect.

```
__CHANNELMAPALL: dict[str, dict[str, str]] = {'8406-2BIO': {'A': 'Bio1', 'B': 'Bio2', 'C': 'NC', 'D': 'NC'}, '8406-BIO': {'A': 'Bio', 'B': 'NC', 'C': 'NC', 'D': 'NC'}, '8406-EEG2BIO': {'A': 'Bio1', 'B': 'EEG1', 'C': 'EMG', 'D': 'Bio2'}, '8406-SE': {'A': 'Bio', 'B': 'EEG1', 'C': 'EMG', 'D': 'EEG2'}, '8406-SE3': {'A': 'Bio', 'B': 'EEG1', 'C': 'EEG3', 'D': 'EEG2'}, '8406-SE31M': {'A': 'EMG', 'B': 'EEG1', 'C': 'EEG3', 'D': 'EEG2'}, '8406-SE4': {'A': 'EEG4', 'B': 'EEG1', 'C': 'EEG3', 'D': 'EEG2'}, '8406-SL': {'A': 'Bio', 'B': 'EEG1', 'C': 'EMG', 'D': 'EEG2'}, '8407-SE': {'A': 'Bio', 'B': 'EEG1', 'C': 'EMG', 'D': 'EEG2'}, '8407-SE-2BIO': {'A': 'Bio1', 'B': 'Bio2', 'C': 'EMG', 'D': 'EEG2'}, '8407-SE3': {'A': 'Bio', 'B': 'EEG1', 'C': 'EEG3', 'D': 'EEG2'}, '8407-SE31M': {'A': 'EEG3', 'B': 'EEG1', 'C': 'EMG', 'D': 'EEG2'}, '8407-SE4': {'A': 'EEG4', 'B': 'EEG1', 'C': 'EEG3', 'D': 'EEG2'}, '8407-SL': {'A': 'Bio', 'B': 'EEG1', 'C': 'EMG', 'D': 'EEG2'}, '8407-SL-2BIO': {'A': 'Bio1', 'B': 'Bio2', 'C': 'EMG', 'D': 'EEG2'}}
```

Class-level dictionary containing the channel map for all preamplifier devices.

**PodApi.Devices.PodDevice\_8480SC module**

**class** `PodApi.Devices.PodDevice_8480SC.Pod8480SC(port: str | int, baudrate: int = 9600)`

Bases: `Pod`

POD\_8480SC handles communication using an 8480-SC POD device.

**static** `DecodeStimulusConfigBits(config: int) → dict`

Converts an integer into 3 values, representing 3 individual bits of the Stimulus Config Bits.

**Parameters**

**config** (*int*) – an integer is passed in, and it represents the Config Flag byte.

**Returns**

Keys as the names of the bits, the values representing values at each bit.

**Return type**

dict

**static DecodeSyncConfigBits**(*config: int*) → dict

Converts an integer into 3 values, representing 3 individual bits of the Sync Config Bits.

**Parameters**

**config** (*int*) – an integer is passed in, and it represents the Sync Config Flag byte.

**Returns**

Keys as the names of the bits, the values representing values at each bit.

**Return type**

dict

**static DecodeTTLConfigBits**(*config: int*) → dict

Converts an integer into 3 values, representing 3 individual bits of the TTL Config Bits.

**Parameters**

**config** (*int*) – an integer is passed in, and it represents the TTL Setup Config Flag Byte.

**Returns**

Keys as the names of the bits, the values representing values at each bit.

**Return type**

dict

**ReadPODpacket**(*validateChecksum: bool = True, timeout\_sec: int | float = 5*) → *PacketStandard*

Reads a complete POD packet, either in standard or binary format, beginning with STX and ending with ETX. Reads first STX and then starts recursion.

**Parameters**

- **validateChecksum** (*bool, optional*) – Set to True to validate the checksum. Set to False to skip validation. Defaults to True.
- **timeout\_sec** (*int | float, optional*) – Time in seconds to wait for serial data. Defaults to 5.

**Returns**

POD packet beginning with STX and ending with ETX. This may be a standard packet, binary packet, or an unformatted packet (STX+something+ETX).

**Return type**

*Packet*

**static StimulusConfigBits**(*optoElec: bool, monoBiphasic: bool, Simul: bool*) → int

Incoming inputs are bitmasked into an integer value. This value is later given as part of a payload to command #102 ‘SET STIMULUS’.

**Parameters**

- **optoElec** (*bool*) – Bit is Opto/Electrical.
- **monoBiphasic** (*bool*) – Bit 1 is Monophasic/Biphasic.
- **Simul** (*bool*) – Bit 2 is Simultaneous.

**Returns**

which represents the Config flag byte in the Stimulus Command. The return value is the seventh item in the payload for command ‘SET STIMULUS’.

**Return type**

int

**static SyncConfigBits**(*sync\_level: bool, sync\_idle: bool, signal\_trigger: bool*) → int

Incoming inputs are bitmasked into an integer value. This value is later given as the payload to command #127 'SET SYNC CONFIG'.

**Parameters**

- **sync\_level** (*bool*) – Bit 0 is Sync Level.
- **sync\_idle** (*bool*) – Bit 1 is Stimulus Triggering.
- **signal\_trigger** (*bool*) – Bit 2 is Signal/Trigger.

**Returns**

which represents the Sync Config Bits format value.

**Return type**

int

**static TtlConfigBits**(*trigger: bool, stimtrig: bool, input\_sync: bool*) → int

Incoming inputs are bitmasked into an integer value. This value is later given as part of the payload to command #109 'SET TTL SETUP'. This commands accepts 3 items in the payload, and the return value of this function is given as the second item.

**Parameters**

- **trigger** (*bool*) – Bit 0 is 0 for rising edge, 1 for falling edge.
- **stimtrig** (*bool*) – Bit 1 is 0 for TTL event notifications, 1 for TTL inputs as triggers.
- **input\_sync** (*bool*) – Bit 7 is 0 for normal TTL operation, 1 for TTL pin operates as a sync output.

**Returns**

which represents the TTL Config Bits Format value.

**Return type**

int

**WritePacket**(*cmd: str | int, payload: int | bytes | tuple[int | bytes] | None = None*) → [PacketStandard](#)

Builds a POD packet and writes it to the POD device.

**Parameters**

- **cmd** (*str | int*) – Command number.
- **payload** (*int | bytes | tuple[int | bytes], optional*) – None when there is no payload. If there is a payload, set to an integer value, bytes string, or tuple. Defaults to None.

**Returns**

Packet that was written to the POD device.

**Return type**

Packet\_Standard

**static \_Custom108GETTTLSETUP**(*payload: bytes*) → tuple[int | dict]

Custom function to translate the TTL setup for command #108 GET TTL SETUP.

**Parameters**

**payload** (*bytes*) – Bytes string of the POD packet payload.

**Returns**

Tuple of the TTL setup.

**Return type**

tuple[int|dict]

**static** `_Custom109SETTTLSETUP(payload: bytes) → tuple[int | dict]`

Custom function to translate the TTL setup for command #109 SET TTL SETUP.

**Parameters****payload** (bytes) – Bytes string of the POD packet payload.**Returns**

Tuple of the TTL setup.

**Return type**

tuple[int|dict]

**static** `_CustomSTIMULUS(payload: bytes, defaultPayload: tuple) → tuple``_summary_`**Parameters**

- **payload** (bytes) – Bytes string of the POD packet payload.
- **defaultPayload** (tuple) – Default translated payload.

**Returns**

Tuple of the translated stimulus payload.

**Return type**

tuple

**static** `_CustomSYNCONFIG(payload: bytes) → dict`

Custom function to translate the sync config.

**Parameters****payload** (bytes) – Bytes string of the POD packet payload.**Returns**

Keys as the names of the bits, the values representing values at each bit.

**Return type**

dict

## Module contents

### 1.1.1.3 PodApi.Packets package

#### Submodules

#### PodApi.Packets.Binary module

**class** `PodApi.Packets.Binary.PacketBinary(pkt: bytes, commands: CommandSet | None = None)`Bases: [Packet](#)

Container class that stores a standard binary command packet for a POD device. The format is STX (1 byte) + command number (4 bytes) + length of binary (4 bytes) + checksum (2 bytes) + ETX (1 bytes) + binary (LENGTH bytes) + checksum (2 bytes) + ETX (1 bytes)

**binaryLength**

Number of bytes of binary data from the packet.

**Type**

bytes

**binaryData**

Variable length binary data from the packet.

**Type**

bytes

**BinaryLength()** → int

Translate the binary ASCII encoding of the binary data length into a readable integer

**Returns**

Integer of the binary data length.

**Return type**

int

**static CheckIfPacketIsValid(msg: bytes)**

Raises an Exception if the packet is incorrectly formatted.

**Parameters**

**msg** (*bytes*) – Bytes string containing a POD packet. Should begin with STX and end with ETX.

**Raises**

- **Exception** – Packet is too small to be a standard packet.
- **Exception** – A standard binary packet must have an ETX before the binary bytes.

**GetBinaryData()** → bytes

Gets the binary data from a POD packet.

**Parameters**

**pkt** (*bytes*) – Bytes string containing a POD packet.

**Returns**

Bytes string containing binary data.

**Return type**

bytes

**static GetBinaryLength(pkt: bytes)** → bytes

Gets the length, or number of bytes, of the binary data in a POD packet.

**Parameters**

**pkt** (*bytes*) – Bytes string containing a POD packet.

**Returns**

Bytes string of the length of the binary data.

**Return type**

bytes

**static GetMinimumLength()** → int

Gets the number of bytes in the smallest possible packet; STX (1 byte) + something + ETX (1 byte).

**Returns**

integer representing the minimum length of a binary POD command packet. Format is STX (1 byte) + command number (4 bytes) + length of binary (4 bytes) + checksum (2 bytes) + ETX (1 bytes) + binary (LENGTH bytes) + checksum (2 bytes) + ETX (1 bytes)

**Return type**

int

**TranslateAll()** → dict[str, Any]

Builds a dictionary containing all parts of the POD packet in readable values.

**Returns**

Dictionary with the command number, binary packet length, and binary data.

**Return type**

dict[str,Any]

**UnpackAll()** → dict[str, bytes]

Builds a dictionary containing all parts of the POD packet in bytes.

**Returns**

Dictionary with the command number, binary packet length, and binary data.

**Return type**

dict[str,bytes]

## PodApi.Packets.Binary4 module

**class** PodApi.Packets.Binary4.**PacketBinary4**(*pkt: bytes, preampGain: int, commands: CommandSet | None = None*)

Bases: [Packet](#)

Container class that stores a binary4 command packet for a POD device. The format is STX (1 byte) + command (4 bytes) + packet number (1 byte) + TTL (1 byte) + CH0 (2 bytes) + CH1 (2 bytes) + CH2 (2 bytes) + checksum (2 bytes) + ETX (1 byte).

**preampGain**

Preamplifier gain. This should be 10 or 100 for an 8206-HR device.

**Type**

int

**packetNumber**

Packet number for this POD packet.

**Type**

bytes

**ttl**

TTL data for this packet.

**Type**

bytes

**ch0**

channel 0 data for this packet.

**Type**

bytes



**ch1**

channel 1 data for this packet.

**Type**

bytes

**ch2**

channel 2 data for this packet.

**Type**

bytes

**static BinaryBytesToVoltage**(*value: bytes, preampGain: int*) → float

Converts a binary bytes value read from POD device and converts it to the real voltage value at the preamplifier input.

**Parameters**

**value** (*bytes*) – Bytes string containing voltage measurement.

**Returns**

A number containing the voltage in Volts [V].

**Return type**

float

**Ch**(*n: int*) → float

Translates the binary channel n bytes into a voltage.

**Parameters**

**n** (*int*) – Channel number. Should be 0, 1, or 2.

**Raises**

**Exception** – Channel does not exist.

**Returns**

Voltage of channel n in Volts.

**Return type**

float

**static CheckIfPacketIsValid**(*msg: bytes*)

Raises an Exception if the packet is incorrectly formatted.

**Parameters**

**msg** (*bytes*) – Bytes string containing a POD packet. Should begin with STX and end with ETX.

**Raises**

**Exception** – Packet the wrong size to be a binary4 packet.

**static GetBinaryLength**() → int

Gets the number of bytes of binary data in a binary4 packet.

**Returns**

Integer representing the number of binary encoded bytes in a binary4 packet.

**Return type**

int

**static GetCh**(*n: int, pkt: bytes*) → bytes

Gets the channel n bytes from a POD packet.

**Parameters**

- **n** (*int*) – Channel number. Should be 0, 1, or 2.
- **pkt** (*bytes*) – Bytes string containing a POD packet. Should begin with STX and end with ETX.

**Returns**

Bytes string of the channel 0 data.

**Return type**

bytes

**Returns**

Channel does not exist.

**Return type**

bytes

**static GetMinimumLength()** → int

Gets the number of bytes in the smallest possible binary4 packet; STX (1 byte) + command (4 bytes) + packet number (1 byte) + TTL (1 byte) + CH0 (2 bytes) + CH1 (2 bytes) + CH2 (2 bytes) + checksum (2 bytes) + ETX (1 byte).

**Returns**

Integer representing the minimum length of a binary4 POD packet.

**Return type**

int

**static GetPacketNumber**(*pkt: bytes*) → bytes

Gets the packet number in bytes from a POD packet.

**Parameters**

**pkt** (*bytes*) – Bytes string containing a POD packet. Should begin with STX and end with ETX.

**Returns**

Bytes string of the packet number.

**Return type**

bytes

**static GetTTL**(*pkt: bytes*) → bytes

Gets the TTL bytes from a POD packet

**Parameters**

**pkt** (*bytes*) – Bytes string containing a POD packet. Should begin with STX and end with ETX.

**Returns**

Bytes string of the TTL data.

**Return type**

bytes

**PacketNumber**() → int

Translates the binary packet number into a readable integer.

**Returns**

Integer of the packet number.

**Return type**

int

**TranslateAll()** → dict[str, Any]

Builds a dictionary containing all parts of the POD packet in readable values.

**Returns**

Dictionary with the command number, packet number, TTL and channels 0, 1, and 2.

**Return type**

dict[str, Any]

**static TranslateBinaryTTLbyte(ttlByte: bytes)** → dict[str, int]

Separates the bits of each TTL (0-3) from a binary encoded byte.

**Parameters****ttlByte** (bytes) – One byte string for the TTL (binary encoded).**Returns**

Dictionary of the TTLs. Values are 1 when input, 0 when output.

**Return type**

dict[str,int]

**Ttl()** → dict[str, int]

Translates the binary TTL bytes into a dictionary containing each TTL value.

**Returns**

Dictionary with TTL name keys and TTL data as values.

**Return type**

dict[str,int]

**UnpackAll()** → dict[str, bytes]

Builds a dictionary containing all parts of the POD packet in bytes.

**Returns**

Dictionary with the command number, packet number, TTL and channels 0, 1, and 2.

**Return type**

dict[str,bytes]

**PodApi.Packets.Binary5 module**

**class PodApi.Packets.Binary5.PacketBinary5**(pkt: bytes, ssGain: dict[str, int | None] = {'A': None, 'B': None, 'C': None, 'D': None}, preampGain: dict[str, int | None] = {'A': None, 'B': None, 'C': None, 'D': None}, commands: CommandSet | None = None)

Bases: [Packet](#)

Container class that stores a binary5 command packet for a POD device. The format is STX (1 byte) + command (4 bytes) + packet number (1 byte) + status (1 byte) + channels (9 bytes) + AEXT0 (2 bytes) + AEXT1 (2 bytes) + ATTL1 (2 bytes) + ATTL2 (2 bytes) + ATTL3 (2 bytes) + ATTL4 (2 bytes) + checksum (2 bytes) + EXT (1 byte)

**\_ssGain**

Dictionary with A, B, C, D keys and second stage gain values (1, 5, or None).

**Type**

dict[str,int|None]

**\_preampGain**

Dictionary with A, B, C, D keys and preamplifier gain values (10, 100, or None).

**Type**

dict[str,int|None]

**packetNumber**

Packet number for this POD packet.

**Type**

bytes

**status**

Status for this POD packet.

**Type**

bytes

**channels**

channel A, B, C, and D data for this POD packet.

**Type**

bytes

**aEXT0**

Analog EXT0 data for this POD packet.

**Type**

bytes

**aEXT1**

Analog EXT1 data for this POD packet.

**Type**

bytes

**aTTL1**

Analog TTL1 data for this POD packet.

**Type**

bytes

**aTTL2**

Analog TTL2 data for this POD packet.

**Type**

bytes

**aTTL3**

Analog TTL3 data for this POD packet.

**Type**

bytes

**aTTL4**

Analog TTL4 data for this POD packet.

**Type**

bytes

**AnalogEXT**(*n: int*) → float

Translates the analog EXT value into a voltage.

**Parameters**

**n** (*int*) – Analog EXT number. Should be 0 or 1.

**Raises**

**Exception** – AEXT does not exist.

**Returns**

Analog EXT voltage in volts (V).

**Return type**

float

**AnalogTTL**(*n: int*) → float

Translates the analog TTL value into a voltage.

**Parameters**

**n** (*int*) – Analog TTL number. Should be 1, 2, 3, or 4.

**Raises**

**Exception** – ATTL does not exist.

**Returns**

Analog TTL voltage in volts (v).

**Return type**

float

**Channel**(*c: str*) → float

Translates the channel data into a voltage.

**Parameters**

**c** (*str*) – Channel character. Should be A, B, C, or D.

**Raises**

**Exception** – Channel does not exist.

**Returns**

Voltage of the channel in volts (V).

**Return type**

float

**static CheckIfPacketIsValid**(*msg: bytes*)

Raises an Exception if the packet is incorrectly formatted.

**Parameters**

**msg** (*bytes*) – Bytes string containing a POD packet. Should begin with STX and end with ETX.

**Raises**

**Exception** – Packet the wrong size to be a binary5 packet.

**static GetAnalogEXT**(*n: int, pkt: bytes*) → bytes

Gets the analog EXT from a POD packet.

**Parameters**

- **n** (*int*) – Analog EXT number. Should be 0 or 1.

- **pkt** (*bytes*) – Bytes string containing a POD packet. Should begin with STX and end with ETX.

**Raises**

**Exception** – AEXT does not exist.

**Returns**

Bytes string of the AEXT.

**Return type**

bytes

**static GetAnalogTTL**(*n: int, pkt: bytes*) → bytes

Gets the analog TTL from a POD packet.

**Parameters**

- **n** (*int*) – Analog TTL number. Should be 1, 2, 3, or 4.
- **pkt** (*bytes*) – Bytes string containing a POD packet. Should begin with STX and end with ETX.

**Raises**

**Exception** – ATTL does not exist.

**Returns**

Bytes string of the ATTL.

**Return type**

bytes

**static GetBinaryLength**() → int

Gets the number of bytes of binary data in a binary5 packet.

**Returns**

Integer representing the number of binary encoded bytes in a binary5 packet.

**Return type**

int

**static GetChannels**(*pkt: bytes*) → bytes

Gets the channel bytes for channels A, B, C, and D together from a POD packet.

**Parameters**

**pkt** (*bytes*) – Bytes string containing a POD packet. Should begin with STX and end with ETX.

**Returns**

Bytes string of the channels A, B, C, and D together.

**Return type**

bytes

**static GetMinimumLength**() → int

Gets the number of bytes in the smallest possible binary4 packet; STX (1 byte) + command (4 bytes) + packet number (1 byte) + status (1 byte) + channels (9 bytes) + AEXT0 (2 bytes) + AEXT1 (2 bytes) + ATTL1 (2 bytes) + ATTL2 (2 bytes) + ATTL3 (2 bytes) + ATTL4 (2 bytes) + checksum (2 bytes) + EXT (1 byte)

**Returns**

Integer representing the minimum length of a binary5 POD packet.

**Return type**

int

**static GetPacketNumber**(*pkt: bytes*) → bytes

Gets the packet number in bytes from a POD packet.

**Parameters****pkt** (*bytes*) – Bytes string containing a POD packet. Should begin with STX and end with ETX.**Returns**

Bytes string of the packet number.

**Return type**

bytes

**static GetStatus**(*pkt: bytes*) → bytes

Gets the status value in bytes from a POD packet.

**Parameters****pkt** (*bytes*) – Bytes string containing a POD packet. Should begin with STX and end with ETX.**Returns**

Bytes string of the status.

**Return type**

bytes

**PacketNumber**() → int

Translates the binary packet number into a readable integer.

**Returns**

Integer of the packet number.

**Return type**

int

**Status**() → int

Translates the binary status value into a readable integer

**Returns**

Integer status value.

**Return type**

int

**TranslateAll**() → dict[str, Any]

Builds a dictionary containing all parts of the POD packet in readable values.

**Returns**

Dictionary with the command number, packet number, status, channels, analog EXT, and analog TTL.

**Return type**

dict[str, bytes]

**UnpackAll**() → dict[str, bytes]

Builds a dictionary containing all parts of the POD packet in bytes.

**Returns**

Dictionary with the command number, packet number, status, channels, analog EXT, and analog TTL.

**Return type**

dict[str, bytes]

**static** **\_Voltage\_PrimaryChannels**(*value: int, ssGain: int | None = None, PreampGain: int | None = None*) → float

Converts a value to a voltage for a primary channel.

**Parameters**

- **value** (*int*) – Value to be converted to voltage.
- **ssGain** (*int | None, optional*) – Second stage gain. Defaults to None.
- **PreampGain** (*int | None, optional*) – Preamplifier gain. Defaults to None.

**Returns**

Number of the voltage in volts [V]. Returns value if no gain is given (no-connect).

**Return type**

float

**static** **\_Voltage\_PrimaryChannels\_Biosensor**(*value: int, ssGain: int*) → float

Converts a value to a voltage for a biosensor primary channel.

**Parameters**

- **value** (*int*) – Value to be converted to voltage.
- **ssGain** (*int*) – Second stage gain.

**Returns**

Number of the voltage in volts [V].

**Return type**

float

**static** **\_Voltage\_PrimaryChannels\_EEGEMG**(*value: int, ssGain: int, PreampGain: int*) → float

Converts a value to a voltage for an EEG/EMG primary channel.

**Parameters**

- **value** (*int*) – Value to be converted to voltage.
- **ssGain** (*int*) – Second stage gain.
- **PreampGain** (*int*) – Preamplifier gain.

**Returns**

Number of the voltage in volts [V].

**Return type**

float

**static** **\_Voltage\_SecondaryChannels**(*value: int*) → float

Converts a value to a voltage for a secondary channel.

**Parameters**

**value** (*int*) – Value to be converted to voltage.

**Returns**

Number of the voltage in volts [V].



**Return type**  
float

## PodApi.Packets.Packet module

**class** PodApi.Packets.Packet.**Packet**(*pkt: bytes, commands: CommandSet | None = None*)

Bases: object

Container class that stores a command packet for a POD device. The format is STX (1 byte) + command number (4 bytes) + data (? bytes) + ETX (1 byte). This class also collection of methods for creating and interpreting POD packets.

**\_commands**

Available commands for a POD device.

**Type**

POD\_Commands | None

**rawPacket**

Bytes string containing a POD packet. Should begin with STX and end with ETX.

**Type**

bytes

**commandNumber**

Command number from the Pod packet.

**Type**

bytes | None

**static** ASCIIbytesToInt\_Split(*msg: bytes, keepTopBits: int, cutBottomBits: int*) → int

Converts a specific bit range in an ASCII-encoded bytes object to an integer.

**Parameters**

- **msg** (*bytes*) – Bytes message holding binary information to be converted into an integer.
- **keepTopBits** (*int*) – Integer position of the msb of desired bit range.
- **cutBottomBits** (*int*) – Integer number of lsb to remove.

**Returns**

Integer result from the ASCII-encoded bytes message in a given bit range.

**Return type**

int

**static** AsciiBytesToInt(*msg\_b: bytes, signed: bool = False*) → int

Converts a ASCII-encoded bytes message into an integer. It does this using a base-16 conversion. If the message is signed and the msb is '1', the integer will be converted to it's negative 2's complement.

**Parameters**

- **msg\_b** (*bytes*) – Bytes message to be converted to an integer. The bytes must be base-16 or the conversion will fail.
- **signed** (*bool, optional*) – True if the message is signed, false if unsigned. Defaults to False.

**Returns**

Integer result from the ASCII-encoded byte conversion.

**Return type**

int

**static BinaryBytesToInt**(*msg: bytes, byteorder: str = 'big', signed: bool = False*) → int

Converts binary-encoded bytes into an integer.

**Parameters**

- **msg** (*bytes*) – Bytes message holding binary information to be converted into an integer.
- **byteorder** (*str, optional*) – Ordering of bytes. ‘big’ for big endian and ‘little’ for little endian. Defaults to ‘big’.
- **signed** (*bool, optional*) – Boolean flag to mark if the msg is signed (True) or unsigned (False). Defaults to False.

**Returns**

Integer result from the binary-encoded bytes message.

**Return type**

int

**static BinaryBytesToInt\_Split**(*msg: bytes, keepTopBits: int, cutBottomBits: int, byteorder: str = 'big', signed: bool = False*) → int

Converts a specific bit range in a binary-encoded bytes object to an integer.

**Parameters**

- **msg** (*bytes*) – Bytes message holding binary information to be converted into an integer.
- **keepTopBits** (*int*) – Integer position of the msb of desired bit range.
- **cutBottomBits** (*int*) – Integer number of lsb to remove.
- **byteorder** (*str, optional*) – Ordering of bytes. ‘big’ for big endian and ‘little’ for little endian. Defaults to ‘big’.
- **signed** (*bool, optional*) – Boolean flag to mark if the msg is signed (True) or unsigned (False). Defaults to False.

**Returns**

Integer result from the binary-encoded bytes message in a given bit range.

**Return type**

int

**static BuildPODpacket\_Standard**(*commandNumber: int, payload: bytes | None = None*) → bytes

Builds a standard POD packet as bytes: STX (1 byte) + command number (4 bytes) + optional packet (? bytes) + checksum (2 bytes)+ ETX (1 bytes).

**Parameters**

- **commandNumber** (*int*) – Integer representing the command number. This will be converted into a 4 byte long ASCII-encoded bytes string.
- **payload** (*bytes | None, optional*) – bytes string containing the payload. Defaults to None.

**Returns**

Bytes string of a complete standard POD packet.

**Return type**

bytes

**static CheckIfPacketIsValid**(*msg: bytes*)

Raises an Exception if the packet is incorrectly formatted.

**Parameters**

**msg** (*bytes*) – Bytes string containing a POD packet. Should begin with STX and end with ETX.

**Raises**

- **Exception** – Packet must begin with STX.
- **Exception** – Packet must end in ETX

**static Checksum**(*bytesIn: bytes*) → bytes

Calculates the checksum of a given bytes message. This is achieved by summing each byte in the message, inverting, and taking the last byte.

**Parameters**

**bytesIn** (*bytes*) – Bytes message containing POD packet data.

**Returns**

Two ASCII-encoded bytes containing the checksum for bytesIn.

**Return type**

bytes

**CommandNumber**() → int

Translate the binary ASCII encoding into a readable integer

**Returns**

Integer of the command number.

**Return type**

int

**static ETX**() → bytes

Get end-of-transmission (ETX) character in bytes. ETX marks the end byte of a POD Packet.

**Returns**

Bytes for ETX(0x03).

**Return type**

bytes

**static GetCommandNumber**(*pkt: bytes*) → bytes | None

Gets the command number bytes from a POD packet.

**Parameters**

**pkt** (*bytes*) – Bytes string containing a POD packet. Should begin with STX and end with ETX.

**Returns**

Bytes string of the command number, if available.

**Return type**

bytes|None

**static GetMinimumLength**() → int

Gets the number of bytes in the smallest possible packet; STX (1 byte) + something + ETX (1 byte).

**Returns**

integer representing the minimum length of a generic bytes string.

**Return type**

int

**HasCommandNumber()** → bool

Checks if the packet has a command number.

**Returns**

True if the packet has a command number, False otherwise.

**Return type**

bool

**HasCommands()** → bool

Checks if the Packet instance has commands set.

**Returns**

True if the commands have been set, false otherwise.

**Return type**

bool

**static IntToAsciiBytes**(*value: int, numChars: int*) → bytes

Converts an integer value into ASCII-encoded bytes.

First, it converts the integer value into a usable uppercase hexadecimal string. Then it converts the ASCII code for each character into bytes. Lastly, it ensures that the final message is the desired length.

Example: if value=2 and numBytes=4, the returned ASCII will show b'0002', which is '0x30 0x30 0x30 0x32' in bytes. Uses the 2's complement if the val is negative.

**Parameters**

- **value** (*int*) – Integer value to be converted into ASCII-encoded bytes.
- **numChars** (*int*) – Number characters to be the length of the ASCII-encoded message.

**Returns**

Bytes that are ASCII-encoded conversions of the value parameter.

**Return type**

bytes

**static PayloadToBytes**(*payload: int | bytes | tuple[int | bytes], argSizes: tuple[int]*) → bytes

Converts a payload into a bytes string.

**Parameters**

- **payload** (*int | bytes | tuple[int | bytes]*) – Integer, bytes, or tuple containing the payload.
- **argSizes** (*tuple[int]*) – Tuple of the argument sizes.

**Raises**

- **Exception** – Payload requires multiple arguments, use a tuple.
- **Exception** – Payload is the wrong size.
- **Exception** – Payload has an incorrect number of items.
- **Exception** – Payload has invalid values.
- **Exception** – Payload is an invalid type.

**Returns**

Bytes string of the payload.

**Return type**

bytes

**static STX()** → bytes

Get start-of-transmission (STX) character in bytes. STX marks the starting byte of a POD Packet.

**Returns**

Bytes for STX (0x02).

**Return type**

bytes

**TranslateAll()** → dict[str, Any]

Builds a dictionary containing all parts of the POD packet in readable values.

**Raises****Exception** – Nothing to translate.**Returns**

Dictionary with the command number.

**Return type**

dict[str,Any]

**static TwosComplement**(*val: int, nbits: int*) → int

Gets the 2's complement of the argument value.

**Parameters**

- **val** (*int*) – Value to be complemented.
- **nbits** (*int*) – Number of bits in the value.

**Returns**

Integer of the 2's complement for the val.

**Return type**

int

**UnpackAll()** → dict[str, bytes]

Builds a dictionary containing all parts of the POD packet in bytes.

**Raises****Exception** – Nothing to unpack.**Returns**

Dictionary with the command number.

**Return type**

dict[str,bytes]

**PodApi.Packets.Standard module****class** PodApi.Packets.Standard.**PacketStandard**(*pkt: bytes, commands: CommandSet*)Bases: [Packet](#)

Container class that stores a standard command packet for a POD device. The format is STX (1 byte) + command number (4 bytes) + optional payload (? bytes) + checksum (2 bytes) + ETX (1 bytes)

**\_customPayload**

Optional function to translate the payload.

**Type**

Callable[[Any],tuple]|None

**\_customPayloadArgs**

Optional arguments for the \_customPayload.

**Type**

tuple[Any]|None

**payload**

Optional payload from the packet.

**Type**

bytes

**static CheckIfPacketIsValid(msg: bytes)**

Raises an Exception if the packet is incorrectly formatted.

**Parameters**

**msg** (bytes) – Bytes string containing a POD packet. Should begin with STX and end with ETX.

**Raises**

**Exception** – Packet is too small to be a standard packet.

**DefaultPayload()** → tuple[int]

Splits the payload up into its components and translates the binary ASCII encoding into a readable integer.

**Returns**

Tuple with integer values for each component of the payload.

**Return type**

tuple[int]

**static GetMinimumLength()** → int

Gets the number of bytes in the smallest possible packet.

**Returns**

integer representing the minimum length of a standard POD command packet. Format is STX (1 byte) + command number (4 bytes) + optional packet (? bytes) + checksum (2 bytes) + ETX (1 bytes)

**Return type**

int

**static GetPayload(pkt: bytes)** → bytes | None

Gets the payload from a POD packet, if available.

**Parameters**

**pkt** (bytes) – Bytes string containing a POD packet.

**Returns**

Bytes string of the payload, if available.

**Return type**

bytes|None

**HasCustomPayload()** → bool

Checks if a custom payload has been set.

**Returns**

True if there is a custom payload, False otherwise.

**Return type**

bool

**HasPayload()** → bool

Checks if this Packet\_Standard instance has a payload.

**Returns**

True if there is a payload, false otherwise.

**Return type**

bool

**Payload()** → tuple | None

Gets the payload as a readable tuple of values.

**Returns**

Translated payload, if available.

**Return type**

tuple|None

**SetCustomPayload(func: Callable[[Any], tuple], args: tuple[Any] | None = None) → None**

Sets a custom function with optional arguments to translate the payload.

**Parameters**

- **func** (*Callable[[Any], tuple]*) – Function to translate the payload.
- **args** (*tuple[Any], optional*) – Arguments . Defaults to None.

**TranslateAll()** → dict[str, Any]

Builds a dictionary containing all parts of the POD packet in readable values.

**Returns**

Dictionary with the command number and payload.

**Return type**

dict[str,Any]

**UnpackAll()** → dict[str, bytes]

Builds a dictionary containing all parts of the POD packet in bytes.

**Returns**

Dictionary with the command number and payload.

**Return type**

dict[str,bytes]

## Module contents

### 1.1.1.4 PodApi.Parameters package

#### Submodules

#### PodApi.Parameters.Params8206HR module

```
class PodApi.Parameters.Params8206HR.Params8206HR(port: str, sampleRate: int, preamplifierGain: int,  
                                                  lowPass: tuple[int], checkForValidParams: bool =  
                                                  True)
```

Bases: *Params*

Container class that stores parameters for an 8206-HR POD device.

#### **port**

Name of the COM port.

#### **Type**

str

#### **sampleRate**

Sample rate in 100-2000 Hz range.

#### **Type**

int

#### **preamplifierGain**

Preamplifier gain. Should be 10x or 100x.

#### **Type**

int

#### **lowPass**

Low-pass for EEG/EMG in 11-500 Hz range.

#### **Type**

tuple[int]

#### **EEG1()** → int

Gets the filter value of EEG1 in Hz from the low-pass.

#### **Returns**

EEG1 low-pass filter in Hz.

#### **Return type**

int

#### **EEG2()** → int

Gets the filter value of EEG2 in Hz from the low-pass.

#### **Returns**

EEG2 low-pass filter in Hz.

#### **Return type**

int



**EEG3\_EMG()** → int

Gets the filter value of EEG3/EMG in Hz from the low-pass.

**Returns**

EEG3/EMG low-pass filter in Hz.

**Return type**

int

**GetInit()** → str

Builds a string that represents the Params\_8206HR constructor with the arguments set to the values of this class instance.

**Returns**

String that represents the Params\_8206HR constructor.

**Return type**

str

**\_CheckParams()** → None

Throws an exception if Params\_8206HR instance variable is an invalid value.

**Raises**

- **Exception** – Sample rate must be between 100-2000 Hz.
- **Exception** – Preamplidier gain must be 10x or 100x.
- **Exception** – Low-pass EEG/EMG must be between 11-500 Hz.

**lowPassLabels:** tuple[str] = ('EEG1', 'EEG2', 'EEG/EMG')

Tuple describing the items in the lowPass.

## PodApi.Parameters.Params8229 module

```
class PodApi.Parameters.Params8229.Params8229(port: str, systemID: int, motorDirection: bool,
        motorSpeed: int, randomReverse: bool, mode: int,
        reverseBaseTime: int | None = None, reverseVarTime:
        int | None = None, schedule: dict[str, tuple[bool]] |
        None = None, checkForValidParams: bool = True)
```

Bases: [Params](#)

Container class that stores parameters for an 8229 POD device.

**port**

Name of the COM port.

**Type**

str

**systemID**

ID of this 8229 POD system. Must be a positive integer.

**Type**

int

**motorDirection**

False for clockwise and true for counterclockwise.

**Type**  
bool

**motorSpeed**

Motor speed as a percentage 0-100%.

**Type**  
int

**randomReverse**

True to enable random reverse, False otherwise. The random reverse time will be reverseBaseTime + random value in reverseVarTime range.

**Type**  
bool

**reverseBaseTime**

Base time for a random reverse in seconds. Must be a positive integer.

**Type**  
int

**reverseVarTime**

Variable time for a random reverse in seconds. Must be a positive integer.

**Type**  
int

**mode**

System mode; 0 = Manual, 1 = PC Control, and 2 = Internal Schedule.

**Type**  
int

**schedule**

Schedule for a week. The keys are the weekdays (Sunday-Saturday). The values are a tuple of 24 bools that are either 1 for motor on or 0 for motor off

**Type**  
dict[str, tuple[int]]

**static BuildEmptySchedule()** → dict[str, tuple[bool]]

Creates a schedule where the motor is off for all hours of every day.

**Returns**  
Dictionary of the empty schedule. The keys are the days of the week. The values are tuples of 24 zeros.

**Return type**  
dict[str, tuple[bool]]

**GetInit()** → str

Builds a string that represents the Params\_8229 constructor with the arguments set to the values of this class instance.

**Returns**  
String that represents the Params\_Interface constructor.

**Return type**  
str

**\_CheckParams()** → None

Throws an exception if Params\_8229 instance variable is an invalid value.

#### Raises

- **Exception** – The system ID must be a positive integer.
- **Exception** – The motor speed must be between 0-100%.
- **Exception** – The reverse base time (sec) must be a positive integer.
- **Exception** – The reverse variable time (sec) must be a positive integer.
- **Exception** – The mode must be 0, 1, or 2.
- **Exception** – The schedule must have exactly ('Sunday','Monday','Tuesday', 'Wednesday','Thursday','Friday','Saturday') as keys.
- **Exception** – There must be 24 items in the schedule for each day.

**hoursPerDay:** int = 24

Integer storing the number of hours in a day.

**week:** tuple[str] = ('Sunday', 'Monday', 'Tuesday', 'Wednesday', 'Thursday', 'Friday', 'Saturday')

Tuple containing strings of the 7 days of the week.

### PodApi.Parameters.Params8401HR module

```
class PodApi.Parameters.Params8401HR.Params8401HR(port: str, preampDevice: str, sampleRate: int,
                                                    muxMode: bool, preampGain: tuple[int], ssGain:
                                                    tuple[int], highPass: tuple[float], lowPass:
                                                    tuple[int], bias: tuple[float], dcMode: tuple[str],
                                                    checkForValidParams: bool = True)
```

Bases: [Params](#)

Container class that stores parameters for an 8401-HR POD device.

#### port

Name of the COM port.

#### Type

str

#### preampDevice

Name of the mouse/rat preamplifier device.

#### Type

str

#### sampleRate

Sample rate (2000-20000 Hz).

#### Type

int

#### muxMode

Using mux mode when True, false otherwise.

#### Type

bool

**preampGain**

Preamplifier gain (1, 10, or 100) for all channels.

**Type**

tuple[int]

**ssGain**

Second stage gain (1 or 5) for all channels.

**Type**

tuple[int]

**highPass**

High-pass filter (0, 0.5, 1, or 10 Hz) for all channels.

**Type**

tuple[float]

**lowPass**

Low-pass filter (21-15000 Hz) for all channels.

**Type**

tuple[int]

**bias**

Bias voltage (+/- 2.048 V) for all channels.

**Type**

tuple[float]

**dcMode**

DC mode (VBIAS or AGND) for all channels.

**Type**

tuple[str]

**GetInit()** → str

Builds a string that represents the Params\_8401HR constructor with the arguments set to the values of this class instance.

**Returns**

String that represents the Params\_8401HR constructor.

**Return type**

str

**\_CheckParams()** → None

Throws an exception if Params\_8401HR instance variable is an invalid value.

**Raises**

- **Exception** – Mouse/rat preamplifier does not exist.
- **Exception** – Sample rate must be between 2000-20000 Hz.
- **Exception** – EEG/EMG preamplifier gain must be 10x or 100x. For biosensors, the preampGain is None.
- **Exception** – The second stage gain must be 1x or 5x; set gain to None if no-connect.
- **Exception** – The high-pass filter must be 0.5, 1, or 10 Hz. If the channel is DC, input 0.
- **Exception** – The low-pass filter must be between 21-15000 Hz.

- **Exception** – The bias voltage must be +/- 2.048 V.
- **Exception** – The DC mode must be VBIAS or AGND.

**channelLabels:** `tuple[str] = ('A', 'B', 'C', 'D')`

Tuple listing the four channel characters in order.

### PodApi.Parameters.Params8480SC module

```
class PodApi.Parameters.Params8480SC.Params8480SC(port: str, stimulus=tuple[int], preamp=<class
'int'>, ledCurrent=tuple[int], ttlPullups=<class
'int'>, estimCurrent=tuple[int], syncConfig=<class
'int'>, ttlSetup=tuple[int], checkForValidParams:
bool = True)
```

Bases: [Params](#)

Container class that stores parameters for an 8401-HR POD device.

#### **port**

Name of the COM port.

**Type**

str

#### **stimulus**

Stimulus configuration on selected channel.

**Type**

tuple[int]

#### **preamp**

Preamp value (0-1023).

**Type**

int

#### **ledCurrent**

Led-Current (0-600 mA) for both channels.

**Type**

tuple[int]

#### **ttlPullups**

TTL Pullups disabled for value 0, pullups enabled for values that are non-zero.

**Type**

int

#### **estimCurrent**

Estim-Current (0-100 %) for both channels.

**Type**

tuple[int]

#### **syncConfig**

Sets Sync-Config byte.

**Type**

int

**t1Setup**

TTL-Setup for selected channel.

**Type**

tuple[int]

**GetInit()** → str

Builds a string that represents the Params\_8480SC constructor with the arguments set to the values of this class instance.

**Returns**

String that represents the Params\_8480SC constructor.

**Return type**

str

**\_CheckParams()** → None

Throws an exception if Params\_8206HR member variable is an invalid value.

**Raises**

- **Exception** – The preamp must be between 0-1023.
- **Exception** – Led-Curent must be between 0-600.
- **Exception** – Estim-Current must be between 0-100.

**estimCurrent\_CH0()** → int

Gets the estimCurrent value for Channel 0.

**Returns**

Channel 0 estimCurrent in percentage.

**Return type**

int

**estimCurrent\_CH1()** → int

Gets the estimCurrent value for Channel 1.

**Returns**

Channel 1 estimCurrent in percentage.

**Return type**

int

**ledCurrent\_CH0()** → int

Gets the ledCurrent value for Channel 0.

**Returns**

Channel 0 ledCurrent in mA.

**Return type**

int

**ledCurrent\_CH1()** → int

Gets the ledCurrent value for Channel 1.

**Returns**

Channel 1 ledCurrent in mA.

**Return type**

int

**PodApi.Parameters.ParamsBasic module**

**class** PodApi.Parameters.ParamsBasic.**Params**(*port: str, checkForValidParams: bool = True*)

Bases: object

Interface for a container class that stores parameters for a POD device.

**port**

Name of the COM port.

**Type**

str

**GetInit**() → str

Builds a string that represents the Params\_Interface constructor with the arguments set to the values of this class instance.

**Returns**

String that represents the Params\_Interface constructor.

**Return type**

str

**\_CheckParams**() → None

Throws an exception if Params\_Interface instance variable is an invalid value.

**Raises**

**Exception** – The port name must begin with COM.

**static** **\_FixTypeInTuple**(*arr: tuple, itemType: type*) → tuple['type']

Retypes each item of the arr arguemnt to itemType.

**Parameters**

- **arr** (*tuple*) – Tuple of items to be re-typed.
- **itemType** (*type*) – Type to be casted to each tuple item.

**Returns**

Tuple with values of all itemType types.

**Return type**

tuple[type]

**Module contents****1.1.2 Module contents****1.2 Setup package****1.2.1 Subpackages****1.2.1.1 Setup.Inputs package****Submodules**

## Setup.Inputs.GetUserInput module

**class** Setup.Inputs.GetUserInput.**UserInput**

Bases: object

UserInput contains several methods for getting user input for POD device setup.

**static AskForBool**(*prompt: str*) → bool

Asks user for bool type input.

**Parameters**

**prompt** (*str*) – Statement requesting input from the user.

**Returns**

Boolean type input from user.

**Return type**

bool

**static AskForFloat**(*prompt: str*) → float

Asks user for float type input.

**Parameters**

**prompt** (*str*) – Statement requesting input from the user.

**Returns**

Float type input from user.

**Return type**

float

**static AskForFloatInList**(*prompt: str, goodInputs: list, badInputMessage: str | None = None*) → float

Asks the user for a float that exists in the list of valid options.

**Parameters**

- **prompt** (*str*) – Statement requesting input from the user.
- **goodInputs** (*list*) – List of valid input options.
- **badInputMessage** (*str | None, optional*) – Error message to be printed if invalid input is given. Defaults to None.

**Returns**

User's choice from the options list as a float.

**Return type**

float

**static AskForFloatInRange**(*prompt: str, minimum: float, maximum: float, thisIs: str = 'Input', unit: str = ''*) → float

Asks the user for an float value that falls in a given range.

**Parameters**

- **prompt** (*str*) – Statement requesting input from the user.
- **minimum** (*float*) – Minimum value of range.
- **maximum** (*float*) – Maximum value of range.
- **thisIs** (*str, optional*) – Description of the input/what is being asked for. Used when printing the error message. Defaults to 'Input'.



- **unit** (*str*, *optional*) – Unit of the requested value. Use when printing the error message. Defaults to ‘.’.

**Returns**

Float value given by the user that falls in the given range.

**Return type**

float

**static AskForInput**(*prompt: str*, *append: str = ':'*) → *str*

Asks user for input given a prompt. Will append a colon ‘:’ to the end of prompt by default

**Parameters**

- **prompt** (*str*) – Statement requesting input from the user
- **append** (*str*, *optional*) – Appended to the end of the prompt. Defaults to ‘:’.

**Returns**

String of the user input

**Return type**

str

**static AskForInt**(*prompt: str*) → *int*

Asks user for int type input.

**Parameters**

**prompt** (*str*) – Statement requesting input from the user.

**Returns**

Integer type input from user.

**Return type**

int

**static AskForIntInList**(*prompt: str*, *goodInputs: list*, *badInputMessage: str | None = None*) → *int*

Asks the user for an integer that exists in the list of valid options.

**Parameters**

- **prompt** (*str*) – Statement requesting input from the user
- **goodInputs** (*list*) – List of valid input options.
- **badInputMessage** (*str | None*, *optional*) – Error message to be printed if invalid input is given. Defaults to None.

**Returns**

User’s choice from the options list as an integer.

**Return type**

int

**static AskForIntInRange**(*prompt: str*, *minimum: int*, *maximum: int*, *thisIs: str = 'Input'*, *unit: str = ''*) → *int*

Asks the user for an integer value that falls in a given range.

**Parameters**

- **prompt** (*str*) – Statement requesting input from the user.
- **minimum** (*int*) – Minimum value of range.
- **maximum** (*int*) – Maximum value of range.

- **thisIs** (*str*, *optional*) – Description of the input/what is being asked for. Used when printing the error message. Defaults to ‘Input’.
- **unit** (*str*, *optional*) – Unit of the requested value. Use when printing the error message. Defaults to ‘’.

**Returns**

Integer value given by the user that falls in the given range.

**Return type**

int

**static AskForStrInList** (*prompt: str*, *goodInputs: list*, *badInputMessage: str | None = None*) → str

Asks the user for a string that exists in the list of valid options.

**Parameters**

- **prompt** (*str*) – Statement requesting input from the user.
- **goodInputs** (*list*) – List of valid input options
- **badInputMessage** (*str | None*, *optional*) – Error message to be printed if invalid input is given. Defaults to None.

**Returns**

User’s choice from the options list as a string.

**Return type**

str

**static AskForType** (*typecast: function*, *prompt: str*) → int | float | str

Ask user for input of a specific data type. If invalid input is given, an error message will print and the user will be prompted again.

**Parameters**

- **typecast** (*function*) – Datatype to cast the user input (ex. `_CastInt`, `_CastFloat`, `_CastStr`)
- **prompt** (*str*) – Statement requesting input from the user

**Returns**

Input from user as the requested type.

**Return type**

int|float|str

**static AskForTypeInList** (*typecast: function*, *prompt: str*, *goodInputs: list*, *badInputMessage: str | None = None*) → int | float | str

Asks the user for a value of a given type that exists in the list of valid options. If invalid input is given, an error message will print and the user will be prompted again.

**Parameters**

- **typecast** (*function*) – Datatype to cast the user input (ex. `_CastInt`, `_CastFloat`, `_CastStr`).
- **prompt** (*str*) – Statement requesting input from the user.
- **goodInputs** (*list*) – List of valid input options.
- **badInputMessage** (*str | None*, *optional*) – Error message to be printed if invalid input is given. Defaults to None.

**Returns**

User's choice from the goodInputs list as the given datatype.

**Return type**

int|float|str

**static AskForTypeInRange**(*typecast: function, prompt: str, minimum: int | float, maximum: int | float, thisIs: str = 'Input', unit: str = ''*) → int | float

Asks user for a numerical value that falls between two numbers. If invalid input is given, an error message will print and the user will be prompted again.

**Parameters**

- **typecast** (*function*) – Datatype to cast the user input (ex. `_CastInt`, `_CastFloat`, `_CastStr`).
- **prompt** (*str*) – Statement requesting input from the user.
- **minimum** (*int | float*) – Minimum value of range.
- **maximum** (*int | float*) – Maximum value of range.
- **thisIs** (*str, optional*) – Description of the input/what is being asked for. Used when printing the error message. Defaults to 'Input'.
- **unit** (*str, optional*) – Unit of the requested value. Use when printing the error message. Defaults to ''.

**Returns**

Numerical value given by the user that falls in the given range.

**Return type**

int|float

**static AskYN**(*question: str, append: str = '(y/n): '*) → bool

Asks the user a yes or no question. If invalid input is given, an error message will print and the user will be prompted again.

**Parameters**

- **question** (*str*) – Statement requesting input from the user.
- **append** (*str, optional*) – Appended to the end of the question. Defaults to '(y/n): '.

**Returns**

True for yes, false for no.

**Return type**

bool

**static CastFloat**(*value*) → float

Casts the argument as an float.

**Parameters**

**value** – Value to type casted.

**Returns**

Value type casted as a float.

**Return type**

float

**static CastInt**(*value*) → int

Casts the argument as an integer.

**Parameters**

**value** – Value to type casted.

**Returns**

Value type casted as an integer.

**Return type**

int

**static CastStr**(*value*) → str

Casts the argument as a string.

**Parameters**

**value** – Value to type casted.

**Returns**

Value type casted as a string.

**Return type**

str

**static CheckFileExt**(*f: str, flsExt: bool = True, goodExt: list[str] = ['.txt'], printErr: bool = True*) → bool

Checks if a file name has a valid extension.

**Parameters**

- **f** (*str*) – file name or extension
- **flsExt** (*bool, optional*) – Boolean flag that is true if f is an extension, false otherwise. Defaults to True.
- **goodExt** (*list[str], optional*) – List of valid file extensions. Defaults to ['.txt',].
- **printErr** (*bool, optional*) – Boolean flag that, when true, will print an error statement. Defaults to True.

**Returns**

True if extension is in goodExt list, False otherwise.

**Return type**

bool

**static GetFileName**(*goodExt: list[str] = ['.txt']*) → str

Asks the user for a filename. :param goodExt: List of valid file extensions. Defaults to ['.txt']. :type goodExt: list[str], optional

**Returns**

String of the file name and extension.

**Return type**

str

**static GetFilePath**(*prompt: str | None = None, goodExt: list[str] = ['.txt']*) → str

Asks the user for a file path and file name.

**Parameters**

- **prompt** (*str | None, optional*) – Text to print to the user before requesting the path. Defaults to None.

- **goodExt** (*list[str], optional*) – List of valid file extensions. Defaults to ['.txt'].

**Returns**

File path and name.

**Return type**

str

**Module contents****1.2.1.2 Setup.SetupOneDevice package****Submodules****Setup.SetupOneDevice.Setup\_8206HR module**

**class** Setup.SetupOneDevice.Setup\_8206HR.Setup\_8206HR

Bases: [Setup\\_Interface](#)

Setup\_8206HR provides the setup functions for an 8206-HR POD device.

**static** **GetDeviceName()** → str

Returns the name of the 8206-HR POD device.

**Returns**

8206-HR.

**Return type**

str

**StopStream()** → None

Write a command to stop streaming data to all POD devices.

**static** **\_ChoosePreampGain()** → int

Asks user for the preamplifier gain of their POD device.

**Returns**

Integer 10 or 100 for the preamplifier gain.

**Return type**

int

**\_ConnectPODdevice**(*deviceNum: int, deviceParams: [Params8206HR](#)*) → bool

Creates a POD\_8206HR object and write the setup parameters to it.

**Parameters**

- **deviceNum** (*int*) – Integer of the device's number.
- **deviceParams** (*Params\_8206HR*) – Device's parameters.

**Returns**

True if connection was successful, false otherwise.

**Return type**

bool

**\_GetPODdeviceParameterTable()** → Texttable

Builds a table containing the parameters for all POD devices.

**Returns**

Texttable containing all parameters.

**Return type**

Texttable

**\_GetParam\_onePODdevice**(*forbiddenNames: list[str] = []*) → *Params8206HR*

Asks the user to input all the device parameters.

**Parameters**

**forbiddenNames** (*list[str]*) – List of port names already used by other devices.

**Returns**

Dictionary of device parameters.

**Return type**

dict[str,(str|int|dict[str,int])]

**\_OpenSaveFile\_EDF**(*fname: str, devNum: int*) → EdfWriter

Opens EDF file and write header.

**Parameters**

- **fname** (*str*) – String filename.
- **devNum** (*int*) – Integer device number.

**Returns**

Opened file.

**Return type**

EdfWriter

**\_OpenSaveFile\_TXT**(*fname: str*) → IOBase

Opens a text file and writes the column names. Writes the current date/time at the top of the txt file.

**Parameters**

**fname** (*str*) – String filename.

**Returns**

Opened file.

**Return type**

IOBase

**\_PHYSICAL\_BOUND\_uV: int = 2046**

Class-level integer representing the max/-min physical value in uV. Used for EDF files.

**\_StreamThreading**() → dict[int, threading.Thread]

Opens a save file, then creates a thread for each device to stream and write data from.

**Returns**

Dictionary with keys as the device number and values as the started Thread.

**Return type**

dict[int,Thread]

**\_StreamUntilStop**(*pod: Pod8206HR, file: IOBase | EdfWriter, sampleRate: int*) → None

Streams data from a POD device and saves data to file. Stops looking when a stop stream command is read. Calculates average time difference across multiple packets to collect a continuous time series data.

**Parameters**

- **pod** (*POD\_8206HR*) – POD device to read from.
- **file** (*IOBase* | *EdfWriter*) – open file.
- **sampleRate** (*int*) – Integer sample rate in Hz.

**static \_WriteDataToFile\_EDF**(*file: EdfWriter, data: list[*numpy.ndarray*]*)

Writes data to an open EDF file.

**Parameters**

- **file** (*EdfWriter*) – opened EDF file.
- **data** (*list[*np.ndarray*]*) – List of 3 items, one for each channel.

**static \_WriteDataToFile\_TXT**(*file: IOBase, data: list[*numpy.ndarray*], t: *ndarray**)

Writes data to an open text file.

**Parameters**

- **file** (*IOBase*) – opened write file.
- **data** (*list[*np.ndarray*]*) – List of 3 items, one for each channel.
- **t** (*np.ndarray*) – list with the time stamps (in seconds).

**Setup.SetupOneDevice.Setup\_8229 module**

**class** Setup.SetupOneDevice.Setup\_8229.**Setup\_8229**

Bases: [Setup\\_Interface](#)

Setup\_8229 provides the setup functions for an 8229 POD device.

**\_streamMode**

True when the user wants to stream data from an 8229 POD device, False otherwise.

**Type**

bool

**static GetDeviceName**() → str

Returns the name of the POD device.

**Returns**

8229.

**Return type**

str

**static GetSupportedFileExtensions**() → list[str]

Returns a list containing valid file extensions.

**Returns**

List of string file extensions.

**Return type**

list[str]

**StopStream**() → None

Update the state flag to signal to stop streaming data.

**\_ConnectPODdevice**(*deviceNum: int, deviceParams: Params8229*) → bool

Creates a 8229 POD device object and write the setup parameters to it.

**Parameters**

- **deviceNum** (*int*) – Integer of the device’s number.
- **deviceParams** (*Params\_8229*) – Device parameters.

**Returns**

True if connection was successful, false otherwise.

**Return type**

bool

**\_GetPODdeviceParameterTable**() → Texttable

Builds a table containing the parameters for all POD devices.

**Returns**

Table containing all parameters.

**Return type**

Texttable

**\_GetParam\_onePODdevice**(*forbiddenNames: list[str] = []*) → *Params8229*

Asks the user to input all the device parameters.

**Parameters**

**forbiddenNames** (*list[str], optional*) – List of port names already used by other devices. Defaults to [].

**Returns**

Device parameters.

**Return type**

Params\_8229

**static \_GetScheduleForWeek**() → dict[str, tuple[int]]

Asks the user to input if the motor is on/off for each hour of each day of the week.

**Returns**

Dictionary with the schedule. The keys are the days of the week (Sunday, Monday, ...). The values are a tuple of 24 items for each hour; the items are 1 if the motor is on or 0 if the motor is off

**Return type**

dict[str, tuple[int]]

**\_OpenSaveFile\_EDF**(*fname: str, devNum: int*)

EDF files are not supported for 8229 POD devices. Overwrites the parent’s method, which would open an EDF file and write the header.

**Parameters**

- **fname** (*str*) – String filename. Not used.
- **devNum** (*int*) – Integer device number. Not used.

**Raises**

**Exception** – EDF filetype is not supported for 8229 POD devices.



**\_OpenSaveFile\_TXT**(*fname: str*) → IOBase

Opens a text file and writes the column names. Writes the current date/time at the top of the txt file.

**Parameters**

**fname** (*str*) – String filename.

**Returns**

Opened file.

**Return type**

IOBase

**\_StreamThreading**() → dict[int, threading.Thread]

Opens a save file, then creates a thread for each device to stream and write data from.

**Returns**

Dictionary with keys as the device number and values as the started Thread.

**Return type**

dict[int, Thread]

**\_StreamUntilStop**(*pod: Pod8229, file: IOBase*) → None

Saves a log of all packets recieved from the 8229 POD device until the user decides to stop streaming.

**Parameters**

- **pod** (*POD\_8229*) – POD device to read from.
- **file** (*IOBase*) – Opened text file to save data to.

## Setup.SetupOneDevice.Setup\_8401HR module

**class** Setup.SetupOneDevice.Setup\_8401HR.**Setup\_8401HR**

Bases: [Setup\\_Interface](#)

Setup\_8401HR provides the setup functions for an 8206-HR POD device. REQUIRES FIRMWARE 1.0.2 OR HIGHER.

**static** **GetDeviceName**() → str

returns the name of the 8401-HR POD device.

**Returns**

8401-HR.

**Return type**

str

**StopStream**() → None

Write a command to stop streaming data to all POD devices.

**static** **\_CodeDCmode**(*dcMode: str*) → int

gets the integer payload to use for ‘SET DC MODE’ commands given the mode.

**Parameters**

**dcMode** (*str*) – DC mode VBIAS or AGND.

**Raises**

**Exception** – DC Mode value is not supported.

**Returns**

Integer code representing the DC mode.

**Return type**

int

**static** `_CodeHighpass`(*highpass*: float) → int

Gets the integer payload to use for 'SET HIGHPASS' given a highpass value.

**Parameters****highpass** (float) – Highpass value in Hz.**Raises****Exception** – High-pass value is not supported.**Returns**

Integer code representing the highpass value.

**Return type**

int

**\_ConnectPODdevice**(*deviceNum*: int, *deviceParams*: Params8401HR) → bool

Creates a POD\_8206HR object and write the setup parameters to it.

**Parameters**

- **deviceNum** (int) – Integer of the device's number.
- **deviceParams** (Params\_8401HR) – Device parameters.

**Returns**

True if connection was successful, false otherwise.

**Return type**

bool

**\_GetPODdeviceParameterTable**() → Texttable

Builds a table containing the parameters for all POD devices.

**Returns**

Texttable containing all parameters.

**Return type**

Texttable

**\_GetParam\_onePODdevice**(*forbiddenNames*: list[str] = []) → Params8401HR

Asks the user to input all the device parameters.

**Parameters****forbiddenNames** (list[str]) – List of port names already used by other devices.**Returns**

Device parameters.

**Return type**

Params\_8401HR

**\_GetPreampDeviceName**() → str

Asks the user to select a mouse/rat preamplifier.

**Returns**

String of the chosen preamplifier.

**Return type**

str

**`_NiceABCDtableText`**(*abcdValues: list[int | str | None], channelMap: dict[str, str]*) → str

Builds a string that formats the channel values to be input into the parameter table.

**Parameters**

- **`abcdValueDict`** (*dict[str, int | str | None]*) – Dictionary with ABCD keys.
- **`channelMap`** (*dict[str, str]*) – Maps the ABCD channels to the sensor’s channel name.

**Returns**

String with “channel name: value newline...” for each channel.

**Return type**

str

**`_OpenSaveFile_EDF`**(*fname: str, devNum: int*) → EdfWriter

Opens EDF file and write header.

**Parameters**

- **`fname`** (*str*) – String filename.
- **`devNum`** (*int*) – Integer device number.

**Returns**

Opened file.

**Return type**

EdfWriter

**`_OpenSaveFile_TXT`**(*fname: str*) → IOBase

Opens a text file and writes the column names. Writes the current date/time at the top of the txt file.

**Parameters**

**`fname`** (*str*) – String filename.

**Returns**

Opened file.

**Return type**

IOBase

**`_PHYSICAL_BOUND_uV: int = 2046`**

Class-level integer representing the max/-min physical value in uV. Used for EDF files.

**`static _SetBias`**(*channelName: str*) → float

Asks the user for the bias voltage in V (+/-2.048V).

**Parameters**

**`channelName`** (*str*) – Name of the channel.

**Returns**

A float for the bias voltage in V.

**Return type**

float

**`static _SetDCMode`**(*channelName: str*) → str

Asks the user for the DC mode (VBIAS or AGND).

**Parameters**

**`channelName`** (*str*) – Name of the channel.

**Returns**

String for the DC mode.

**Return type**

str

**\_SetForMappedChannels**(*message: str, channelMap: dict[str, str], func: function*) → tuple[int | None]

Asks the user to input values for all channels (excluding no-connects).

**Parameters**

- **message** (*str*) – Message to ask the user.
- **channelMap** (*dict[str, str]*) – Maps the ABCD channels to the sensor's channel name.
- **func** (*function*) – a function that asks the user for an input. takes one string parameter and returns one value.

**Returns**

Tuple with user inputs for values for the ABCD channels

**Return type**

tuple[int|None]

**static \_SetHighpass**(*channelName: str*) → float | None

Asks the user for the high-pass in Hz (0.5,1,10Hz, or DC).

**Parameters**

**channelName** (*str*) – Name of the channel.

**Returns**

A float for the high-pass frequency in Hz, or None if DC.

**Return type**

float|None

**static \_SetLowpass**(*channelName: str*) → int | None

Asks the user for the low-pass in Hz (21-15000Hz).

**Parameters**

**channelName** (*str*) – Name of the channel.

**Returns**

An integer for the low-pass frequency in Hz.

**Return type**

int|None

**static \_SetPreampGain**(*channelName: str*) → int | None

Asks the user for the preamplifier gain.

**Parameters**

**channelName** (*str*) – Name of the channel.

**Returns**

An integer for the gain, or None if no gain.

**Return type**

int|None

**static \_SetSSGain**(*channelName: str*) → int

Asks the user for the second stage gain.

**Parameters**

**channelName** (*str*) – Name of the channel.

**Returns**

An integer for the gain.

**Return type**

int

**\_StreamThreading()** → dict[int, threading.Thread]

Opens a save file, then creates a thread for each device to stream and write data from.

**Returns**

Dictionary with keys as the device number and values as the started Thread.

**Return type**

dict[int, Thread]

**\_StreamUntilStop**(*pod: Pod8401HR, file: IOBase | EdfWriter, sampleRate: int, devNum: int*) → None

Streams data from a POD device and saves data to file. Stops looking when a stop stream command is read. Calculates average time difference across multiple packets to collect a continuous time series data.

**Parameters**

- **pod** (*POD\_8206HR*) – POD device to read from.
- **file** (*IOBase | EdfWriter*) – open file.
- **sampleRate** (*int*) – Integer sample rate in Hz.

**static \_WriteDataToFile\_EDF**(*file: EdfWriter, data: list[numpy.ndarray]*)

Writes data to an open EDF file.

**Parameters**

- **file** (*EdfWriter*) – opened EDF file.
- **data** (*list[np.ndarray]*) – List with one item for each channel.

**static \_WriteDataToFile\_TXT**(*file: IOBase, data: list[numpy.ndarray], t: ndarray*)

Writes data to an open text file.

**Parameters**

- **file** (*IOBase*) – opened write file.
- **data** (*list[np.ndarray]*) – List with one item for each channel.
- **t** (*np.ndarray*) – list with the time stamps (in seconds).

## Setup.SetupOneDevice.Setup\_8480SC module

**class** Setup.SetupOneDevice.Setup\_8480SC.**Setup\_8480SC**

Bases: [Setup\\_Interface](#)

Setup\_8480SC provides the setup functions for an 8480-SC POD device.

**\_streamMode**

Set to True when the user wants to start streaming data from an 8480 POD device, False otherwise.

**Type**

bool

**static GetDeviceName()** → str

Returns the name of the POD device.

**Returns**

String of \_NAME.

**Return type**

str

**static GetSupportedFileExtensions()** → list[str]

Returns a list containing valid file extensions.

**Returns**

List of string file extensions.

**Return type**

list[str]

**StopStream()** → None

Update the state flag to signal to stop streaming data.

**static \_ChoosePeriod()** → tuple[int]

Asks the user an input value for Stimulus Period, which is then seperated into Period\_ms and Period\_us. Seperation is required because the 'SET STIMULUS' requires 7 items in payload, and the second and third items is the Period\_ms and Period\_us.

**Returns**

Formatted period into millisecs and microsecs.

**Return type**

tuple[int]

**static \_ChooseRepeat()** → int

Asks the user to input a value for the number of stimulus pulses to perform.

**Returns**

representing repeat count for command 'SET STIMULUS'.

**Return type**

int

**static \_ChooseStimulusConfig()** → int

Asks the user to input values for Config format of Stimulus

**Returns**

Formatted Stimulus Config value.

**Return type**

(int)

**static \_ChooseSyncConfig()** → int

Asks the user to input values for Sync Config bits.

**Returns**

Value calculated from the input bits, this value would be given as payload.

**Return type**

int

**static** `_ChooseWidth()` → tuple[int]

Asks the user an input value for Stimulus width, which is then separated into width\_ms and width\_us. Separation is required because the 'SET STIMULUS' requires 7 items in payload, and the fourth and fifth items is the width\_ms and width\_us.

**Returns**

Formatted given width into millisecs and microsecs

**Return type**

tuple[int]

**\_ConnectPODdevice**(*deviceNum*: int, *deviceParams*: Params8480SC) → bool

Creates a POD\_8206HR object and write the setup parameters to it.

**Parameters**

- **deviceNum** (int) – Integer of the device's number.
- **deviceParams** (Params\_8480SC) – Device parameters.

**Returns**

True if connection was successful, false otherwise.

**Return type**

bool

**\_GetPODdeviceParameterTable**() → Texttable

Builds a table containing the parameters for all POD devices.

**Returns**

Texttable containing all parameters.

**Return type**

Texttable

**\_GetParam\_onePODdevice**(*forbiddenNames*: list[str] = []) → Params8480SC

Asks the user to input all the device parameters.

**Parameters**

**forbiddenNames** (list[str]) – List of port names already used by other devices.

**Returns**

Device parameters.

**Return type**

Params\_8480SC

**\_OpenSaveFile\_EDF**(*fname*: str, *devNum*: int)

EDF files are not supported for 8480 POD devices. Overwrites the parent's method, which would open an EDF file and write the header.

**Parameters**

- **fname** (str) – String filename. Not used.
- **devNum** (int) – Integer device number. Not used.

**Raises**

**Exception** – EDF filetype is not supported for 8480 POD devices.

**\_OpenSaveFile\_TXT**(*fname*: str) → IOBase

Opens a text file and writes the column names. Writes the current date/time at the top of the txt file.

**Parameters**

**fname** (*str*) – String filename.

**Returns**

Opened file.

**Return type**

IOBase

**\_StreamThreading()** → dict[int, threading.Thread]

Opens a save file, then creates a thread for each device to stream and write data from.

**Returns**

Dictionary with keys as the device number and values as the started Thread.

**Return type**

dict[int, Thread]

**\_StreamUntilStop**(*pod*: Pod8480SC, *file*: IOBase) → None

Saves a log of all packets recieved from the 8480 POD device until the user decides to stop streaming.

**Parameters**

- **pod** (*POD\_8480*) – POD device to read from.
- **file** (*IOBase*) – Opened text file to save data to.

**static \_TtlSetup()** → int

Asks the user to input values for Config format values of TTL Setup.

**Returns**

Formatted TTL Setup config value.

## Setup.SetupOneDevice.Setup\_PodInterface module

**class Setup.SetupOneDevice.Setup\_PodInterface.Setup\_Interface**

Bases: object

Setup\_Interface provides the basic interface of required methods for subclasses to implement. SetupPodDevices.py is designed to handle any of these children.

**\_podDevices**

Instance-level dictionary of pod device objects. MUST have keys as device number.

**Type**

dict[int, POD\_Basics]

**\_podParametersDict**

Instance-level dictionary of device information. MUST have keys as device number.

**Type**

dict[int, Params\_Interface]

**\_saveFileName**

Instance-level string filename: <path>/file.ext. The device name and number will be appended to the filename.

**Type**

str



**AddPODdevice()** → None

Asks the user for the parameters for the new device. A new device number is generated.

**AreDeviceParamsValid**(*paramDict*: None | dict[int, PodApi.Parameters.ParamsBasic.Params/])

Checks if the parameters dictionary is valid.

**Parameters**

**paramDict** (None | dict[int, Params\_Interface]) – Dictionary of parameters for all POD devices.

**Raises**

- **Exception** – Parameters must be contained in a dictionary.
- **Exception** – Device keys must be integer type.
- **Exception** – Device parameters must be dictionary type.
- **Exception** – Device parameters dictionary is empty.

**ConnectAllPODdevices()** → bool

Connects all setup POD devices.

**Returns**

True if all devices are successfully connected, false otherwise.

**Return type**

bool

**DisplayPODdeviceParameters()** → None

Display all the pod device parameters in a table.

**static GetDeviceName()** → str

returns the name of the POD device.

**Returns**

GENERIC.

**Return type**

str

**GetPODparametersInit()** → str

Gets a dictionary whose keys are the device number and the value is the device parameters dict. :returns: String representation of a dictionary of POD device parameters. The keys are the device number. :rtype: str

**GetSaveFileName()** → str

Gets the path and filename where streaming data is saved to.

**Returns**

String of the save file name and path (\_saveFileName).

**Return type**

str

**static GetSupportedFileExtensions()** → list[str]

Returns a list containing valid file extensions.

**Returns**

List of string file extensions.

**Return type**

list[str]

**PrintSaveFile()** → None

Prints the file path and name that data is saved to. Note that the device name and number will be appended to the end of the filename.

**SetupFileName**(*fileName: str | None = None*) → None

Gets the file path and name to save streaming data to. Note that the device name and number will be appended to the end of the filename.

**Parameters**

**fileName** (*str | None, optional*) – Name and path of the file, if known. Defaults to None.

**SetupPODparameters**(*podParametersDict: dict[int, PodApi.Parameters.ParamsBasic.Params] | None = None*) → None

Sets the parameters for the POD devices.

**Parameters**

**podParametersDict** (*dict[int, Params\_Interface] | None, optional*) – dictionary of the device parameters for all devices. Defaults to None.

**StopStream()** → None

Write a command to stop streaming data to all POD devices.

**Stream()** → dict[int, threading.Thread]

Tests that all devices are connected then starts streaming data.

**Raises**

**Exception** – Test connection failed.

**Returns**

Dictionary with integer device number keys and Thread values.

**Return type**

dict[int, Thread]

**ValidateParams()** → None

Displays a table of the parameters of all devices, then asks the user if everything is correct. The user can then edit the parameters of a device.

**\_BuildFileName**(*devNum: int*) → str

Appends the device name and number to the end of the file name.

**Parameters**

**devNum** (*int*) – Integer of the device number.

**Returns**

String file name.

**Return type**

str

**static \_ChoosePort**(*forbidden: list[str] = []*) → str

Systems checks user's Operating System, and chooses ports accordingly.

**Parameters**

**forbidden** (*list[str], optional*) – List of port names that are already used. Defaults to [].

**Returns**

String name of the port.

**Return type**

str

**\_ConnectPODdevice**(*deviceNum*: int, *deviceParams*: Params) → bool

Creates a POD device object and write the setup parameters to it.

**Parameters**

- **deviceNum** (int) – Integer of the device's number.
- **deviceParams** (Params\_Interface) – Device parameters.

**Returns**

True if connection was successful, false otherwise.

**Return type**

bool

**\_DisconnectAllPODdevices**() → None

Disconnects all POD devices by deleted all POD obejcts.

**\_EditParams**() → None

Asks the user which device to edit, and then asks them to re-input the device parameters.

**\_GetForbiddenNames**(*exclude*: str | None = None) → list[str]

Generates a list of port names used by the active pod devices. There is an option to exclude an additional name from the list.

**Parameters**

**exclude** (str | None, optional) – String port name to exclude from the returned list.  
Defaults to None.

**Returns**

List of string names of ports in use.

**Return type**

list[str]

**\_GetPODdeviceParameterTable**() → Texttable

Builds a table containing the parameters for all POD devices.

**Returns**

Table containing all parameters.

**Return type**

Texttable

**\_GetParam\_allPODdevices**() → dict[int, PodApi.Parameters.ParamsBasic.Params]

First gets the number of POD devices, then asks the user for the information for each device.

**Returns**

Dictionary with device numbers for keys and parameters for values.

**Return type**

dict[int, Params\_Interface]

**\_GetParam\_onePODdevice**(*forbiddenNames*: list[str] = []) → Params

Asks the user to input all the device parameters.

**Parameters**

**forbiddenNames** (list[str]) – List of port names already used by other devices. Defaults to [].

**Returns**

Device parameters.

**Return type**

Params\_Interface

**static** `_GetTimeHeader_forTXT()` → str

Builds a string containing the current date and time to be written to the text file header.

**Returns**

String containing the date and time. Each line begins with '#' and ends with a newline.

**Return type**

str

**\_OpenSaveFile**(*devNum: int*) → IOBase | EdfWriter

Opens a save file for a given device.

**Parameters**

**devNum** (*int*) – Integer of the device number.

**Returns**

Opened file. IOBase for a text file, or EdfWriter for EDF file.

**Return type**

IOBase | EdfWriter

**\_OpenSaveFile\_EDF**(*fname: str, devNum: int*) → EdfWriter

Opens EDF file and write header.

**Parameters**

- **fname** (*str*) – String filename.
- **devNum** (*int*) – Integer device number.

**Returns**

Opened file.

**Return type**

EdfWriter

**\_OpenSaveFile\_TXT**(*fname: str*) → IOBase

Opens a text file and writes the column names. Writes the current date/time at the top of the txt file.

**Parameters**

**fname** (*str*) – String filename.

**Returns**

Opened file.

**Return type**

IOBase

**static** `_PrintDeviceNumber`(*num: int*) → None

Prints a title with the device number.

**Parameters**

**num** (*int*) – Integer of the device number.

**\_RemoveDevice**() → None

Asks the user for a device number to remove, then deletes that device. This will only remove a device if there are more than one options.

**\_SelectDeviceFromDict**(*action: str*) → int

Asks the user to select a valid device number. The input must be an integer number of an existing device.

**Parameters**

**action** (*str*) – Description of the action to be performed on the device.

**Returns**

Integer for the device number.

**Return type**

int

**static \_SetNumberOfDevices**(*name: str*) → int

Asks the user for how many devices they want to setup.

**Parameters**

**name** (*str*) – Name of the POD device type.

**Returns**

Integer number of POD devices desired by the user.

**Return type**

int

**\_StreamThreading**() → dict[int, threading.Thread]

Opens a save file, then creates a thread for each device to stream and write data from.

**Returns**

Dictionary with keys as the device number and values as the started Thread.

**Return type**

dict[int, Thread]

**static \_TestDeviceConnection**(*pod: Pod*) → bool

Tests if a POD device can be read from or written. Sends a PING command.

**Parameters**

**pod** (*POD\_Basics*) – POD device to read to and write from.

**Returns**

True for successful connection, false otherwise.

**Return type**

bool

**\_TestDeviceConnection\_All**() → bool

Tests the connection of all setup POD devices.

**Returns**

True when all devices are successfully connected, false otherwise

**Return type**

bool

**static \_uV**(*voltage: float | int*) → float

Converts volts to microVolts, rounded to 6 decimal places.

**Parameters**

**voltage** (*float | int*) – number of volts.

**Returns**

voltage in of uV.

**Return type**  
float

## Module contents

### 1.2.2 Submodules

#### 1.2.3 Setup.Setup\_PodDevices module

```
class Setup.Setup_PodDevices.Setup_PodDevices(saveFileDict: dict[str, str] | None = None,  
                                              podParametersDict: dict[str, dict[int,  
PodApi.Parameters.ParamsBasic.Params] | None] | None  
                                              = None)
```

Bases: object

Setup\_PodDevices allows a user to set up and stream from any number of POD devices. The streamed data is saved to a file.

REQUIRES FIRMWARE 1.0.2 OR HIGHER.

##### **\_Setup\_PodDevices**

Dictionary containing the Setup\_Interface subclasses for each POD device.

**Type**  
dict[str,dict[int,Params\_Interface]

##### **\_options**

Dictionary listing the different options for the user to complete.

**Type**  
dict[int,str]

**GetOptions()** → dict[int, str]

Gets the dictionary of setup options.

**Returns**  
Dictionary listing the different options for the user to complete (\_options).

**Return type**  
dict[int,str]

**GetPODparametersInit()** → str

Sets up each POD device type. Used in initialization. :returns: String representing a dictionary of all POD devices initialization. The keys are the device name and the entries are the initialization dictionaries. :rtype: str

**GetSaveFileNames()** → str

**Run()** → None

Prints the options and asks the user what to do. Loops until ‘Quit’ is chosen.

**SetupPODparameters**(podParametersDict: dict[str, dict[int, PodApi.Parameters.ParamsBasic.Params] | None]) → None

Sets up each POD device type. Used in initialization.

**Parameters**

**podParametersDict** (*dict[str, dict[int, Params\_Interface] | None]*) – Dictionary of all POD devices initialization. The keys are the device name and the entries are the initialization dictionaries.

**SetupSaveFile**(*saveFileDict: dict[str, str] | None = None*) → None

Gets the path/file name from the user and stores it. Used in initialization.

**Parameters**

**saveFile** (*dict[str, str | None] | None, optional*) – String of the save file, which includes the directory path, filename, and file extension. Defaults to None.

**\_AskOption**() → int

Asks user which option to do.

**Returns**

Integer number representing an option key.

**Return type**

int

**\_AskToStopStream**() → None

Asks user to press enter to stop streaming. The program will then prompt all POD devices to end stream.

**static \_AskUserForDevices**()

Asks the user what POD devices they want to use.

**\_CheckForValidParams**(*podParametersDict: dict[str, None | dict]*) → bool

Checks if the parameters are correctly formatted.

**Parameters**

**podParametersDict** (*dict[str, None | dict]*) – Dictionary with keys as the device names and values as None or the respective parameter dictionary.

**Raises**

- **Exception** – Parameters must be dictionary type.
- **Exception** – Parameters dictionary is empty.
- **Exception** – Invalid device name in parameter dictionary.

**Returns**

True if the parameters are correctly formatted.

**Return type**

bool

**\_ConnectNewDevice**() → None

Asks the user for the POD device type, then it sets up that device.

**\_DoOption**(*choice: int*) → None

Performs the methods associated with the user selected option.

**Parameters**

**choice** (*int*) – Integer number representing an option key.

**\_EditCheckConnect**() → None

Displays the POD devices parameters, asks the user to edit the device, and then reconnects the device for each POD device type.

**\_EditSaveFilePath()** → None

Asks the user for the POD device type, then asks the user for a new file name and path, then sets the value to the POD devices.

**\_GetChosenDeviceType**(*question: str*) → str

Asks the user which type of POD device they want.

**Parameters**

**question** (*str*) – Question to ask the user.

**Returns**

String of the user input (may be invalid POD device type).

**Return type**

str

**\_GetParams**(*podParametersDict: None | dict[str, None]*) → dict[str, dict | None]

If no parameters are give, this asks user which types of POD devices they want to use. Then it checks if the parameters are valid.

**Parameters**

**podParametersDict** (*None | dict[str, None]*) – Dictionary of all POD devices initialization. The keys are the device name and the entries are the initialization dictionaries.

**Returns**

Dictionary whose keys are the POD device name, and value the setup dictionary.

**Return type**

dict[str, dict | None]

**\_PrintInitCode()** → None

Prints code that can be used to initialize and run SetupPodDevices with the current parameters.

**\_PrintOptions()** → None

Prints options available for user.

**\_Reconnect()** → bool

Reconnects all POD devices.

**Returns**

Bool that is true if all devices were successfully connected. False otherwise.

**Return type**

bool

**\_RemoveDevice()** → None

Displays the POD devices parameters, asks the user which device ro remove, and then deletes that POD device.

**\_Set\_Setup\_PodDevices**(*podParametersDict: dict[str, dict | None]*) → None

Sets the \_Setup\_PodDevices varaible to have keys as the POD device name and values as the setup class.

**Parameters**

**podParametersDict** (*dict[str, dict | None]*) – Dictionary with keys as the device names and values as None or the respective parameter dictionary.

**\_ShowCurrentSettings()** → None

Displays the POD device settings for all devices, and then prints the save file name.



**\_Stream()** → float

Streams data from all POD devices and prints the execution time.

**Returns**

Float of the execution time in seconds.

**Return type**

float

**\_StreamAllDevices()** → None

Streams data from all the devices. User is asked to click enter to stop streaming. Data is saved to file. Uses threading.

**static \_TimeFunc(func: function)** → float

Runs a function and gets the calculated execution time.

**Parameters**

**func** (*function*) – Function/method name.

**Returns**

Float of the execution time in seconds rounded to 3 decimal places.

**Return type**

float

## 1.2.4 Module contents



## INDICES AND TABLES

- `genindex`
- `modindex`
- `search`



## PYTHON MODULE INDEX

### p

- PodApi, 51
- PodApi.Commands, 5
- PodApi.Commands.PodCommands, 1
- PodApi.Devices, 26
- PodApi.Devices.BasicPodProtocol, 9
- PodApi.Devices.PodDevice\_8206HR, 13
- PodApi.Devices.PodDevice\_8229, 14
- PodApi.Devices.PodDevice\_8401HR, 19
- PodApi.Devices.PodDevice\_8480SC, 23
- PodApi.Devices.SerialPorts, 9
- PodApi.Devices.SerialPorts.PortAccess, 5
- PodApi.Devices.SerialPorts.SerialComm, 6
- PodApi.Packets, 44
- PodApi.Packets.Binary, 26
- PodApi.Packets.Binary4, 28
- PodApi.Packets.Binary5, 31
- PodApi.Packets.Packet, 37
- PodApi.Packets.Standard, 41
- PodApi.Parameters, 51
- PodApi.Parameters.Params8206HR, 44
- PodApi.Parameters.Params8229, 45
- PodApi.Parameters.Params8401HR, 47
- PodApi.Parameters.Params8480SC, 49
- PodApi.Parameters.ParamsBasic, 51

### S

- Setup, 77
- Setup.Inputs, 57
- Setup.Inputs.GetUserInput, 52
- Setup.Setup\_PodDevices, 74
- Setup.SetupOneDevice, 74
- Setup.SetupOneDevice.Setup\_8206HR, 57
- Setup.SetupOneDevice.Setup\_8229, 59
- Setup.SetupOneDevice.Setup\_8401HR, 61
- Setup.SetupOneDevice.Setup\_8480SC, 65
- Setup.SetupOneDevice.Setup\_PodInterface, 68



# INDEX

## Symbols

|   |   |
|---|---|
| <code>_AskOption()</code> (Setup.Setup_PodDevices.Setup_PodDevices method), 75                          | <code>_ChooseWidth()</code> (Setup.SetupOneDevice.Setup_8480SC.Setup_8480SC static method), 66        |
| <code>_AskToStopStream()</code> (Setup.Setup_PodDevices.Setup_PodDevices method), 75                    | <code>_CodeDCmode()</code> (Setup.SetupOneDevice.Setup_8401HR.Setup_8401HR static method), 61         |
| <code>_AskUserForDevices()</code> (Setup.Setup_PodDevices.Setup_PodDevices static method), 75           | <code>_CodeDecimalAsHex()</code> (PodApi.Devices.PodDevice_8229.Pod8229 static method), 16            |
| <code>_BuildFileName()</code> (Setup.SetupOneDevice.Setup_PodInterface.Setup_Interface method), 70      | <code>_CodeHighpass()</code> (Setup.SetupOneDevice.Setup_8401HR.Setup_8401HR static method), 62       |
| <code>_CheckForValidParams()</code> (Setup.Setup_PodDevices.Setup_PodDevices method), 75                | <code>_ConnectNewDevice()</code> (Setup.Setup_PodDevices.Setup_PodDevices method), 75                 |
| <code>_CheckParams()</code> (PodApi.Parameters.Params8206HR.Params8206HR method), 45                    | <code>_ConnectPODdevice()</code> (Setup.SetupOneDevice.Setup_8206HR.Setup_8206HR method), 57          |
| <code>_CheckParams()</code> (PodApi.Parameters.Params8229.Params8229 method), 46                        | <code>_ConnectPODdevice()</code> (Setup.SetupOneDevice.Setup_8229.Setup_8229 method), 59              |
| <code>_CheckParams()</code> (PodApi.Parameters.Params8401HR.Params8401HR method), 48                    | <code>_ConnectPODdevice()</code> (Setup.SetupOneDevice.Setup_8401HR.Setup_8401HR method), 62          |
| <code>_CheckParams()</code> (PodApi.Parameters.Params8480SC.Params8480SC method), 50                    | <code>_ConnectPODdevice()</code> (Setup.SetupOneDevice.Setup_8480SC.Setup_8480SC method), 67          |
| <code>_CheckParams()</code> (PodApi.Parameters.ParamsBasic.Params method), 51                           | <code>_ConnectPODdevice()</code> (Setup.SetupOneDevice.Setup_PodInterface.Setup_Interface method), 71 |
| <code>_ChoosePeriod()</code> (Setup.SetupOneDevice.Setup_8480SC.Setup_8480SC static method), 66         | <code>_Custom108GETTTLSETUP()</code> (PodApi.Devices.PodDevice_8480SC.Pod8480SC static method), 25    |
| <code>_ChoosePort()</code> (Setup.SetupOneDevice.Setup_PodInterface.Setup_Interface static method), 70  | <code>_Custom109SETTTLSETUP()</code> (PodApi.Devices.PodDevice_8480SC.Pod8480SC static method), 26    |
| <code>_ChoosePortLinux()</code> (PodApi.Devices.SerialPorts.PortAccess.FindPorts static method), 6      | <code>_Custom140SETTIME()</code> (PodApi.Devices.PodDevice_8229.Pod8229 static method), 17            |
| <code>_ChoosePortWindows()</code> (PodApi.Devices.SerialPorts.PortAccess.FindPorts static method), 6    | <code>_CustomSTIMULUS()</code> (PodApi.Devices.PodDevice_8480SC.Pod8480SC static method), 26          |
| <code>_ChoosePreampGain()</code> (Setup.SetupOneDevice.Setup_8206HR.Setup_8206HR static method), 57     | <code>_CustomSYNCCONFIG()</code> (PodApi.Devices.PodDevice_8480SC.Pod8480SC static method), 26        |
| <code>_ChooseRepeat()</code> (Setup.SetupOneDevice.Setup_8480SC.Setup_8480SC static method), 66         |   |
| <code>_ChooseStimulusConfig()</code> (Setup.SetupOneDevice.Setup_8480SC.Setup_8480SC static method), 66 |   |
| <code>_ChooseSyncConfig()</code>  |   |

|   |   |
|---|---|
| <code>dApi.Devices.PodDevice_8480SC.Pod8480SC</code><br><code>static method</code> ), 26  | <code>(Setup.SetupOneDevice.Setup_8401HR.Setup_8401HR</code><br><code>method</code> ), 62   |
| <code>_DecodeDecimalAsHex()</code> (Po-<br><code>dApi.Devices.PodDevice_8229.Pod8229</code><br><code>static method</code> ), 17                     | <code>_GetParam_onePODdevice()</code><br><code>(Setup.SetupOneDevice.Setup_8480SC.Setup_8480SC</code><br><code>method</code> ), 67                |
| <code>_DisconnectAllPODdevices()</code><br><code>(Setup.SetupOneDevice.Setup_PodInterface.Setup_Interface</code><br><code>method</code> ), 71       | <code>_GetParam_onePODdevice()</code><br><code>(Setup.SetupOneDevice.Setup_PodInterface.Setup_Interface</code><br><code>method</code> ), 71       |
| <code>_DoOption()</code> (Setup.Setup_PodDevices.Setup_PodDevices<br><code>method</code> ), 75  | <code>_GetParams()</code> (Setup.Setup_PodDevices.Setup_PodDevices<br><code>method</code> ), 76   |
| <code>_EditCheckConnect()</code><br><code>(Setup.Setup_PodDevices.Setup_PodDevices</code><br><code>method</code> ), 75                              | <code>_GetPreamplifierDeviceName()</code><br><code>(Setup.SetupOneDevice.Setup_8401HR.Setup_8401HR</code><br><code>method</code> ), 62            |
| <code>_EditParams()</code> (Setup.SetupOneDevice.Setup_PodInterface.Setup_PodInterface<br><code>method</code> ), 71                                 | <code>_GetScheduleForWeek()</code><br><code>(Setup.SetupOneDevice.Setup_8229.Setup_8229</code><br><code>static method</code> ), 60                |
| <code>_EditSaveFilePath()</code><br><code>(Setup.Setup_PodDevices.Setup_PodDevices</code><br><code>method</code> ), 75                              | <code>_GetTimeHeader_forTXT()</code><br><code>(Setup.SetupOneDevice.Setup_PodInterface.Setup_Interface</code><br><code>static method</code> ), 72 |
| <code>_FixABCDtype()</code> (PodApi.Devices.PodDevice_8401HR.Pod8401HR<br><code>static method</code> ), 22  | <code>_NiceABCDtableText()</code><br><code>(Setup.SetupOneDevice.Setup_8401HR.Setup_8401HR</code><br><code>method</code> ), 62                    |
| <code>_FixTypeInTuple()</code> (Po-<br><code>dApi.Parameters.ParamsBasic.Params</code> <code>static</code><br><code>method</code> ), 51             | <code>_OpenSaveFile()</code> (Setup.SetupOneDevice.Setup_PodInterface.Setup_PodInterface<br><code>method</code> ), 72                             |
| <code>_GetChosenDeviceType()</code><br><code>(Setup.Setup_PodDevices.Setup_PodDevices</code><br><code>method</code> ), 76                           | <code>_OpenSaveFile_EDF()</code><br><code>(Setup.SetupOneDevice.Setup_8206HR.Setup_8206HR</code><br><code>method</code> ), 58                     |
| <code>_GetForbiddenNames()</code><br><code>(Setup.SetupOneDevice.Setup_PodInterface.Setup_PodInterface</code><br><code>method</code> ), 71          | <code>_OpenSaveFile_EDF()</code><br><code>(Setup.SetupOneDevice.Setup_8229.Setup_8229</code><br><code>method</code> ), 60                         |
| <code>_GetPODdeviceParameterTable()</code><br><code>(Setup.SetupOneDevice.Setup_8206HR.Setup_8206HR</code><br><code>method</code> ), 57             | <code>_OpenSaveFile_EDF()</code><br><code>(Setup.SetupOneDevice.Setup_8401HR.Setup_8401HR</code><br><code>method</code> ), 63                     |
| <code>_GetPODdeviceParameterTable()</code><br><code>(Setup.SetupOneDevice.Setup_8229.Setup_8229</code><br><code>method</code> ), 60                 | <code>_OpenSaveFile_EDF()</code><br><code>(Setup.SetupOneDevice.Setup_8480SC.Setup_8480SC</code><br><code>method</code> ), 67                     |
| <code>_GetPODdeviceParameterTable()</code><br><code>(Setup.SetupOneDevice.Setup_8401HR.Setup_8401HR</code><br><code>method</code> ), 62             | <code>_OpenSaveFile_EDF()</code><br><code>(Setup.SetupOneDevice.Setup_PodInterface.Setup_PodInterface</code><br><code>method</code> ), 72         |
| <code>_GetPODdeviceParameterTable()</code><br><code>(Setup.SetupOneDevice.Setup_8480SC.Setup_8480SC</code><br><code>method</code> ), 67             | <code>_OpenSaveFile_TXT()</code><br><code>(Setup.SetupOneDevice.Setup_8206HR.Setup_8206HR</code><br><code>method</code> ), 58                     |
| <code>_GetPODdeviceParameterTable()</code><br><code>(Setup.SetupOneDevice.Setup_PodInterface.Setup_PodInterface</code><br><code>method</code> ), 71 | <code>_OpenSaveFile_TXT()</code><br><code>(Setup.SetupOneDevice.Setup_8229.Setup_8229</code><br><code>method</code> ), 60                         |
| <code>_GetParam_allPODdevices()</code><br><code>(Setup.SetupOneDevice.Setup_PodInterface.Setup_PodInterface</code><br><code>method</code> ), 71     | <code>_OpenSaveFile_TXT()</code><br><code>(Setup.SetupOneDevice.Setup_8401HR.Setup_8401HR</code><br><code>method</code> ), 63                     |
| <code>_GetParam_onePODdevice()</code><br><code>(Setup.SetupOneDevice.Setup_8206HR.Setup_8206HR</code><br><code>method</code> ), 58                  | <code>_OpenSaveFile_TXT()</code><br><code>(Setup.SetupOneDevice.Setup_8480SC.Setup_8480SC</code><br><code>method</code> ), 67                     |
| <code>_GetParam_onePODdevice()</code><br><code>(Setup.SetupOneDevice.Setup_8229.Setup_8229</code><br><code>method</code> ), 60                      | <code>_OpenSaveFile_TXT()</code><br><code>(Setup.SetupOneDevice.Setup_PodInterface.Setup_PodInterface</code><br><code>method</code> ), 72         |
| <code>_GetParam_onePODdevice()</code>   |   |



---

```

_PHYSICAL_BOUND_uV (Setup.SetupOneDevice.Setup_8206HR.Setup_8206HR attribute), 64
attribute), 58
_Set_Setup_PodDevices()
_PHYSICAL_BOUND_uV (Setup.SetupOneDevice.Setup_8401HR.Setup_8401HR attribute), 63
method), 76
_PrintDeviceNumber()
_Setup_PodDevices (Setup.Setup_PodDevices.Setup_PodDevices
(Setup.SetupOneDevice.Setup_PodInterface.Setup_Interface attribute), 74
static method), 72
_ShowCurrentSettings()
_PrintInitCode() (Setup.Setup_PodDevices.Setup_PodDevices (Setup.Setup_PodDevices.Setup_PodDevices
method), 76 method), 76
_PrintOptions() (Setup.Setup_PodDevices.Setup_PodDevices (Setup.Setup_PodDevices.Setup_PodDevices
method), 76 method), 76
_ReadPODpacket_Recursive() (Po- _StreamAllDevices()
dApi.Devices.BasicPodProtocol.Pod method), (Setup.Setup_PodDevices.Setup_PodDevices
11 method), 77
_Read_Binary() (PodApi.Devices.BasicPodProtocol.Pod _StreamThreading() (Setup.SetupOneDevice.Setup_8206HR.Setup_8206
method), 11 method), 58
_Read_Binary() (PodApi.Devices.PodDevice_8206HR.Pod8206HR _StreamThreading() (Setup.SetupOneDevice.Setup_8229.Setup_8229
method), 13 method), 61
_Read_Binary() (PodApi.Devices.PodDevice_8401HR.Pod8401HR _StreamThreading() (Setup.SetupOneDevice.Setup_8401HR.Setup_8401
method), 22 method), 65
_Read_GetCommand() (Po- _StreamThreading() (Setup.SetupOneDevice.Setup_8480SC.Setup_8480
dApi.Devices.BasicPodProtocol.Pod method), method), 68
12 _StreamThreading() (Setup.SetupOneDevice.Setup_PodInterface.Setup_
method), 73
_Read_Standard() (Po- _StreamUntilStop() (Setup.SetupOneDevice.Setup_8206HR.Setup_8206
dApi.Devices.BasicPodProtocol.Pod method), method), 58
12 _StreamUntilStop() (Setup.SetupOneDevice.Setup_8229.Setup_8229
method), 61
_Read_ToETX() (PodApi.Devices.BasicPodProtocol.Pod _StreamUntilStop() (Setup.SetupOneDevice.Setup_8401HR.Setup_8401
method), 12 method), 65
_Reconnect() (Setup.Setup_PodDevices.Setup_PodDevices _StreamUntilStop() (Setup.SetupOneDevice.Setup_8480SC.Setup_8480
method), 76 method), 68
_RemoveDevice() (Setup.SetupOneDevice.Setup_PodInterface.Setup_PodInterface _StreamUntilStop() (Setup.SetupOneDevice.Setup_8480SC.Setup_8480
method), 72 method), 68
_RemoveDevice() (Setup.Setup_PodDevices.Setup_PodDevices _TestDeviceConnection()
method), 76 (Setup.SetupOneDevice.Setup_PodInterface.Setup_Interface
static method), 73
_SelectDeviceFromDict()
(Setup.SetupOneDevice.Setup_PodInterface.Setup_PodInterface _TestDeviceConnection_All()
method), 72 (Setup.SetupOneDevice.Setup_PodInterface.Setup_Interface
static method), 73
_SetBias() (Setup.SetupOneDevice.Setup_8401HR.Setup_8401HR method), 73
static method), 63
_TimeFunc() (Setup.Setup_PodDevices.Setup_PodDevices
static method), 77
_SetDCMode() (Setup.SetupOneDevice.Setup_8401HR.Setup_8401HR static method), 63
_TranslateTTLbyte_ASCII() (Po-
dApi.Devices.PodDevice_8206HR.Pod8206HR
static method), 14
_SetForMappedChannels()
(Setup.SetupOneDevice.Setup_8401HR.Setup_8401HR static method), 64
_TtlSetup() (Setup.SetupOneDevice.Setup_8480SC.Setup_8480SC
static method), 68
_SetHighpass() (Setup.SetupOneDevice.Setup_8401HR.Setup_8401HR static method), 64
_ValidateChecksum() (Po-
static method), 64 dApi.Devices.BasicPodProtocol.Pod static
method), 13
_SetLowpass() (Setup.SetupOneDevice.Setup_8401HR.Setup_8401HR static method), 23
_SetNumberOfDevices()
_ValidatePreampGain() (Po-
(Setup.SetupOneDevice.Setup_PodInterface.Setup_Interface dApi.Devices.PodDevice_8401HR.Pod8401HR
static method), 73 static method), 23
_SetPreampGain() (Setup.SetupOneDevice.Setup_8401HR.Setup_8401HR static method), 64
_SetSSGain() (Setup.SetupOneDevice.Setup_8401HR.Setup_8401HR static method), 23

```

---

|  |      |   |
|--|------|---|
| <code>_Validate_Day()</code><br><i>dApi.Devices.PodDevice_8229.Pod8229 static method), 17</i>                    | (Po- | <code>__U16</code> ( <i>PodApi.Commands.PodCommands.CommandSet attribute), 5</i>                                |
| <code>_Validate_Hours()</code><br><i>dApi.Devices.PodDevice_8229.Pod8229 static method), 17</i>                  | (Po- | <code>__U32</code> ( <i>PodApi.Commands.PodCommands.CommandSet attribute), 5</i>                                |
| <code>_Validate_Schedule()</code><br><i>dApi.Devices.PodDevice_8229.Pod8229 static method), 18</i>               | (Po- | <code>__U8</code> ( <i>PodApi.Commands.PodCommands.CommandSet attribute), 5</i>                                 |
| <code>_Validate_Speed()</code><br><i>dApi.Devices.PodDevice_8229.Pod8229 static method), 18</i>                  | (Po- | <code>__commands</code> ( <i>PodApi.Commands.PodCommands.CommandSet attribute), 1</i>                           |
| <code>_Voltage_PrimaryChannels()</code><br><i>dApi.Packets.Binary5.PacketBinary5 method), 36</i>                 | (Po- | <code>__serialInst</code> ( <i>PodApi.Devices.SerialPorts.SerialComm.PortIO attribute), 6</i>                   |
| <code>_Voltage_PrimaryChannels_Biosensor()</code><br><i>dApi.Packets.Binary5.PacketBinary5 method), 36</i>       | (Po- | <code>_commands</code> ( <i>PodApi.Devices.BasicPodProtocol.Pod attribute), 9</i>                               |
| <code>_Voltage_PrimaryChannels_EEGEMG()</code><br><i>dApi.Packets.Binary5.PacketBinary5 method), 36</i>          | (Po- | <code>_commands</code> ( <i>PodApi.Packets.Packet.Packet attribute), 37</i>                                     |
| <code>_Voltage_SecondaryChannels()</code><br><i>dApi.Packets.Binary5.PacketBinary5 method), 36</i>               | (Po- | <code>_customPayload</code> ( <i>PodApi.Packets.Standard.PacketStandard attribute), 41</i>                      |
| <code>_WriteDataToFile_EDF()</code><br><i>(Setup.SetupOneDevice.Setup_8206HR.Setup_8206HR static method), 59</i> | (Po- | <code>_customPayloadArgs</code> ( <i>PodApi.Packets.Standard.PacketStandard attribute), 42</i>                  |
| <code>_WriteDataToFile_EDF()</code><br><i>(Setup.SetupOneDevice.Setup_8401HR.Setup_8401HR static method), 65</i> | (Po- | <code>_options</code> ( <i>Setup.Setup_PodDevices.Setup_PodDevices attribute), 74</i>                           |
| <code>_WriteDataToFile_TXT()</code><br><i>(Setup.SetupOneDevice.Setup_8206HR.Setup_8206HR static method), 59</i> | (Po- | <code>_podDevices</code> ( <i>Setup.SetupOneDevice.Setup_PodInterface.Setup_Interface attribute), 68</i>        |
| <code>_WriteDataToFile_TXT()</code><br><i>(Setup.SetupOneDevice.Setup_8401HR.Setup_8401HR static method), 65</i> | (Po- | <code>_podParametersDict</code> ( <i>Setup.SetupOneDevice.Setup_PodInterface.Setup_Interface attribute), 68</i> |
| <code>__ARGUMENTS</code> ( <i>PodApi.Commands.PodCommands.CommandSet attribute), 4</i>                           |      | <code>_port</code> ( <i>PodApi.Devices.BasicPodProtocol.Pod attribute), 9</i>                                   |
| <code>__BINARY</code> ( <i>PodApi.Commands.PodCommands.CommandSet attribute), 4</i>                              |      | <code>_preampGain</code> ( <i>PodApi.Devices.PodDevice_8206HR.Pod8206HR attribute), 13</i>                      |
| <code>__BuildPortName()</code><br><i>dApi.Devices.SerialPorts.SerialComm.PortIO method), 8</i>                   |      | <code>_preampGain</code> ( <i>PodApi.Devices.PodDevice_8401HR.Pod8401HR attribute), 19</i>                      |
| <code>__CHANNELMAPALL</code><br><i>dApi.Devices.PodDevice_8401HR.Pod8401HR attribute), 23</i>                    |      | <code>_preampGain</code> ( <i>PodApi.Packets.Binary4.PacketBinary4 attribute), 28</i>                           |
| <code>__DESCRIPTION</code> ( <i>PodApi.Commands.PodCommands.CommandSet attribute), 4</i>                         |      | <code>_preampGain</code> ( <i>PodApi.Packets.Binary5.PacketBinary5 attribute), 32</i>                           |
| <code>__NAME</code> ( <i>PodApi.Commands.PodCommands.CommandSet attribute), 4</i>                                |      | <code>_saveFileName</code> ( <i>Setup.SetupOneDevice.Setup_PodInterface.Setup_Interface attribute), 68</i>      |
| <code>__NOVALUE</code> ( <i>PodApi.Commands.PodCommands.CommandSet attribute), 5</i>                             |      | <code>_ssGain</code> ( <i>PodApi.Devices.PodDevice_8401HR.Pod8401HR attribute), 19</i>                          |
| <code>__RETURNS</code> ( <i>PodApi.Commands.PodCommands.CommandSet attribute), 5</i>                             |      | <code>_ssGain</code> ( <i>PodApi.Packets.Binary5.PacketBinary5 attribute), 31</i>                               |
|  |      | <code>_streamMode</code> ( <i>Setup.SetupOneDevice.Setup_8229.Setup_8229 attribute), 59</i>                     |
|  |      | <code>_streamMode</code> ( <i>Setup.SetupOneDevice.Setup_8480SC.Setup_8480SC attribute), 65</i>                 |
|  |      | <code>_uV()</code> ( <i>Setup.SetupOneDevice.Setup_PodInterface.Setup_Interface static method), 73</i>          |

## A

|  |
|--|
| <code>AddCommand()</code> ( <i>PodApi.Commands.PodCommands.CommandSet method), 1</i>                     |
| <code>AddPODdevice()</code> ( <i>Setup.SetupOneDevice.Setup_PodInterface.Setup_Interface method), 68</i> |
| <code>GetEXT0</code> ( <i>PodApi.Packets.Binary5.PacketBinary5 attribute), 32</i>                        |

- aEXT1 (*PodApi.Packets.Binary5.PacketBinary5* attribute), 32
- AnalogEXT() (*PodApi.Packets.Binary5.PacketBinary5* method), 32
- AnalogTTL() (*PodApi.Packets.Binary5.PacketBinary5* method), 33
- AreDeviceParamsValid() (*Setup.SetupOneDevice.Setup\_PodInterface.Setup\_Interface* method), 69
- ArgumentHexChar() (*PodApi.Commands.PodCommands.CommandSet* method), 2
- AsciiBytesToInt() (*PodApi.Packets.Packet.Packet* static method), 37
- ASCIIbytesToInt\_Split() (*PodApi.Packets.Packet.Packet* static method), 37
- AskForBool() (*Setup.Inputs.GetUserInput.UserInput* static method), 52
- AskForFloat() (*Setup.Inputs.GetUserInput.UserInput* static method), 52
- AskForFloatInList() (*Setup.Inputs.GetUserInput.UserInput* static method), 52
- AskForFloatInRange() (*Setup.Inputs.GetUserInput.UserInput* static method), 52
- AskForInput() (*Setup.Inputs.GetUserInput.UserInput* static method), 53
- AskForInt() (*Setup.Inputs.GetUserInput.UserInput* static method), 53
- AskForIntInList() (*Setup.Inputs.GetUserInput.UserInput* static method), 53
- AskForIntInRange() (*Setup.Inputs.GetUserInput.UserInput* static method), 53
- AskForStrInList() (*Setup.Inputs.GetUserInput.UserInput* static method), 54
- AskForType() (*Setup.Inputs.GetUserInput.UserInput* static method), 54
- AskForTypeInList() (*Setup.Inputs.GetUserInput.UserInput* static method), 54
- AskForTypeInRange() (*Setup.Inputs.GetUserInput.UserInput* static method), 55
- AskYN() (*Setup.Inputs.GetUserInput.UserInput* static method), 55
- aTTL1 (*PodApi.Packets.Binary5.PacketBinary5* attribute), 32
- aTTL2 (*PodApi.Packets.Binary5.PacketBinary5* attribute), 32
- aTTL3 (*PodApi.Packets.Binary5.PacketBinary5* attribute), 32
- aTTL4 (*PodApi.Packets.Binary5.PacketBinary5* attribute), 32
- ## B
- bias (*PodApi.Parameters.Params8401HR.Params8401HR* attribute), 48
- BinaryBytesToInt() (*PodApi.Packets.Packet.Packet* static method), 38
- BinaryBytesToInt\_Split() (*PodApi.Packets.Packet.Packet* static method), 38
- BinaryBytesToVoltage() (*PodApi.Packets.Binary4.PacketBinary4* static method), 29
- binaryData (*PodApi.Packets.Binary.PacketBinary* attribute), 27
- binaryLength (*PodApi.Packets.Binary.PacketBinary* attribute), 26
- BinaryLength() (*PodApi.Packets.Binary.PacketBinary* method), 27
- BuildEmptySchedule() (*PodApi.Parameters.Params8229.Params8229* static method), 46
- BuildPODpacket\_Standard() (*PodApi.Packets.Packet.Packet* static method), 38
- BuildSetDayScheduleArgument() (*PodApi.Devices.PodDevice\_8229.Pod8229* static method), 14
- ## C
- CalculateBiasDAC\_GetDACValue() (*PodApi.Devices.PodDevice\_8401HR.Pod8401HR* static method), 19
- CalculateBiasDAC\_GetVout() (*PodApi.Devices.PodDevice\_8401HR.Pod8401HR* static method), 19
- CastFloat() (*Setup.Inputs.GetUserInput.UserInput* static method), 55
- CastInt() (*Setup.Inputs.GetUserInput.UserInput* static method), 55
- CastStr() (*Setup.Inputs.GetUserInput.UserInput* static method), 56
- Ch() (*PodApi.Packets.Binary4.PacketBinary4* method), 29
- ch0 (*PodApi.Packets.Binary4.PacketBinary4* attribute), 28
- ch1 (*PodApi.Packets.Binary4.PacketBinary4* attribute), 28
- ch2 (*PodApi.Packets.Binary4.PacketBinary4* attribute), 29
- Channel() (*PodApi.Packets.Binary5.PacketBinary5* method), 33
- channelLabels (*PodApi.Parameters.Params8401HR.Params8401HR* attribute), 49
- channels (*PodApi.Packets.Binary5.PacketBinary5* attribute), 32

|                         |   |                              |   |
|-------------------------|---|------------------------------|---|
| CheckFileExt()          | (Setup.Inputs.GetUserInput.UserInput static method), 56                     | DecodeDayOfWeek()            | (PodApi.Devices.PodDevice_8229.Pod8229 static method), 15                   |
| CheckIfPacketIsValid()  | (PodApi.Packets.Binary.PacketBinary static method), 27                      | DecodeDaySchedule()          | (PodApi.Devices.PodDevice_8229.Pod8229 static method), 15                   |
| CheckIfPacketIsValid()  | (PodApi.Packets.Binary4.PacketBinary4 static method), 29                    | DecodeLcdSchedule()          | (PodApi.Devices.PodDevice_8229.Pod8229 static method), 16                   |
| CheckIfPacketIsValid()  | (PodApi.Packets.Binary5.PacketBinary5 static method), 33                    | DecodeSSConfigBitmask()      | (PodApi.Devices.PodDevice_8401HR.Pod8401HR static method), 19               |
| CheckIfPacketIsValid()  | (PodApi.Packets.Packet.Packet static method), 38                            | DecodeStimulusConfigBits()   | (PodApi.Devices.PodDevice_8480SC.Pod8480SC static method), 23               |
| CheckIfPacketIsValid()  | (PodApi.Packets.Standard.PacketStandard static method), 42                  | DecodeSyncConfigBits()       | (PodApi.Devices.PodDevice_8480SC.Pod8480SC static method), 23               |
| Checksum()              | (PodApi.Packets.Packet.Packet static method), 39                            | DecodeTTLByte()              | (PodApi.Devices.PodDevice_8401HR.Pod8401HR static method), 20               |
| ChoosePort()            | (PodApi.Devices.BasicPodProtocol.Pod static method), 9                      | DecodeTTLConfigBits()        | (PodApi.Devices.PodDevice_8480SC.Pod8480SC static method), 24               |
| ChoosePort()            | (PodApi.Devices.SerialPorts.PortAccess.Firmware static method), 5           | DecodeTTLPayload()           | (PodApi.Devices.PodDevice_8401HR.Pod8401HR static method), 20               |
| CloseSerialPort()       | (PodApi.Devices.SerialPorts.SerialComm.PortIO static method), 6             | DefaultPayload()             | (PodApi.Packets.Standard.PacketStandard static method), 42                  |
| CodeDayOfWeek()         | (PodApi.Devices.PodDevice_8229.Pod8229 static method), 14                   | Description()                | (PodApi.Commands.PodCommands.CommandSet static method), 2                   |
| CodeDaySchedule()       | (PodApi.Devices.PodDevice_8229.Pod8229 static method), 15                   | DisplayPODdeviceParameters() | (Setup.SetupOneDevice.Setup_PodInterface.Setup_Interface static method), 69 |
| commandNumber           | (PodApi.Packets.Packet.Packet attribute), 37                                | DoesCommandExist()           | (PodApi.Commands.PodCommands.CommandSet static method), 2                   |
| CommandNumber()         | (PodApi.Packets.Packet.Packet static method), 39                            | <b>E</b>                     |   |
| CommandNumberFromName() | (PodApi.Commands.PodCommands.CommandSet static method), 2                   | EEG1()                       | (PodApi.Parameters.Params8206HR.Params8206HR static method), 44             |
| CommandSet              | (class in PodApi.Commands.PodCommands), 1                                   | EEG2()                       | (PodApi.Parameters.Params8206HR.Params8206HR static method), 44             |
| ConnectAllPODdevices()  | (Setup.SetupOneDevice.Setup_PodInterface.Setup_Interface static method), 69 | EEG3_EMG()                   | (PodApi.Parameters.Params8206HR.Params8206HR static method), 44             |
| <b>D</b>                |   | estimCurrent                 | (PodApi.Parameters.Params8480SC.Params8480SC attribute), 49                 |
| dcMode                  | (PodApi.Parameters.Params8401HR.Params8401HR attribute), 48                 | estimCurrent_CH0()           | (PodApi.Parameters.Params8480SC.Params8480SC static method), 50             |
| DecodeChannelBitmask()  | (PodApi.Devices.PodDevice_8401HR.Pod8401HR static method), 19               | estimCurrent_CH1()           | (PodApi.Parameters.Params8480SC.Params8480SC static method), 50             |
| DecodeDayAndSchedule()  | (PodApi.Devices.PodDevice_8229.Pod8229 static method), 15                   |                              |   |



ETX() (*PodApi.Packets.Packet.Packet* static method), 39

## F

FindPorts (class in *PodApi.Devices.SerialPorts.PortAccess*), 5

Flush() (*PodApi.Devices.SerialPorts.SerialComm.PortIO* method), 6

FlushPort() (*PodApi.Devices.BasicPodProtocol.Pod* method), 9

## G

GetAllPortNames() (*PodApi.Devices.SerialPorts.PortAccess.FindPorts* static method), 5

GetAnalogEXT() (*PodApi.Packets.Binary5.PacketBinary5* static method), 33

GetAnalogTTL() (*PodApi.Packets.Binary5.PacketBinary5* static method), 34

GetBasicCommands() (*PodApi.Commands.PodCommands.CommandSet* static method), 2

GetBinaryData() (*PodApi.Packets.Binary.PacketBinary* method), 27

GetBinaryLength() (*PodApi.Packets.Binary.PacketBinary* static method), 27

GetBinaryLength() (*PodApi.Packets.Binary4.PacketBinary4* static method), 29

GetBinaryLength() (*PodApi.Packets.Binary5.PacketBinary5* static method), 34

GetCh() (*PodApi.Packets.Binary4.PacketBinary4* static method), 29

GetChannelBitmask() (*PodApi.Devices.PodDevice\_8401HR.Pod8401HR* static method), 20

GetChannelMapForPreampDevice() (*PodApi.Devices.PodDevice\_8401HR.Pod8401HR* static method), 20

GetChannels() (*PodApi.Packets.Binary5.PacketBinary5* static method), 34

GetCommandNumber() (*PodApi.Packets.Packet.Packet* static method), 39

GetCommands() (*PodApi.Commands.PodCommands.CommandSet* method), 3

GetCurrentTime() (*PodApi.Devices.PodDevice\_8229.Pod8229* static method), 16

GetDeviceCommands() (*PodApi.Devices.BasicPodProtocol.Pod* method), 9

GetDeviceName() (*Setup.SetupOneDevice.Setup\_8206HR.Setup\_8206HR* static method), 57

GetDeviceName() (*Setup.SetupOneDevice.Setup\_8229.Setup\_8229* static method), 59

GetDeviceName() (*Setup.SetupOneDevice.Setup\_8401HR.Setup\_8401HR* static method), 61

GetDeviceName() (*Setup.SetupOneDevice.Setup\_8480SC.Setup\_8480SC* static method), 65

GetDeviceName() (*Setup.SetupOneDevice.Setup\_PodInterface.Setup\_Interface* static method), 69

GetFileName() (*Setup.Inputs.GetUserInput.UserInput* static method), 56

GetFilePath() (*Setup.Inputs.GetUserInput.UserInput* static method), 56

GetInit() (*PodApi.Parameters.Params8206HR.Params8206HR* method), 45

GetInit() (*PodApi.Parameters.Params8229.Params8229* method), 46

GetInit() (*PodApi.Parameters.Params8401HR.Params8401HR* method), 48

GetInit() (*PodApi.Parameters.Params8480SC.Params8480SC* method), 50

GetInit() (*PodApi.Parameters.ParamsBasic.Params* method), 51

GetMinimumLength() (*PodApi.Packets.Binary.PacketBinary* static method), 27

GetMinimumLength() (*PodApi.Packets.Binary4.PacketBinary4* static method), 30

GetMinimumLength() (*PodApi.Packets.Binary5.PacketBinary5* static method), 34

GetMinimumLength() (*PodApi.Packets.Packet.Packet* static method), 39

GetMinimumLength() (*PodApi.Packets.Standard.PacketStandard* static method), 42

GetOptions() (*Setup.Setup\_PodDevices.Setup\_PodDevices* method), 74

GetPacketNumber() (*PodApi.Packets.Binary4.PacketBinary4* static method), 30

GetPacketNumber() (*PodApi.Packets.Binary5.PacketBinary5* static method), 35

GetPayload() (*PodApi.Packets.Standard.PacketStandard* static method), 42

GetPODpacket() (*PodApi.Devices.BasicPodProtocol.Pod* method), 9

GetPODparametersInit() (*Setup.Setup\_PodDevices.Setup\_PodDevices* method), 74

GetPODparametersInit()

(Setup.SetupOneDevice.Setup\_PodInterface.SetupLibCommandBinary() (Pod-  
 method), 69  
 GetPortName() (PodApi.Devices.SerialPorts.SerialComm.PortIO method), 3  
 method), 7  
 IsPreampDeviceSupported() (Po-  
 GetSaveFileName() (Setup.SetupOneDevice.Setup\_PodInterface.SetupPodDevice\_8401HR.Pod8401HR  
 method), 69  
 static method), 21  
 GetSaveFileNames() (Setup.Setup\_PodDevices.Setup\_PodDevice\_8401HR.Pod8401HR  
 method), 74  
 IsSerialClosed() (Po-  
 method), 7  
 GetSelectPortNames() (Pod-  
 dApi.Devices.SerialPorts.PortAccess.FindPorts static method), 6  
 IsSerialOpen() (PodApi.Devices.SerialPorts.SerialComm.PortIO  
 method), 7  
 GetSSConfigBitmask() (Po-  
 dApi.Devices.PodDevice\_8401HR.Pod8401HR static method), 20  
 L  
 GetStatus() (PodApi.Packets.Binary5.PacketBinary5 static method), 35  
 ledCurrent (PodApi.Parameters.Params8480SC.Params8480SC  
 attribute), 49  
 GetSupportedFileExtensions() (Setup.SetupOneDevice.Setup\_8229.Setup\_8229  
 static method), 59  
 ledCurrent\_CH0() (Po-  
 dApi.Parameters.Params8480SC.Params8480SC  
 method), 50  
 GetSupportedFileExtensions() (Setup.SetupOneDevice.Setup\_8480SC.Setup\_8480SC  
 static method), 66  
 ledCurrent\_CH1() (Po-  
 dApi.Parameters.Params8480SC.Params8480SC  
 method), 50  
 GetSupportedFileExtensions() (Setup.SetupOneDevice.Setup\_PodInterface.SetupPodDevice\_8401HR.Pod8401HR  
 static method), 69  
 lowPass (PodApi.Parameters.Params8206HR.Params8206HR  
 attribute), 44  
 GetSupportedPreampDevices() (Pod-  
 dApi.Devices.PodDevice\_8401HR.Pod8401HR static method), 21  
 lowPass (PodApi.Parameters.Params8401HR.Params8401HR  
 attribute), 48  
 GetTTL() (PodApi.Packets.Binary4.PacketBinary4 static  
 method), 30  
 lowPassLabels (PodApi.Parameters.Params8206HR.Params8206HR  
 attribute), 45  
 GetTTLbitmask() (Pod-  
 dApi.Devices.PodDevice\_8401HR.Pod8401HR static method), 21  
 M  
 GetU() (PodApi.Devices.BasicPodProtocol.Pod static  
 method), 10  
 mode (PodApi.Parameters.Params8229.Params8229 at-  
 tribute), 46  
 module  
 H  
 HasCommandNumber() (PodApi.Packets.Packet.Packet  
 method), 40  
 PodApi, 51  
 HasCommands() (PodApi.Packets.Packet.Packet method),  
 40  
 PodApi.Commands, 5  
 HasCustomPayload() (Pod-  
 dApi.Packets.Standard.PacketStandard  
 method), 42  
 PodApi.Commands.PodCommands, 1  
 HasPayload() (PodApi.Packets.Standard.PacketStandard  
 method), 43  
 PodApi.Devices, 26  
 highPass (PodApi.Parameters.Params8401HR.Params8401HR  
 attribute), 48  
 PodApi.Devices.BasicPodProtocol, 9  
 hoursPerDay (PodApi.Parameters.Params8229.Params8229  
 attribute), 47  
 PodApi.Devices.PodDevice\_8206HR, 13  
 PodApi.Devices.PodDevice\_8229, 14  
 PodApi.Devices.PodDevice\_8401HR, 19  
 PodApi.Devices.PodDevice\_8480SC, 23  
 PodApi.Devices.SerialPorts, 9  
 PodApi.Devices.SerialPorts.PortAccess, 5  
 PodApi.Devices.SerialPorts.SerialComm, 6  
 PodApi.Packets, 44  
 PodApi.Packets.Binary, 26  
 PodApi.Packets.Binary4, 28  
 PodApi.Packets.Binary5, 31  
 PodApi.Packets.Packet, 37  
 PodApi.Packets.Standard, 41  
 PodApi.Parameters, 51  
 PodApi.Parameters.Params8206HR, 44  
 PodApi.Parameters.Params8229, 45  
 PodApi.Parameters.Params8401HR, 47  
 I  
 IntToAsciiBytes() (PodApi.Packets.Packet.Packet  
 static method), 40

PodApi.Parameters.Params8480SC, 49  
 PodApi.Parameters.ParamsBasic, 51  
 Setup, 77  
 Setup.Inputs, 57  
 Setup.Inputs.GetUserInput, 52  
 Setup.Setup\_PodDevices, 74  
 Setup.SetupOneDevice, 74  
 Setup.SetupOneDevice.Setup\_8206HR, 57  
 Setup.SetupOneDevice.Setup\_8229, 59  
 Setup.SetupOneDevice.Setup\_8401HR, 61  
 Setup.SetupOneDevice.Setup\_8480SC, 65  
 Setup.SetupOneDevice.Setup\_PodInterface, 68  
 motorDirection (PodApi.Parameters.Params8229.Params8229 attribute), 45  
 motorSpeed (PodApi.Parameters.Params8229.Params8229 attribute), 46  
 muxMode (PodApi.Parameters.Params8401HR.Params8401HR attribute), 47

## N

NoValue() (PodApi.Commands.PodCommands.CommandSet static method), 3

## O

OpenSerialPort() (PodApi.Devices.SerialPorts.SerialComm.PortIO method), 7

## P

Packet (class in PodApi.Packets.Packet), 37  
 PacketBinary (class in PodApi.Packets.Binary), 26  
 PacketBinary4 (class in PodApi.Packets.Binary4), 28  
 PacketBinary5 (class in PodApi.Packets.Binary5), 31  
 packetNumber (PodApi.Packets.Binary4.PacketBinary4 attribute), 28  
 packetNumber (PodApi.Packets.Binary5.PacketBinary5 attribute), 32  
 PacketNumber() (PodApi.Packets.Binary4.PacketBinary4 method), 30  
 PacketNumber() (PodApi.Packets.Binary5.PacketBinary5 method), 35  
 PacketStandard (class in PodApi.Packets.Standard), 41  
 Params (class in PodApi.Parameters.ParamsBasic), 51  
 Params8206HR (class in PodApi.Parameters.Params8206HR), 44  
 Params8229 (class in PodApi.Parameters.Params8229), 45  
 Params8401HR (class in PodApi.Parameters.Params8401HR), 47  
 Params8480SC (class in PodApi.Parameters.Params8480SC), 49  
 payload (PodApi.Packets.Standard.PacketStandard attribute), 42  
 Payload() (PodApi.Packets.Standard.PacketStandard method), 43  
 PayloadToBytes() (PodApi.Packets.Packet.Packet static method), 40  
 Pod (class in PodApi.Devices.BasicPodProtocol), 9  
 Pod8206HR (class in PodApi.Devices.PodDevice\_8206HR), 13  
 Pod8229 (class in PodApi.Devices.PodDevice\_8229), 14  
 Pod8401HR (class in PodApi.Devices.PodDevice\_8401HR), 19  
 Pod8480SC (class in PodApi.Devices.PodDevice\_8480SC), 23  
 PodApi  
   module, 51  
   PodApi.Commands  
     module, 5  
     PodApi.Commands.PodCommands  
       module, 1  
   PodApi.Devices  
     module, 26  
     PodApi.Devices.BasicPodProtocol  
       module, 9  
     PodApi.Devices.PodDevice\_8206HR  
       module, 13  
     PodApi.Devices.PodDevice\_8229  
       module, 14  
     PodApi.Devices.PodDevice\_8401HR  
       module, 19  
     PodApi.Devices.PodDevice\_8480SC  
       module, 23  
     PodApi.Devices.SerialPorts  
       module, 9  
     PodApi.Devices.SerialPorts.PortAccess  
       module, 5  
     PodApi.Devices.SerialPorts.SerialComm  
       module, 6  
   PodApi.Packets  
     module, 44  
     PodApi.Packets.Binary  
       module, 26  
     PodApi.Packets.Binary4  
       module, 28  
     PodApi.Packets.Binary5  
       module, 31  
     PodApi.Packets.Packet  
       module, 37  
     PodApi.Packets.Standard  
       module, 41  
   PodApi.Parameters  
     module, 51  
     PodApi.Parameters.Params8206HR  
       module, 44  
     PodApi.Parameters.Params8229  
       module, 45

PodApi.Parameters.Params8401HR module, 47

PodApi.Parameters.Params8480SC module, 49

PodApi.Parameters.ParamsBasic module, 51

port (PodApi.Parameters.Params8206HR.Params8206HR attribute), 44

port (PodApi.Parameters.Params8229.Params8229 attribute), 45

port (PodApi.Parameters.Params8401HR.Params8401HR attribute), 47

port (PodApi.Parameters.Params8480SC.Params8480SC attribute), 49

port (PodApi.Parameters.ParamsBasic.Params attribute), 51

PortIO (class in PodApi.Devices.SerialPorts.SerialComm), 6

preamp (PodApi.Parameters.Params8480SC.Params8480SC attribute), 49

preampDevice (PodApi.Parameters.Params8401HR.Params8401HR attribute), 47

preampGain (PodApi.Parameters.Params8401HR.Params8401HR attribute), 47

preamplifierGain (PodApi.Parameters.Params8206HR.Params8206HR attribute), 44

PrintSaveFile() (Setup.SetupOneDevice.Setup\_PodInterface.Setup\_Interface method), 70

**R**

randomReverse (PodApi.Parameters.Params8229.Params8229 attribute), 46

rawPacket (PodApi.Packets.Packet.Packet attribute), 37

Read() (PodApi.Devices.SerialPorts.SerialComm.PortIO method), 7

ReadLine() (PodApi.Devices.SerialPorts.SerialComm.PortIO method), 8

ReadPODpacket() (PodApi.Devices.BasicPodProtocol.Pod method), 10

ReadPODpacket() (PodApi.Devices.PodDevice\_8206HR.Pod8206HR method), 13

ReadPODpacket() (PodApi.Devices.PodDevice\_8229.Pod8229 method), 16

ReadPODpacket() (PodApi.Devices.PodDevice\_8401HR.Pod8401HR method), 21

ReadPODpacket() (PodApi.Devices.PodDevice\_8480SC.Pod8480SC method), 24

ReadUntil() (PodApi.Devices.SerialPorts.SerialComm.PortIO method), 8

RemoveCommand() (PodApi.Commands.PodCommands.CommandSet method), 3

RestoreBasicCommands() (PodApi.Commands.PodCommands.CommandSet method), 3

ReturnHexChar() (PodApi.Commands.PodCommands.CommandSet method), 3

reverseBaseTime (PodApi.Parameters.Params8229.Params8229 attribute), 46

reverseVarTime (PodApi.Parameters.Params8229.Params8229 attribute), 46

Run() (Setup.Setup\_PodDevices.Setup\_PodDevices method), 74

**S**

sampleRate (PodApi.Parameters.Params8206HR.Params8206HR attribute), 44

sampleRate (PodApi.Parameters.Params8401HR.Params8401HR attribute), 47

schedule (PodApi.Parameters.Params8229.Params8229 attribute), 46

Search() (PodApi.Commands.PodCommands.CommandSet method), 4

SetBaudrate() (PodApi.Devices.SerialPorts.SerialComm.PortIO method), 8

SetBaudrateOfDevice() (PodApi.Devices.BasicPodProtocol.Pod method), 10

SetCustomPayload() (PodApi.Packets.Standard.PacketStandard method), 43

Setup module, 77

Setup.Inputs module, 57

Setup.Inputs.GetUserInput module, 52

Setup.Setup\_PodDevices module, 74

Setup.SetupOneDevice module, 74

Setup.SetupOneDevice.Setup\_8206HR module, 57

Setup.SetupOneDevice.Setup\_8229 module, 59

Setup.SetupOneDevice.Setup\_8401HR module, 61

Setup.SetupOneDevice.Setup\_8480SC module, 65



Setup.SetupOneDevice.Setup\_PodInterface module, 68  
 Setup\_8206HR (class in Setup.SetupOneDevice.Setup\_8206HR), 57  
 Setup\_8229 (class in Setup.SetupOneDevice.Setup\_8229), 59  
 Setup\_8401HR (class in Setup.SetupOneDevice.Setup\_8401HR), 61  
 Setup\_8480SC (class in Setup.SetupOneDevice.Setup\_8480SC), 65  
 Setup\_Interface (class in Setup.SetupOneDevice.Setup\_PodInterface), 68  
 Setup\_PodDevices (class in Setup.Setup\_PodDevices), 74  
 SetupFileName() (Setup.SetupOneDevice.Setup\_PodInterface.Setup\_Interface method), 70  
 SetupPODparameters() (Setup.Setup\_PodDevices.Setup\_PodDevices method), 74  
 SetupPODparameters() (Setup.SetupOneDevice.Setup\_PodInterface.Setup\_Interface method), 70  
 SetupSaveFile() (Setup.Setup\_PodDevices.Setup\_PodDevices method), 75  
 ssGain (PodApi.Parameters.Params8401HR.Params8401HR attribute), 48  
 status (PodApi.Packets.Binary5.PacketBinary5 attribute), 32  
 Status() (PodApi.Packets.Binary5.PacketBinary5 method), 35  
 stimulus (PodApi.Parameters.Params8480SC.Params8480SC attribute), 49  
 StimulusConfigBits() (PodApi.Devices.PodDevice\_8480SC.Pod8480SC static method), 24  
 StopStream() (Setup.SetupOneDevice.Setup\_8206HR.Setup\_8206HR method), 57  
 StopStream() (Setup.SetupOneDevice.Setup\_8229.Setup\_8229 method), 59  
 StopStream() (Setup.SetupOneDevice.Setup\_8401HR.Setup\_8401HR method), 61  
 StopStream() (Setup.SetupOneDevice.Setup\_8480SC.Setup\_8480SC method), 66  
 StopStream() (Setup.SetupOneDevice.Setup\_PodInterface.Setup\_Interface method), 70  
 Stream() (Setup.SetupOneDevice.Setup\_PodInterface.Setup\_Interface method), 70  
 STX() (PodApi.Packets.Packet.Packet static method), 41  
 syncConfig (PodApi.Parameters.Params8480SC.Params8480SC attribute), 49  
 SyncConfigBits() (PodApi.Devices.PodDevice\_8480SC.Pod8480SC static method), 24  
 systemID (PodApi.Parameters.Params8229.Params8229 attribute), 45  
 T  
 TranslateAll() (PodApi.Packets.Binary.PacketBinary method), 28  
 TranslateAll() (PodApi.Packets.Binary4.PacketBinary4 method), 31  
 TranslateAll() (PodApi.Packets.Binary5.PacketBinary5 method), 35  
 TranslateAll() (PodApi.Packets.Packet.Packet method), 41  
 TranslateAll() (PodApi.Packets.Standard.PacketStandard method), 43  
 TranslateBinaryTTLbyte() (PodApi.Packets.Binary4.PacketBinary4 static method), 31  
 ttl (PodApi.Packets.Binary4.PacketBinary4 attribute), 28  
 Ttl() (PodApi.Packets.Binary4.PacketBinary4 method), 31  
 TTLConfigBits() (PodApi.Devices.PodDevice\_8480SC.Pod8480SC static method), 25  
 ttlPullups (PodApi.Parameters.Params8480SC.Params8480SC attribute), 49  
 ttlSetup (PodApi.Parameters.Params8480SC.Params8480SC attribute), 49  
 TwosComplement() (PodApi.Packets.Packet.Packet static method), 41  
 U  
 U16() (PodApi.Commands.PodCommands.CommandSet static method), 4  
 U32() (PodApi.Commands.PodCommands.CommandSet static method), 4  
 U8() (PodApi.Commands.PodCommands.CommandSet static method), 4  
 UnpackAll() (PodApi.Packets.Binary.PacketBinary method), 28  
 UnpackAll() (PodApi.Packets.Binary4.PacketBinary4 method), 31  
 UnpackAll() (PodApi.Packets.Binary5.PacketBinary5 method), 35  
 UnpackAll() (PodApi.Packets.Packet.Packet method), 41  
 UnpackAll() (PodApi.Packets.Standard.PacketStandard method), 43  
 V  
 ValidateParams() (Setup.SetupOneDevice.Setup\_PodInterface.Setup\_Interface method), 70  
 V

## W

`week` (*PodApi.Parameters.Params8229.Params8229* attribute), [47](#)

`Write()` (*PodApi.Devices.SerialPorts.SerialComm.PortIO* method), [8](#)

`WritePacket()` (*PodApi.Devices.BasicPodProtocol.Pod* method), [11](#)

`WritePacket()` (*PodApi.Devices.PodDevice\_8229.Pod8229* method), [16](#)

`WritePacket()` (*PodApi.Devices.PodDevice\_8401HR.Pod8401HR* method), [22](#)

`WritePacket()` (*PodApi.Devices.PodDevice\_8480SC.Pod8480SC* method), [25](#)

`WriteRead()` (*PodApi.Devices.BasicPodProtocol.Pod* method), [11](#)